
# U.S. ARMY MEDICAL DEPARTMENT CENTER AND SCHOOL FORT SAM HOUSTON, TEXAS 78234-6100



**Hematology I** 

#### DEVELOPMENT

This subcourse is approved for resident and correspondence course instruction. It reflects the current thought of the Academy of Health Sciences and conforms to printed Department of the Army doctrine as closely as currently possible. Development and progress render such doctrine continuously subject to change.

#### **ADMINISTRATION**

For comments or questions regarding enrollment, student records, or shipments, contact the Nonresident Instruction Branch at DSN 471-5877, commercial (210) 221-5877, toll-free 1-800-344-2380; fax: 210-221-4012 or DSN 471-4012, e-mail accp@amedd.army.mil, or write to:

COMMANDER
AMEDDC&S
ATTN MCCS HSN
2105 11TH STREET SUITE 4192
FORT SAM HOUSTON TX 78234-5064

Approved students whose enrollments remain in good standing may apply to the Nonresident Instruction Branch for subsequent courses by telephone, letter, or e-mail.

Be sure your social security number is on all correspondence sent to the Academy of Health Sciences.

#### CLARIFICATION OF TRAINING LITERATURE TERMINOLOGY

When used in this publication, words such as "he," "him," "his," and "men" are intended to include both the masculine and feminine genders, unless specifically stated otherwise or when obvious in context.

.

## **TABLE OF CONTENTS**

Lessons			<u>Paragraphs</u>
	INTRODUCT	TION	
1	BLOOD		
	Section I. Section II. Section IV. Section V.	Normal Cell Maturation Abnormal Cell Maturation	1-1-1-4 1-5-1-6 1-7-1-9 1-10-1-12 1-13-1-16
	Exercises		
2	MATERIAL E	EMPLOYED IN HEMOTOLOGY	
	Section I. Section II. Section III.	Laboratory Glassware	2-1-2-3 2-4-2-7 2-8-2-9
	Exercises		
3	COLLECTIO BLOOD SME	N OF BLOOD AND PREPARATION OF EARS	
	Section I.	Collection of Blood Specimens	3-1-3-4
	Section II.	Preparation and Staining of Blood Smears	3-5-3-8
	Exercises		
4	MORPHOLO	GY OF BLOOD CELLS	
	Section I. Section II. Section IV. Section V.	General Information Erythrocytes Leukocytes Thrombocytes The Leukemias	4-1-4-2 4-3-4-5 4-6-4-11 4-12-4-13 4-14-4-18
	Exercises		

MD0853

## **TABLE OF CONTENTS (continued**

Les	<u>ssons</u>		<u>Paragraphs</u>
5	MANUAL CELI	COUNTS	
	Section I. Section II	Manual Counts, Blood Specimens Manual Counts, Other Body Fluids	5-1-5-4 5-5-5-6
	Exercises		
6		RIT, ERYTHOCYTE SEDIMENTATION HEMOGLOBIN	
	Section I. Section II. Section III.	Erythrocyte Sedimentation Rate	6-1-6-2 6-3-6-5 6-6-6-14
	Exercises		
	APPENDIX		

MD0853 ii

## CORRESPONDENCE COURSE OF THE U.S. ARMY MEDICAL DEPARTMENT CENTER AND SCHOOL

#### **SUBCOURSE MD0853**

#### **HEMATOLOGY I**

#### INTRODUCTION

This subcourse is concerned with the blood tests performed in the hematology section of the laboratory. The purpose of these tests is to aid the physician in diagnosis. Thus, these tests are important and often essential to the health and life of the patient. Thorough study of this subcourse should enable you to better fulfill your role in health care.

#### Subcourse Components:

This subcourse consists of 6 lessons and an Appendix. The lessons are:

Lesson 1. Blood

Lesson 2. Material Employed in Hematology.

Lesson 3. Collection of Blood and Preparation of Blood Smears.

Lesson 4. Morphology of Blood Cells.

Lesson 5. Manual Cell Counts.

Lesson 6. Hematocrit, Erythrocyte Sedimentation Rate, and Hemoglobin.

Appendix. Glossary of Terms.

#### Credit Awarded:

To receive credit hours, you must be officially enrolled and complete an examination furnished by the Nonresident Instruction Branch at Fort Sam Houston, Texas. Upon successful completion of the examination for this subcourse, you will be awarded 12 credit hours.

You can enroll by going to the web site <a href="http://atrrs.army.mil">http://atrrs.army.mil</a> and enrolling under "Self Development" (School Code 555).

A listing of correspondence courses and subcourses available through the Nonresident Instruction Section is found in Chapter 4 of DA Pamphlet 350-59, Army Correspondence Course Program Catalog. The DA PAM is available at the following website: http://www.usapa.army.mil/pdffiles/p350-59.pdf.

MD0853 iii

#### **LESSON ASSIGNMENT**

LESSON 1 Blood.

**TEXT ASSIGNMENT** Paragraphs 1-1 through 1-16.

**LESSON OBJECTIVES** After completing this lesson, you should be able to:

1-1. Select the statement that best describes the composition of blood.

1-2. Select the statement that best describes the cellular and liquid constituent characteristics of blood.

1-3. Select the characteristic features of normal and abnormal blood cell formation and maturation.

1-4. Correctly select the functions of blood cells.

**SUGGESTION** After completing the assignment, complete the

exercises at the end of this lesson. These exercises

will help you to achieve the lesson objectives.

#### LESSON 1

#### **BLOOD**

#### Section 1. THE COMPOSITION OF BLOOD

#### 1-1. INTRODUCTION

Blood is a complex and unique fluid of variable composition circulating through the heart, arteries, capillaries, and veins, known as the vascular system of the body. It is a tissue in which cellular constituents are suspended in a liquid medium performing specialized functions. The prime function of blood is to carry oxygen from the lungs to the body tissues and carbon dioxide from the tissues to the lungs. Blood carries fluid to and from the tissues, thus maintaining proper fluid balance throughout the body (explained later in this paragraph). The average pH value of blood is 7.40. Blood also carries nutrients or food supplies from the digestive system to the body cells or tissues and transports waste products for the tissues to the kidneys and bowel for excretion to prevent accumulation. In response to trauma or infection, blood cells and antibodies are carried in the blood to a point protection against the causative agents of disease, or to transport blood-clotting substances to a break in a blood vessel to promote the clotting process when injury is caused by bleeding or hemorrhage. The blood also carries hormones from the endrocrine glands to the target organs, and it assists in the regulation of the body temperature by carrying excess heat from the interior of the body to the surface layers of the skin, where the heat is dissipated to the surrounding air. To perform these complex functions, blood is complex and is composed of two main parts. One part is composed of blood cells, the particles suspended in the plasma, making up approximately 45 percent of total blood volume and including erythrocytes (red blood cells), leukocytes (white blood cells), and thrombocytes (platelets). The red and white blood cells are known as corpuscles. The other part is plasma, the fluid portion of the blood, which consists primarily of water in which are dissolved proteins and many inorganic and organic substances carried by the blood to and from the tissues. Plasma makes up 55 percent of the blood.

#### 1-2. GENERAL CELLULAR STRUCTURE

- a. A typical cell is composed of a single nucleus embedded in cytoplasm. The living substance of the cell is a grayish, viscous liquid called protoplasm. Protoplasm is enclosed in the cell interfaces by a cell membrane that selectively regulates the interchange of materials between the cell and its environment. A typical cell tree is diagramed in figure 1-1, however, all cells and forms or development of cells are not depicted.
- b. The nucleus is a spherical oval body surrounded by a thin membrane (nuclear membrane). Contained in the nucleus is a sphere called the nucleolus. The nucleus is thought to be an organizing center for the cell and can have the capacity for cell production. The absence of nuclei signifies the end of cell development. Also found




Figure 1-1. Typical cell structure.

within the nucleus is a network of nuclear fibrils made up of DNA (deoxyribonucleic acid) and protein called chromatin. It is thought that the decreasing growth activity of a cell during maturation is regulated by chromatin.

c. Surrounding the nucleus is a mass of protoplasm called cytoplasm. Contained within the cytoplasm are numerous granules, filaments, and globules. These structures are divided into two groups known as organoids (organelles) and inclusions. The organoids are thought to perform most of the metabolic functions of the cell. Mitochondria, Golgi apparatus, fibrils, centrioles, and the chromatin substance are classified as organelles. Cytoplasm inclusions are usually seen as granulation. The granulation is an accumulation of proteins, lipids, carbohydrates, pigments, and secretory granules.

#### 1-3. CELLULAR CONSTITUENTS

a. **Erythrocytes**. An erythrocyte (red blood cell) is an elastic, non-nucleated, biconcave disc having a diameter of approximately 7.2 microns. The mature red cell contains about 34 percent hemoglobin (a complex iron-bearing pigment that transports oxygen). Hemoglobin is contained in the interior of the cell, and the outer surface of the cell is surrounded by a cell membrane. When unstained, the cell has a pale, greenish-yellowish appearance. It is buff pink with an accented central zone of pallor when stained with Wright's stain. The production of erythrocytes or erythropoiesis, occur primarily in the red marrow of the spongy bones. Erythrocytes make up the great majority of cells found in the peripheral blood. Their vast surface area is important in the transport of oxygen from the lungs to the tissues because of quick exchange of oxygen in both sites that occurs across the red cell surface. Erythrocytes are subject to

many alterations in shape, size, staining properties, and structure in different disease processes. An adult female has approximately 4.8 million/cu mm red cells, and an adult male has approximately 5.4 million/cu mm red cells. The erythrocytes have an average life span of 80 to 120 days. See Table 1-1 for purposes of the blood cell tree as described in this subcourse.

Plasmocytic series Leukocytes Plasmoblast Granulocytes series Proplasmocyte Myeloblast Promycelocyte Plasmocyte Neutrophils Thrombocytes Neutrophilic myelocyte Megakarycotic series Neutrophilic metayelocyte Neutrophilic band Megakaryoblast Neutrophil segmented cell Megakaryocyte Thrombocyte Eosinophils Basophilic myelocytes Metamyelocyte bands Aragraulocytes Seamented forms Basophils Erythrocyte series Rubriblast Lymphocyte series Lymphoblasts Prorubricyte Prolymphocytes Rubricyte Metarubricyte Lymphocytes Monocyte series Reticulocyte Monoblasts Erythrocyte Promonocytes Monocytes

Table 1-1. Blood cell tree.

- b. **Leukocytes**. Leukocytes are commonly known as white blood cells because of their lack of color in unstained preparations. They are nucleated cells that have an average diameter of 8 to 12 microns. There are 5,000-10,000 leukocytes per cubic millimeter of blood. Leukocytes are divided into three groups: (1) granulocytes (neutrophils, eosinophils, and basohils) that can phagocytose or ingest bacteria or other particles; (2) lymphocytes which participate in humoral (B lymphocyte) and cell mediated (T lymphocyte) immunity and (3) monocytes that phagocytose bacteria cellular debris and interact with lymphocytes in the processing of antigens in the immune reaction. In humans, there are 3 types of phagocytes: macrophages, granulocytes, and monocytes. All three types are attracted to bacteria and certain other particles by substances released by the particles (the process of chemotaxis). They then pin the particle against a surface and engulf it. They are differentiated by the specific nuclear and cytoplasmic staining properties.
- c. **Thrombocytes**. Thrombocytes are found commonly known as platelets. They are the morphologically recognizable precursors of platelets and are detached

fragments of the cytoplasm-megakaryocytes that are found in the bone marrow (others coming from red bone marrow include: promegakarycocytes and megakaryoblasts). Unstained platelets appear as small hyaline structures with a diameter of approximately two microns. When stained with Wright's stain, they have a pale blue cytoplasm with a dark granular center. Platelets are very fragile and live for a period of about 3 to 5 days.

#### 1-4. LIQUID CONSTITUENTS

- a. **Plasma**. Blood plasma is a clear, yellowish fluid that accounts for about 55 percent of the total volume of the blood. The chemical nature of plasma is very complex. It consists of 92 percent water, 7 percent proteins (albumin, globulin, and fibrinogen), carbohydrates (glucose), lipids (fats, lecithin, and cholesterol), dissolved gases (oxygen, carbon dioxide, and nitrogen), non-protein nitrogenous substances, and less than 1 percent of inorganic salts. Blood plasma is the fluid portion of the blood before clotting occurs. Put another way, plasma from which fibrinogen has been removed is called serum.
- b. **Serum**. Serum is the fluid portion of blood after the clotting process is complete. Fibrinogen, the precursor of fibrin, is removed from plasma to form the framework of the blood clot.

#### Section II. FORMATION OF BLOOD CELLS

#### 1-5. EMBRYONIC HEMATOPOIESIS

a. The primary source of blood cells is the mesenchyme connective tissue in the embryo. Three phases of embryonic hematopoiesis merge, resulting in the formation of blood (see figure 1-2).




Figure 1-2. Embryonic hematophoiesis.

b.

- b. During the first two months of embryonic development, the mesoblastic phase occurs. The blood cells are formed in the blood islands of the yolk sac. Immature blood cells develop, having large nuclei containing a vesicular chromatin meshwork. The nuclei are surrounded by a thin rim of cytoplasm. The cells gradually develop hemoglobin to become a nucleated red blood cell. Mitosis occurs and daughter cells contain more hemoglobin. Finally the nucleus is lost and the nature erythrocyte is produced.
- c. At two months, the hepatic phase begins. Blood cell development shifts to the body of the fetus as the organs of the reticuloendothelial system (liver, spleen, thymus, etc.) are formed. During this phase, red cells begin their normal development. Granulocytes and thrombocytes are formed in the liver while lymphocytes and monocytes are formed in the spleen and thymus.
- d. The myeloid phase begins during the fifth month of gestation. Blood cells are formed in the bone marrow and lymphatic system that at the time of birth constitute the total sources of hematopoiesis. The bone marrow is the principle source of production of erythrocytes, granulocytes, and thrombocytes. The lymph nodes are the primary sites of production of lymphocytes, monocytes, and plasmacytes.

#### 1-6. POSTNATAL HEMATOPOIESIS

- a. Blood formation at birth is confined primarily to the bone marrow (central medullary structure of the bone). Blood cells multiply by mitosis and then mature to a specific cell type. The mature cells lose the ability to reproduce and develop a definite life span. Regeneration of blood cells after birth involves multiplication of precursor cells, evolution of the definitive characteristics of each type, and release of mature cells.
- b. Myelopoiesis is the production of blood cells and bone marrow by the bone marrow (medullary site of production). The red bone marrow is the principle source of production of red cells and white cells of the granulocytic series. At birth, the central medullary structure of bones is red bone marrow and it is actively engaged in hematopoiesis. At about 5 years of age non-hematopoietic type of marrow (fatty yellow bone marrow) that is a reserve potential tends to replace most of the red bone marrow. This partial replacement of red marrow is complete when the individual reaches maturity (about 18 years of age) at which time, active hematopoietic centers in bone tissue are limited almost exclusively to the sternum, pelvic area, vertebrae, skull, ribs, clavicle, scapuli, and the epiphyses of the long bones.
- c. Extramedullary hematopoiesis is blood production that occurs in sites other than the bone marrow. Active sites are the spleen, thymus, lymph nodes, and other lymphoid tissues. Cell production is largely limited to lymphopoiesis. The lymph nodes are the primary source of lymphocytes and plasmacytes (see figure 1-3).




Figure 1-3. Hematopoietic system in an adult 18-20 years of age.

#### Section III. NORMAL CELL MATURATION

#### 1-7. GENERAL FEATURES

In the course of blood cell maturation, certain specific features are developed (see Figure 1-4). Each of the component parts of the cell undergoes a transformation during maturation. The immature cell or blast cell contains a large nucleus, a small amount of cytoplasm, and no granules. As the cell ages, the cytoplasm becomes less basophilic and nuclear chromatin becomes heavier (darker stain). Reduction in size and loss of nucleoli occurs as the cell becomes older. The three types of granulation (neutrophilic, basophilic, and eosinophilic) become more specific and smaller as the cell ages. Maturation, in general, involves: (1) cytoplasmic differentiation, (2) nuclear maturation, and (3) reduction in cell size.




Figure 1-4. Development of blood cells.

#### 1-8. CYTOPLASM

The basophilia of a blast cell is proportional to the ribonucleic acid (RNA) content. As the cell matures, the RNA content decreases and the cell becomes a paler blue. In the myeloid cells, a specific type of granulation occurs. When granules appear they are pinkish-red and few in number. The granules increase in number and differentiate into three types upon maturation. As the cell matures, it develops an affinity for the acid or basic portion of the stain (Wright's stain). Basophilic granules stain blue-black, eosinophilic granules stain red-orange, and neutrophilic granules stain pinkish-purple. Lymphocytes are usually devoid of cytoplasmic granulation but they can possess non-specific azurophilic (dark-purple) granules, usually characteristic of monocytes and plasmocytes. Upon maturation, the erythrocyte develops a light orange respiratory pigment called hemoglobin.

#### 1-9. NUCLEUS

The nucleus of the young cell is large, round, and occupies most of the cell. As the cell matures, the size of the nucleus decreases. Nuclei of early or primitive cells usually have one or more nucleoli. The latter are small, round, homogeneous areas that usually stain light blue with a darker boundary. In appearance nucleoli are somewhat like craters in the nucleus. They are surrounded by strands of chromatin. These nucleoli, plus a delicate reticular network of chromatin, are the principle indicators of blood cell immaturity. As the cells mature, the nucleus gradually becomes smaller, stains darker, and chromatin meshwork become "coarse" with the strands of chromatin less fine and lacelike. In the course of cell development, the nucleus changes its shape, particularly in the granulocytic series, where it becomes indented, lobulated, segmented, or fragmented. As maturation or development progresses, the nucleus, if still intact, becomes small, compact, usually dark and structureless, and can completely disappear. The loss or shrinking of the nucleus is accompanied by a decrease in cell size.

#### Section IV. ABNORMAL CELL MATURATION

#### 1-10. GENERAL FEATURES

Abnormal cell maturation is asynchronous development as opposed to normal cell maturation or synchronous development. Since the normal sequence of cell maturation is upset, atypical cells will be present. Abnormal cells can be recognized by: (1) abnormal cytoplasmic maturation, (2) abnormal nuclear maturation, and (3) abnormal size.

#### 1-11. CYTOPLASM

Asynchronism of the cytoplasm is most commonly seen in the granulocytes. Granulation can be primitive or absent. In some instances, the granules fail to differentiate. Erythrocytes show basophilia late in the series and retarded

hemoglobination. Inclusions in the cytoplasm, such as Dohle bodies (infectious diseases), Auer rods (leukemia), and toxic granulation (infection affecting the marrow) are seen in the abnormal white cells.

#### 1-12. NECLEUS

Abnormal cells often show two nuclei in severe disturbances, such as leukemia. Nucleoli have a retarded reduction. The nucleus can have an irregular outline or indentation (Rieder cells). Hypersegmented nuclei occur in the neutrophils in sepsis and in pernicious anemia. Abnormal maturation of the nucleus often results in variation in cell size.

#### Section V. FUNCTIONS OF BLOOD CELLS

#### 1-13. INTRODUCTION

Each component of blood is uniquely capable of performing one or more functions. Together, these components provide the maintenance of a relatively stable environment of the body by a variety of mechanism. This maintenance of a relative biological constancy or integrity is known as homeostasis. Once the blood cells reach full maturity, they enter the bloodstream and begin fulfilling their functions.

#### 1-14. ERYTHROCYTES

Hemaglobin is the main functioning component of the cell. It carries out the transportation of oxygen to the tissues and the removal of carbon dioxide. Hemoglobin also aids in the maintenance of the delicate acid-based buffer system of the body. The erythrocyte must also supply energy to accomplish the active transport of glucose and ions against a gradient across the red cell membrane.

#### 1-15. LEUKOCYTES

Leukocytes remove invading antigens (for example, bacteria) and to some extent transport and distribute antibodies. Monocytes and all of the granulocytes have been shown to demonstrate directional movement. Their movement is subject to chemotaxis or the response of living protoplasm to a chemical stimulus. As mentioned previously, the attraction of the leukocytes to bacteria and certain other particles by substances released by the particles is called chemotaxis. They either transport the particles or engulf them. The chemotaxis process of engulfing and destroying bacteria, or phagocytosis, is a prime function of leukocytes.

a. **Monocytes**. These cells will engulf bacteria and larger materials, including even protozoa and red cells, and are transformed into and/or are called macrophages. In this regard, monocytes are perhaps the most efficient phagocytes of all the cells. Monocytes contain many of the lytic enzymes that are found in microphages (granulocytes). In addition, monocytes contain lipases that enable them to dissolve the lipoid capsules of certain bacteria.

- b. **Neutrophils**. Neutrophilic leukocytes are excellent microphages. That is, they engulf bacteria and other microscopic particles. The particles are first surrounded by cellular pseudopodia and then incorporated into a cell vacuole. There the foreign bodies mix with substances released from the cytoplasm of leukocytes. In this way the leukocyte is not injured by whatever "combat activity" is taking place in the vacuole. Neutrophils are fully developed (mature) cells that are incapable of mitotic division. They carry on active metabolism. Eventually the granulocytes disintegrate and in inflammatory processes are succeeded by monocytes.
- c. **Eosinophils**. The eosinophils are easily distinguishable by their large cytoplasmic granules. They have some similarities with the neutrophil: they arise from a common stem cell and share same morphologic features. Unlike the neutrophil, however, eosinophils initially develop in the bone marrow, then are released into the circulation, and undergo further maturation in the spleen. Their emergence from the marrow seems to depend upon a stimulus provided by lymphocytes. Eosinophils are found in tissue fluid as well as in peripheral blood, especially in areas where there is an allergic reaction. Current thinking holds that eosinophils are involved in antigenantibody reactions, and have been shown to phagocytize antigen-antibody reactants. Eosinophils are also thought to transport, or at least contain, lysins that act on fibrin. It is suggested that eosinophils limit the action of substances such as histamine. How this is accomplished is not yet clear. The mobilization of eosinophils from their reserve in the bone marrow is at least in part under hormonal control. If the adrenal cortex is functioning properly, an injection of adrenocorticotropic hormone (ACTH) results in a marked decrease in the number of circulating eosinophils and in the number of circulating lymphocytes. On the other hand, there is an increase in the number of circulating neutrophils.
- d. **Basophils**. The function of basophils in man has not been ascertained. They quite possibly represent a vestige of evolution. Their granules have been found to contain heparin and these cells frequently appear during the clot dissolution phase of an injury. Hence, it has been suggested they may be involved in clot absorption.
- e. **Lymphocytes**. The lymphocyte is now believed to be directly connected with antibody production. Undoubtedly, the lymphocyte performs important immunologic functions. According to very recent studies, many *of* the activities previously thought to take place in the reticuloendothelial (RE) system actually take place in lymphocytic tissue.

#### 1-16. PLATELETS

Platelets possess metabolic system, expend energy, and respond to stimulus. They contain many enzymes and undergo respiratory activity and glycolysis. They possess coagulation factors usually designated as PF-1, PF-2, and on through PF-7. The cells contain fibrinogen and vasoconstrictor substances, calcium, and many other components that are either known or presumed to participate in the clotting mechanism. Clot-promoting lipoproteins are also found in platelets. In addition, well-defined

antigens have been found in platelets. The role of platelets in the blood coagulation mechanism will be described in more detail in Subcourse MD0857 Lesson 2.

**NOTE**: See Table 1-2 for a summary of blood cell functions.

	보이는 그는 그는 그는 그는 그들은 그는 그들은 사람들이 얼마나 되었다. 그는 그는 그들은 그는 그는 그는 그는 그는 그는 그는 그는 그를 가는 것이 없다면 그는 그를 가는 것이 없다면 하다.
BLOOD CELL NAME	FUNCTION
<u>Erthrocytes</u>	Contain hemoglobin, transport ${\rm O}_2$ and ${\rm CO}_2$ , acid-base buffer, and supply energy and ions.
<u>Leukocytes</u>	Remove invading antigens and transport and distribute antibodies.
Monocytes	Phagocyte cells surround, engulf, and diges bacteria or other particles. Disolve lipoi capsules of certain bacteria.
Neutrophils	Engulf bacteria and other microscopic particles and perform "combat activity" in vacuoles. Carry on active metabolism.
Esoinophils	Involved in and phagocytize antigen-antibody reactants. Possibly transport or contain lysins and may limit action of substances like histamines
Basophils	Involved with possible clot absorption.
Lymphocytes	Directly connected with antibody production and perform immunologic functions.
<u>Platelets</u>	Expend energy, respond to stimuli, undergo respiratory activity and glycolysis, and participate in the clotting mechanism.

Table 1-2. Blood cell functions.

**Continue with Exercises** 

## **EXERCISES, LESSON 1**


**INSTRUCTIONS**: Answer the following exercises by marking the lettered response that best answers the exercise, by completing the incomplete statement, or by writing the answer in the space provided at the end of the exercise. Some exercises may require you to refer to the appendices located at the back of the subcourse.

After you have completed all the exercises, turn to "Solutions to Exercises" at the end of the lesson and check your answers. For each exercise answered incorrectly, reread the material referenced with the solution.

1.	Plasma constitutes approximately percentage of the bland is composed mainly of		
	a.	45 percent; water.	
	b.	50 percent; protoplasm.	
	C.	55 percent; water.	
	d.	70 percent; cytoplasm.	
2.	Wł	nat does chromatin regulate in a cell during maturatio	n?
	a.	Movement.	
	b.	Growth.	
	C.	Granulation.	
	d.	Waste products.	
3.	Су	toplasmic granulation is usually due to cytoplasmic:	
	a.	Fibrils.	
	b.	Inclusions.	
	C.	Organelles.	
	d.	Parachromatin.	

- 4. When stained with Wright's stain, what is the color of the normal erythrocyte and how is it accented?
  - a. Light orange *I* granular zone.
  - b. Pale green, central pale zone.
  - c. Pale, greenish-yellowish, central pale zone.
  - d. Buff pink, central zone of pallor.

INSTRUCTIONS: Exercise items 5 through 10 require you to select the name of the particle, activity or definition of the typical cell structure as shown below.



Typical cell structure.

- 5. A grayish, viscous liquid living substance within the cell that regulates the interchange of materials between the cell and the environment is called:
  - a. Cytoplasm.
  - b. Protoplasm.
  - c. Centriole.
  - d. Mitochondria.

	C.	Chromatin.
	d.	Mitochondria.
7.		nich chromatin substances are contained within cytoplasm and are classified as anelles?
	a.	All granules, filaments, and globules inside the nucleus.
	b.	Mitochondria, Golgi apparatus, fibrils, and centrioles are some of them.
	C.	Some granules, filaments, and globules inside the nucleus.
8.	Wh	ich organoid substance appears like a cork-screw cigar?
	a.	Golgi apparatus.
	b.	Fibrils.
	C.	Centrioles.
	d.	Mitochondria.
9.	Mit	ochondria is defined in this subcourse as:
	a.	A meshwork of lipid containing fibrils within the cytoplasmic portion of a cell.
	b.	A spherical or oval body surrounded by a thin membrane (nuclear membrane)
	C.	Granular components of a cell cytoplasm active in oxidative processes
	d.	A minute cell organoid within the centrosome.

1-15

6. The nuclear fibrils that possibly regulate growth activity during maturation and are made up of DNA (deoxyribonucleic acid) and protein are called:

a. Lysosome.

MD0853

b. Golgi apparatus.

10.	Се	ntriole is defined in this subcourse as:
	a.	A meshwork of lipid containing fibrils within the cytoplasmic portion of a cell.
	b.	A spherical or oval body surrounded by a thin membrane (nuclear membrane)
	C.	Granular components of a cell cytoplasm active in oxidative processes
	d.	A minute cell organoid within the centrosome.
11.	On	ly nucleated red blood cells are formed in the islands of the yolk sac during the of embryonic development.
	a.	First 2 months.
	b.	Third to fifth months.
	C.	Fifth to final months.
	d.	Entire duration.
12.		e bone marrow is the principal organ of red blood cell production during the eloid phase as well as at birth and during the remainder of life.
	a.	First 2 months.
	b.	Third to fifth months.
	C.	Fifth to final months.
	d.	Entire duration.
13.	As	the child matures, myelopoiesis is increasingly confined to the:
	a.	Bones.
	b.	Medullary sites.
	C.	Ligaments.
	d.	Muscles.

14.	Wh	Which definition of myelopoeisis is correct?		
	a.	It is the production of blood cells and bone marrow by the kidney.		
	b.	It is the production of blood cells and bone marrow by the liver.		
	C.	It is the production of blood cells and bone marrow by the bone marrow (medullary site of production).		
	d.	A decrease in the number of neutrophils in the blood.		
15.	Du	ring the formation of blood cells, which tissue is their primary source?		
	a.	Striated tissue.		
	b.	Mesenchyme connective tissue in the embryo.		
	C.	Liver.		
	d.	Tissue thromboplastin.		
16.		e three types of granulation occurring during normal cell maturation that come more specific and smaller as the cell matures are:		
	a.	Basophilic, hematopoietic, and eosinophilic.		
	b.	Non-hematopietic, neutrophilic, and eosinophilic.		
	C.	Netrophilic, basophilic, an eosinophilic.		
	d.	Basophllic, neutrophilic, and idiopathic.		
17.		mature blood cells that are so young that they can hardly be distinguished orphologically from each other are called:		
	a.	First cells.		
	b.	Blast cells.		
	C.	Band cells.		

18.	Wł	nat occurs to the RNA content of the blast cell as it matures?
	a.	Increases.
	b.	Decreases.
	C.	All of the above.
19.	Ge	nerally speaking, as normal cells maturate:
	a.	The nucleus decreases in size.
	b.	The cytoplasm differentiates.
	C.	The cell reduces in size.
	d.	All of the above.
20.	Wł	nat color does the blast cell develop into as it matures?
	a.	Green.
	b.	Dark purple.
	C.	Blue.
	d.	Pinkish-purple.
21.		mphocytes during normal cell maturation are devoid of cytoplasmic granulation t can possess nonspecific granules.
	a.	Azurophilic.
	b.	Dark green.
	C.	Light pink.
	d.	Rich orange.

22.		at colors do basophilic and neutrophilic grant blied?	ules appear when Wright's stain is
	a.	Pink or brown-orange.	
	b.	Red-orange; blue-black.	
	C.	Blue-black; pinkish-purple.	
	d.	Reddish-purple; blue-black.	
23.	Wit	h Wright's stain, eosinophilic granules appea	r large and:
	a.	Pink or brown.	
	b.	Red-orange.	
	C.	Blue-black.	
	d.	Reddish-purple.	
24.	Wit	h Wright's stain, the nonspecific azurophilic g	granule is:
	a.	Pink or brown.	
	b.	Red-orange.	
	C.	Colorless.	
	d.	Red/dark purple.	
25.	The sign	e presence of n of a cell	shaped nucleoli is generally a
	a.	Oblong; maturity.	
	b.	Crater; immaturity.	
	C.	Lacelike; cytoplasmic change.	
	d.	Coarse: granulocytic change.	

26.	Abnormal cytoplasmic maturation is seen most frequently in:		
	a.	Monocytes.	
	b.	Lymphocytes.	
	C.	Granulocytes.	
	d.	Erythrocytes.	
27.	Wh	at are Dohle bodies?	
	a.	Greek food.	
	b.	Atypical cells.	
	C.	Abnormal cells.	
	d.	Infectious diseases.	
28.	Wh	at another term for normal cell maturation?	
	a.	Polychromasic.	
	b.	Synchronous.	
	C.	Asynchronous.	
	d.	Rieder.	
29.		at term is used to describe abnormal cell maturation or uncoordinated cell relopment?	
	a.	Dohle bodies.	
	b.	Asynchronous.	
	C.	Binary fission.	
	d.	Synchronous.	
	e.	Azurophilic.	
	f.	Sepsis.	

30.	30. If asynchronism occurs in the cytoplasm, where would it most commonly be se		
	a.	Granulocytes.	
	b.	Nucleolus.	
	C.	Lysosome.	
	d.	Mitochondria.	
31.	Wh	en there is abnormal cell maturation, where does inclusion occur?	
	a.	Cell membrane.	
	b.	Nucleus.	
	C.	Cytoplasm.	
	d.	Fat droplet.	
32		ich cytoplasmic cells show basophilia late in the series and hemoglobinization ardation during abnormal cells maturation?	
	a.	Monocytes.	
	b.	Lymphocytes.	
	C.	Granulocytes.	
	d.	Erythrocytes.	
	e.	Neutrophils.	
	f.	Basophils.	

33.	Leu	ıkemia is also known as:
	a.	Auer rods.
	b.	Rieder cells.
	C.	Juvenile cells.
	d.	Macrocytosis.
	e.	Aplastic anemia.
34.		en a cell is maturing abnormally, how many nuclei might it have if it has a ere disturbance such as leukemia?
	a.	0.
	b.	1.
	C.	2.
	d.	3.
35.	The	e nuclei of abnormal cells often have:
	a.	Two nuclei, in severe disturbances, as in leukemia.
	b.	Reduced reduction of nucleoli.
	C.	Undifferentiating granules.
	d.	Rieder cell indentation.
	e.	Hypersegmented nuclei.
	f.	Cell variation.
	g.	All of the above.
	h.	None of the above.

36.	-	posegmented nuclei occur in the neutrophils in and in rnicious		
	a.	Sepsis; leukemia.		
	b.	Sepsis; anemia.		
	C.	Sepsis; hemotoma.		
	d.	Sepsis; sickle cell anemia.		
37.	Homeostasis is defined as:			
	a.	Biological constancy or integrity of blood or the checking of the flow of blood especially from a vessel.		
	b.	The spitting of blood; coughing up blood.		
	C.	Infection of the marrow; infectious disease.		
	d.	Blood cells reaching full maturity and entering the mainstream to begin their functions.		
38.	Which statement is correct concerning the function of the blood cell?			
	a.	Individually, each component of blood provides the maintenance for a relatively stable environment of the body by a variety of mechanisms.		
	b.	Together, the blood components provide the maintenance for a relatively stable environment of the body by a variety of mechanisms.		
39.	One function of the erythrocyte is to aid in the maintenance of the delicate acid- base buffer system of the body. Which part of the cell performs this function?			
	a.	Hemoglobin.		
	b.	Membrane.		
	C.	Lysosome.		
	d.	Glycoproteins.		

40.	Wh	nich are the oxygen-carrying cells of the body?
	a.	Monocytes.
	b.	Granulocytes.
	C.	Lymphocytes.
	d.	Erythrocytes.
41.		mainly controls the level of oxygen in the tissues.
	a.	Hemoglobin.
	b.	Chemotaxis.
	C.	Phagocytosis.
	d.	Adrenocorticotropic hormone.
42.	Ery	ythrocytes, leukocytes, and platelets are considered to be the of (in) the blood.
	a.	Coagulation factors.
	b.	Immunologic functions.
	C.	Components.
43.	Wh	nich cells possess the powers of locomotion and phagocytosis?
	a.	Erythrocytes.
	b.	Leukocytes.
	C.	Thrombocytes.
	d.	Platelets.

44.	Which cell is able to phagocytize the largest particles and because of this is called macrophage?			
	a.	Neutrophils.		
	b.	Eosinophils.		
	C.	Basophils.		
	d.	Monocytes.		
45.	Which cell has the main function of transporting O2 to the tissues and removing CO2?			
	a.	Neutrophils.		
	b.	Eosinophils.		
	C.	Basophils.		
	d.	Erythrocyte.		
46.	Which cell is surrounded by cellular pseudopodia and is not injured by "combat activity" taking place in the vacuole?			
	a.	Eosinophils.		
	b.	Neutrophils.		
	C.	Monocytes.		
	d.	Basophils.		
47.	Which cell is incapable of mitotic division?			
	a.	Eosinophils.		
	b.	Neutrophils.		
	C.	Monocytes.		
	d.	Basophils.		

#### 48. Which statement is correct?

- Monocytes are cells containing foreign bodies with substances released from the cytoplasm.
- b. Lymphocytes are now believed to contain heparin.
- c. Easinophils are involved with antigen-antibody reactions and eat antigenantibody reactants.
- d. Basophils may be involved in attracting cells to substances and then either transport or engulf them.

#### 49. Cellular particles associated with the clotting of blood are:

- a. Thrombocytes (platelets).
- b. Hemoglobin particles.
- c. White blood cells.
- d. Red blood cells.

#### 50. Which statement is correct?

- a. Platelets possess metabolic systems, expend energy, respond to stimuli, and then disintegrate.
- b. Platelets possess metabolic systems, mobilize the bone marrow, expend energy, and respond to stimuli.
- c. Platelets possess carbon dioxide and metabolic systems, expend energy, and respond to stimuli.
- d. Platelets possess metabolic systems, expend energy, respond to stimuli, and undergo glycolysis.

**Check Your Answers on Next Page** 

## **SOLUTIONS TO EXERCISES, LESSON 1**

- 1. c (para 1-1)
- 2. b (para 1-2b)
- 3. b (para 1-2c)
- 4. d (para 1-3a)
- 5. b (para 1-2a)
- 6. c (para 1-2b)
- 7. b (para 1-2c)
- 8. d (figure 1-1)
- 9. c (appendix)
- 10. d (appendix)
- 11. a (para 1-5b)
- 12. c (para 1-5d)
- 13. a (para 1-6b)
- 14. c (para 1-6b)
- 15. b (para 1-5a)
- 16. c (para 1-7)
- 17. b (para 1-7)
- 18. b (para 1-8)
- 19. d (para 1-7)
- 20. c (para 1-8)
- 21. a (para 1-8)
- 22. c (para 1-8; figure 1-4; appendix)

- 23. b (para 1-8; figure 1-4; appendix)
- 24. d (para 1-8; appendix)
- 25. b (para 1-9)
- 26. c (para 1-11)
- 27. d (para 1-11)
- 28. b (para 1-10)
- 29. b (para 1-10)
- 30. a (para 1-11)
- 31. c (para 1-11)
- 32. d (para 1-11)
- 33. a (para 1-11)
- 34. c (para 1-12)
- 35. g (para 1-11, 1-12)
- 36. b (para 1-12)
- 37. a (para 1-13)
- 38. b (para 1-13)
- 39. a (para 1-14)
- 40. d (para 1-14)
- 41. a (para 1-14)
- 42. c (paras 1-13 to 1-16)
- 43. b (para 1-15)
- 44. d (para 1-15a)
- 45. d (para 1-14)

- 46. b (para 1-15b)
- 47. b (para 1-15b)
- 48. c (para 1-15c)
- 49. a (para 1-16)
- 50. d (para 1-16)

End of Lesson 1

#### **LESSON ASSIGNMENT**

LESSON 2

Material Employed in Hematology.

**TEXT ASSIGNMENT** 

Paragraphs 2-1 through 2-9.

**LESSON OBJECTIVES** 

After completing this lesson, you should be able to:

- 2-1. Correctly list the procedures in preparing laboratory reagents.
- 2-2. Select the safety precautions necessary to ensure that reagents and equipment are labeled and handled carefully and appropriately.
- 2-3. Select the statement that best describes important characteristics and functions of common hematological glassware.
- 2-4. Correctly identify the parts of a compound microscope.
- 2-5. Correctly select the functions of a compound microscope.
- 2-6. Select the statement that correctly describes the important characteristics and functions of equipment used in the laboratory.
- 2-7. Select the correct procedures for cleaning and maintaining laboratory equipment.

SUGGESTION

After completing the assignment, complete the exercises at the end of this lesson. These exercises will help you to achieve the lesson objectives.

MD0853 2-1

#### LESSON 2

#### MATERIAL EMPLOYED IN HEMATOLOGY

#### Section I. LABORATORY REAGENTS

#### 2-1. PREPARATION

Various stains and solutions are utilized in routine hematological examinations. These stains and solutions must be prepared with the utmost care and precisely according to formulations. Detailed directions for the preparation of all reagents that are required for performing procedures outlined throughout this subcourse are contained in the respective procedure. Careful attention should be given to precise measurements order in which reagents are added, control of temperature where indicated, filtration, and aging. Particular attention must be given to storage of reagents particularly with reference to requirements for refrigeration, incubation, and protection from intense light.

#### 2-2. LABELING REAGENT CONTAINERS

Proper labeling of reagents is an extremely important detail. Labels should be complete, securely attached, and neatly and legibly written or preferably typewritten. Items recorded on the label should include all constituents and quantities utilized, date of preparation, initials of the individual who prepared the reagent, and expiration date if the solution deteriorates with age. Labels should be protected against damage from water or other fluids by covering with a protective coating of cellophane tape over the surface of the label.

#### 2-3. SAFETY PRECAUTIONS

There are various precautions that must be taken in handling reagents in the hematology laboratory. Among the most important are the following:

- a. Once a portion of a reagent has been removed from the original container, it should never be poured back because it can contaminate the remaining reagent.
- b. Reagents are preferably stored in alphabetical order on shelving protected from dust, moisture, and direct sunlight.
- c. Never use a reagent that cannot be clearly identified from the label on the container. Discard all reagents that cannot be accurately identified.
  - d. Always read the label before dispensing a reagent.
- e. When working with newly prepared reagents, especially stains, ascertain whether desired results are being obtained. Unsatisfactory solutions should be discarded and replaced.

MD0853 2-2

- f. All mixing containers, stirring rods, and containers used for storage of reagents should be chemically cleaned prior to use.
- g. During mixing and preparation, as well as in storage, it is good practice to avoid contact of reagents with metals. Many reagents contain substances that will react chemically with metals and produce changes that will render than unusable for laboratory work.
- h. Do not allow inexperienced personnel to prepare reagents without close supervision.
- i. Certain reagents are poisonous and adequate precautions should be taken to prevent accidental poisoning. All highly toxic reagents should be conspicuously labeled "POISON" and should be stored in a separate cabinet in the laboratory.
- j. Commercial reagents should be checked with standards for purity. Record all lot numbers in case a reagent is not pure.
  - k. Test all new reagents to assure that proper results are attainable.

#### Section II. LABORATORY GLASSWARE

#### 2-4. DILUTING PIPETS

- a. The blood cell diluting pipet (see figure 2-1) consists of a graduated capillary tube having an arbitrary volume of one unit and marked in increments of that unit, each designated as 0.1; note that this unit is not a standard measurement but merely an arbitrarily selected unit. Above the capillary tube is a mixing bulb containing a color-coded glass bead, and above the bulb another shorter capillary tube with an engraved mark. The pipet for performing the white blood cell count has a white bead, the mixing bulb is smaller than that *of* the red count pipet, and the marking above the bulb reads "11." The pipet for performing the red blood cell count has a red bead in the mixing bulb and the marking above the bulb is "101."
- b. These pipets are used to take the specimen directly from a capillary puncture or, after careful mixing, from a vial of fluid or of blood treated with an anticoagulant, such as EDTA. The blood or fluid is drawn into the pipet to a predetermined point and diluted to the correct mark with diluting fluid. After proper mixing, the diluted substance is placed in the counting chamber and the cells are counted.
- c. Usually, the technique for diluting the blood specimen with a pipet calls for whole blood to be drawn exactly to the 0.5 mark and diluted only to the "11" or "101" mark with appropriate diluting fluid dependent upon the type of cell count. Since the volume of fluid in the stem does not enter into the dilution, the dilution is calculated on the volume in the bulb, thus with the white blood cell count:




Figure 2-1. Hematological pipets.

$$\frac{\text{Whole blood}}{\text{Dilution = Volume in bulb}} = \frac{0.5}{10} = \frac{1}{20}$$
and in the RBC count

$$\frac{\text{Whole blood}}{\text{Dilution = Volume in bulb}} = \frac{0.5}{100} = \frac{1}{200}$$

The dilution of the blood sample in the white blood count is 1 to 20 and the dilution factor is 20. In the red blood count the blood sample is diluted 1 to 200 and the dilution factor is 200. The permissible error of the red cell pipet is  $\pm 5$  percent and the white cell pipet has a permissible error of  $\pm 3.5$  percent.

- d. Blood cell diluting pipets are delicate pieces of equipment and should have careful treatment. It is also important that clean, dry diluting pipets be used to prevent dilution errors and hemolysis of cells.
- e. If an aspirator is available, pipets can be easily washed by drawing cold tap water (distilled water is preferred) through the pipet. If a deposit remains, it can be dissolved by drawing dilute sodium hypochlorite (household bleach) through the pipet followed by several washings with distilled water. If pipets become plugged through neglect, the hole may be opened with a fine steel wire, and the above procedures may then be carried out. When necessary, pipets may be dried rapidly by drawing a small amount of alcohol through then and then ether or acetone followed by a stream of air until the glass bead removes freely and no moisture remains in the pipet.

f. All pipets should be inspected each time before use so as to prevent the use of equipment that is dirty or that has chipped ends.

## 2-5. UNOPETTE SYSTEM

a. The Unopette System (Becton- Dickinson and Co.) consists of a disposable uniform-bore glass capillary pipet with an attached plastic tab for handling. The pipet (see figure 2-2) is attached to a plastic reservoir in which predetermined amounts of diluting fluids, depending on their purpose, can be placed. The blood aspirated into the diluting fluid can be mixed and dispensed with great ease and speed.




Figure 2-2. Unopette System.

- b. The Unopette System can be used for a variety of hematological procedures. In general, the Unopette System is used as follows:
- (1) Puncture diaphragm. Using the protective shield on the capillary pipette, puncture the diaphragm of the reservoir as follows:
- (a) Place reservoir on a flat surface. Grasping reservoir in one hand, take pipette assembly in other hand and push tip of pipette shield firmly through diaphragm in neck of reservoir, then removed.




Figure 2-3. Unopetted filling procedure (continued)

(b) Remove shield from pipette assembly with a twist.




Figure 2-3. Unopetted filling procedure (continued)

- (2) Add sample. Fill capillary with whole blood and transfer to reservoir as follows:
- (a) Holding pipette almost horizontally, touch tip of pipette to blood. Pipette will fill by capillary action. Filling is complete and will stop automatically when blood reaches end of capillary bore in neck of pipette.




Figure 2-3. Unopetted filling procedure (continued)

- (b) Wipe excess blood from outside of capillary pipette making certain that no sample is removed from capillary bore.
- (c) Squeeze reservoir slightly to force out some air. Do not expel any liquid. Maintain pressure on reservoir.




Figure 2-3. Unopetted filling procedure (continued)

(d) Cover opening of overflow chamber of pipette with index finger and seat pipette securely in reservoir neck.




Figure 2-3. Unopetted filling procedure (continued)

- (e) Release pressure on reservoir. Then remove finger from pipette opening. Negative pressure will draw blood into diluent.
- (f) Squeeze reservoir gently two or three times to raise capillary bore, forcing diluent up into, but not out of, overflow chamber, releasing pressure each time to return mixture of reservoir.



Figure 2-3. Unopetted filling procedure (continued)

(g) Place index finger over upper opening and gently invert several times to thoroughly mix blood with diluent.




Figure 2-3. Unopetted filling procedure (concluded).

(h) Let stand for ten (10) minutes to allow red cells to hemolyze.

#### 2-6. BLOOD CELL COUNTING CHAMBERS

a. The most common type of hemacytometer consists of two counting chambers separated by grooves or canals. On the smooth glass surface of the counting chambers are straight lines etched into glass in a gridwork pattern. The Neubauer ruling, preferred for hematological work, consists of a gridwork with dimensions of 3 mm by 3 mm. It is further divided into 9 smaller squares with dimensions of 1 mm by 1 mm; 4 of these squares are used for the white count. The 8 outer squares are further subdivided into 16 squares 0.25 mm on a side. The central square is divided into 25 squares, 0.20 mm on a side, which are used for the red cell count. Thus the large squares are 1 square mm, the 16 small squares in the outer large squares are 1/16 square mm and the 25 central squares are 1/25 square mm (see figure 2-4).




Figure 2-4. Rulings on a hemacytometer.

b. The cover glass must be free of visible defects and must be optically plane on both sides within  $\pm$  0.002 mm according to the United States (US) Bureau of Standards. When the cover glass is placed on the platform, the space between it and the ruled platform should be 0.1 mm.

#### 2-7. CLEANING

Glassware can be cleaned in hot, soapy water and thoroughly rinsed in distilled water. Blood dilution pipets can be washed by flushing water and acetone through them. Ordinary household bleach can be used to remove blood clots in the bore of pipets. Dilution pipets are dry when the bead moves freely in the bulb. Glassware to be used for coagulation studies must be scrupulously clean. This glassware should be

cleaned in nonorganic detergents and rinsed well with distilled water. The etched surface of the hemacytometer should be rinsed in water and blotted dry with lens paper to avoid marring or further etching of the lines on the surfaces. A method of cleaning small bore tubes and pipets (such as Wintrobe sedimentation rate tubes) is to attach a capillary pipet, by a rubber hose, to a water type suction pump. Attach the tube to the flat end of the pipet and hold the tube under water. Blood is drawn from the tube as water is drawn in. When the tube is clean, invert the tube, remove the residual water by suction, and allow to dry. Occasionally, the tube should be cleaned with dilute sodium hypochlorite (household bleach) to remove deposits of residual blood.

# Section III. LABORATORY EQUIPMENT

#### 2-8. MICROSCOPES

a. **Introduction**. A modern microscope for use in the hematology laboratory is equipped with an illuminator system, a substage condenser system, an objective system, a projector (eyepiece or ocular system), an iris diaphragm, nicol prisms, a tubular barrel (monocular or binocular bodies), and a mechanical stage (see figure 2-5). A compound microscope uses a combination of lenses, the objective lens (lens closer to the object) and the ocular lens (lens closer to the eye) to project the image to the retina of the eye. The objective lens acts much like a small projection lens which projects an enlarged primary image near the top of the tubular barrel. This image, formed in air, is known as an "aerial image". This object is viewed through the projector or eyepiece that acts like a magnifier except that it magnifies an aerial object instead of an actual object. The final image projected on the retina of the eye is called a "virtual image" because the light rays appear to come from the image. The rays are actually created by an increase in magnification by the lens system.

# b. Magnification.

- (1) Magnification in a microscope is limited to the useful magnification that can be achieved, that is, the ability to obtain fine detail of the object being examined. This ability to render visible the fine detail is the resolving power of the microscope. The resolving power of a microscope is dependent on the numerical aperture (N.A.) of the objective lens and condenser lens. Therefore, proper adjustment of these lenses is essential in order to obtain useful magnification.
- (2) Microscopes in general use in medical laboratories are provided with three objectives with focal lengths of 1.9 mm, 4 mm, and 16 mm, respectively. Microscopes are usually provided with 5X and 10X (most common) oculars. Multiplying the power of the ocular by the power of the objective gives the degree of magnification of the object under observation. The degree of magnification is expressed in diameters (refers to an increase in diameter). The ocular magnification, the millimeter length of the objective, its magnification power, and the total apparent increase obtained using oculars and objectives of the powers shown are given:




Figure 2-5. Compound microscope.

<u>Ocular</u>	<u>Objective</u>	<u>Magnification</u>
10X	16 mm (10X)	100 diameters
10X	4 mm (43X)	430 diameters
10X	1.9 mm (97X)	970 diameters

Magnification is increased in practice by using a higher power objective. Most microscopes are equipped with a revolving nosepiece, and selection of an objective lens is done with ease.

#### c. Illumination.

- (1) Correct illumination of the object under study is an extremely important detail. Incorrect lighting of the object can lead to inaccurate results and conclusions. Correct illumination can be obtained from a concave mirror or substage light.
- (2) Illumination entering the microscope can be central or oblique. To obtain central illumination the light from the source must be reflected from the mirror directly into the microscope tube. Oblique lighting must be avoided because an object in the center of the field will sway from side to side when the fine adjustment is rotated. Oblique light is usually no problem with substage illuminators.

- (3) Regulation of the amount of light admitted is accomplished by the iris diaphragm in the substage condenser. The size of the opening in the diaphragm is controlled by a lever on the side of the condenser. The lever of the iris diaphragm should never be forced to the full limit in either direction. Doing so may damage the delicate leaves of the diaphragm. Generally, when observing liquid preparations under low power, the diaphragm opening should be partially closed. Under the high dry objective, the diaphragm is generally opened to a greater degree to allow more light to pass through the material. When observing stained preparations under the oil immersion objective, the iris diaphragm is usually opened wide.
- (4) The substage condenser functions to direct a light beam of the desired numerical aperture (N.A.) and field size onto the specimen. The size of the opening in the condenser together with its position up or down controls the light entering the system. When the condenser is close to the stage, concentration of light is greater; as the condenser is moved downward, less light passes upward through the object under observation.
- (5) Improper illumination is indicated when: (1) dark points or shadows appear in the field; (2) the outline of an object is bright on one side and dark on the other; or (3) the object appears to be in a glare of light. This can usually be corrected by changing the position of the mirror, by reducing the amount of light by adjusting the size of the opening in the iris diaphragm, or by raising or lowering the condenser.
- d. **Focusing**. Focusing can be defined as the adjustment of the relationship between the optical system and the object so that a clear image is obtained. Several important rules to be observed when focusing the microscope on the preparation are:
- (1) After the object is mounted on the stage, the objective to be used is turned into line with the eyepiece.
- (2) Movement of the objective is accomplished by revolving the nosepiece. The nosepiece is provided in order to enable rapid, convenient substitution of one objective for another. This change is effected by grasping two of the objectives between the thumb and forefinger of the right hand and rotating them until the desired objective is brought into line with the axis of the body tube. It is very important that exact alignment be obtained. The correct setting is indicated by a slight "click" as the objective comes into position.
- (3) Whenever the nosepiece is revolved, its movements should be observed to make certain that the objectives do not come into contact with the object. Some microscopes are not parfocal; that is, objects in focus under low power will not be in focus when the nosepiece is rotated to a higher power of magnification. It may, therefore, be necessary to refocus when changing to higher magnification. In microscopes that are parfocal, it is possible to swing other objectives into place without touching the coarse adjustment and with only a slight turn of the fine adjustment knob required to restore perfect focusing.

(4) To bring an object into focus, watch from the side and use the coarse adjustment to lower the objective until it is below the point at which the object would normally be expected to come into view.

**NOTE**: To avoid damage to slide or microscope, view from side for preliminary focusing. Then, using the coarse adjustment and at the same time looking through the ocular, raise the objective very slowly until the field comes into view. Further adjust to the best image, using only the fine adjustment.

- (5) In focusing upward with the fine adjustment, the object will first appear in faint outline, then gradually more distinctly, and finally, sharply defined. If the adjustment goes beyond the point of sharp definition, return to the point of greatest clarity by using the fine adjustment.
- (6) Never move an objective downward while looking through the eyepiece. When the objective is moved downward, always observe the downward motion with the eye held level with the microscope stage. Failure to observe these precautions can result in damage to the lens of the objective or the object under study.
- e. **Care of the Microscope**. The microscope is an instrument of precision with many delicate parts, and it must be handled with the utmost care. Care should not be confined to the optical elements alone. The microscope is a combination of optical and mechanical excellence, one complementing the other. The following precautions should always be observed in the care of the microscope:
  - (1) No unauthorized person should manipulate the microscope.
- (2) Keep the microscope as free from dirt and dust as possible. Dusty lenses produce foggy images, while dust in the focusing mechanisms causes excessive wear of those parts.
  - (3) The microscope should be always covered when not in use.
- (4) Care should be taken to prevent all parts of the microscope from coming into contact with acid, alkali, chloroform, alcohol, or other substances that corrode metal or dissolve the cementing substance by means of which the lenses are secured into the objectives and oculars.
  - (5) Always carry the microscope with two hands by the arm and base.
- (6) Avoid sudden jars *I* such as placing the microscope on the table with undue force.
- (7) No dust should be permitted to settle on the lenses nor should the finger come in contact with any of the surfaces.
- (8) The lens system should never be separated, as the lenses are liable to become decentered and dust can enter.

- (9) Avoid all violent contact of the objective lens and the cover glass.
- (10) Keep eyepieces in the microscope at all tines to keep free of dust.
- (11) To remove dust, brush the lenses with a camel's hairbrush. Avoid hard wiping, as dust is often hard and abrasive.
- (12) Xylene is the only agent that should be used in cleaning lenses or removing oil from objectives. Only small amounts of xylene are necessary and should be used with care.
- (13) When sewing machine oil is used to lubricate moving parts of the microscope, all excess should be wiped off to prevent the collection of grit and dust.
- (14) The microscope should be protected against direct sunlight and moisture.
- (15) In very warm, humid climates, microscopes should be stored in dry cabinets when not in use. Such cabinets should be reasonably airtight, equipped with a light bulb to supply heat, and several cloth bags containing a hygroscopic salt, such as calcium chloride, to absorb moisture. In warm, humid climates, the lenses of unprotected microscopes can be attacked by certain fungi that etch glass and ruin the lenses.
- (16) After use, always turn the nosepiece to a position, which brings the low power objective into direct line with the opening in the substage condenser. If this precaution is not taken, the longer, higher-powered objectives can accidentally come into contact with the condenser lens.
- (17) The entire microscope should be cleaned frequently to remove dust, finger marks, oil, grease, and remnants of specimens. All parts of the microscope should be kept scrupulously clean at all times.
- (18) Never tamper with any of the parts of the microscope. If the instrument does not seem to be functioning properly, immediately call the matter to the attention of the laboratory supervisor.
- (19) Maintenance of the microscope should be done in accordance with the manufacturer's booklet of instruction.
- (20) Immediately after use, the oil immersion objective must be wiped clean of oil with a soft, absorbent lens paper.

#### f. Types.

(1) <u>Binocular</u>. This type of microscope is preferred since it gives a more natural and restful condition of observation. A beam-dividing prism and three mirrors

divide the light equally, sending half to the left eye and half to the right eye. The coating on the mirrors is enhanced by aluminum to increase the reflectivity. The binocular body is protected by cover glass seals to keep dust from entering. The binocular compound microscope is the preferred microscope for routine hematology.

(2) <u>Phase</u>. Phase microscopy is becoming increasingly prevalent in platelet counting. In bright-field illumination, a completely transparent specimen is difficult to see in any detail. By using phase contrast, transparent living objects can be studied. Phase microscopy operates on the principle that if a portion of light is treated differently from the rest, and caused to interfere with the rest, it produces a visible image of an otherwise invisible transparent specimen. Phase contrast accessories are available from the standard optical companies.

#### 2-9. CENTRIFUGES

These are laboratory devices or units that apply a relatively high centrifugal force (up to 25,000 g) to a specimen, causing its separation into different fractions according to their specific gravities. Centrifugation is the process of separating components of a mixture (away from a center as in centrifugal force) on the basis of differences in densities of the different components using a centrifuge.

- a. **Table Top Models**. These units are mounted on rubber feet that absorb vibration. The speed is controlled by means of a rheostat on the front panel. Top speeds of centrifuges will vary and the top speed of a particular instrument should be known in order to use the speed control device. Those centrifuges have adapters to hold 6 tubes and adapters for 12 tubes.
- b. **Floor-Mounted Models**. The heavier floor-mounted models accommodate a large number of tubes at one time. The top speed of these instruments is higher than that of table models. Because of their increased inertia, they are equipped with a brake to facilitate stopping. In these units, the tubes are placed in balanced receptacles that are mounted on spokes emanating from a central hub.
- c. **Microhematocrit Centrifuge**. This centrifuge is a special type of high-speed centrifuge employed to spin capillary tubes (see figure 2-6). The circular tube holder on this centrifuge is flat, surrounded by a rubber ring. It has a capacity of 24 capillary tubes. After a capillary tube is filled with blood, it is closed with a commercial plastic sealing material. During centrifugation the sealed end is always placed in position facing toward the outside of the holder plate. Most centrifuges of this type spin the tubes at 10,000 rpm.
- d. **Precautions**. In all instances where centrifugation is required, careful attention must be given to balancing the units. This means that tubes must be placed exactly opposite each other, they must be of identical weight, and they must contain the same amount of fluid. If at all possible, centrifuges should be equipped with tachometers so that speed nay be checked and controlled. Certain procedures, such as hematocrits, require a critical relative centrifugal force (RCF or g). Relative centrifugal

force is the weight of a particle in a centrifuge relative to its normal weight, the centrifugal force per unit mass in gravities (g). To determine the RCF (or g) for these procedures, consult the serology manual or a monograph. The inside of the centrifuges should occasionally be cleaned to prevent dust particles from being blown into specimens. The lid on the centrifuge should be closed and locked before and during operation. Only open the lid when the centrifuge has stopped rotating.



Figure 2-6. Microhematocrit centrifuge and microhematocrit reader.

**Continue with Exercises** 

## **EXERCISES, LESSON 2**

**INSTRUCTIONS**: Answer the following exercises by marking the lettered response that best answers the exercise, by completing the incomplete statement, or by writing the answer in the space provided at the end of the exercise.

After you have completed all of the exercises, turn to "Solutions to Exercises" at the end of the lesson and check your answers. For each exercise answered incorrectly, reread the material referenced with the solution.

- 1. When using laboratory reagents, what routine hematological care should be taken in their preparation?
  - a. Follow detailed procedural directions.
  - b. Measure regents exactly.
  - c. Follow the order that reagents are to be added.
  - d. Control the temperature where indicated.
  - e. Use filtration of ingredients as stated.
  - f. Age the prepared solution as specified.
  - a. All of the above.
- 2. Particular attention must be given to storage of laboratory reagents, particularly with reference to requirements for:
  - a. Refrigeration.
  - b. Incubation.
  - c. Protection from intense light.
  - d. All of the above.

- 3. A reagent container is correctly labeled if the label contains the:
  - a. Constituents, initials of the individual who prepared the reagent in alphabetical order, expiration date, quantity used, and date prepared.
  - b. Expiration date, initials of the individual who prepared the reagent, constituents quantity used, temperature control, and date prepared.
  - c. Initials of the individual who prepared the reagent, stains, constituents, expiration date, quantity used, and date prepared.
  - d. Constituents, expiration date, initials of the individual who prepared the reagent, quantity used, and date prepared.
- 4. A reagent's container label is properly labeled and protected if the label is:
  - a. Complete and securely attached.
  - b. Neatly and legibly written or preferably typewritten.
  - c. Covered with a protective coating of cellophane tape over the surface of the label.
  - d. All of the above.
- 5. Why must an unused portion of a reagent never be poured back into tile original container?
  - a. Possible explosions.
  - b. Inefficiency.
  - c. Contamination.
  - d. Improper labeling.

	a.	Dust.
	b.	Moisture.
	C.	Direct sunlight.
	d.	All of the above.
7.		at should you do if you cannot clearly read the reagent's label and/or identify contents of the container?
	a.	Use the contents.
	b.	Check with your colleague.
	C.	Discard it properly and use a reagent's container with a label that you can read and the contents you can identify.
	d.	Save it for next time.
8.		en preparing laboratory reagents for storage, which items should be chemically aned prior to use?
	a.	Mixing containers.
	b.	Stirring rods.
	C.	Storage containers.
	d.	All of the above.

6. When storing reagents on shelving, protect them from:

9.	At what point is it is a good, safe, practice to avoid contact of reagents with since metals may become unusable for laboratory work?	metals
	a. Labeling and heating.	
	b. Refrigeration and labeling.	
	c. Preparation and mixing.	
	d. Preparation and heating.	
10.	Highly toxic reagents should be conspicuously labeled:	
	a. Reagent.	
	b. "POISON."	
	c. Poison.	
	d. "REAGENT".	
11.	The white blood cell diluting pipet has a:	
	a. Smaller bulb, a red bead, and an "11."	
	b. Smaller bulb, a white bead, and an "11."	
	c. Larger bulb, a red bead, and a "101."	
	d. Larger bulb, a white bead, and a "101."	
12.	The red blood cell diluting pipet has a:	
	a. Smaller bulb, a red bead, and an "11."	
	b. Smaller bulb, a white bead, and an "11."	

d. Larger bulb, a white read, and a "101."

c. Larger bulb, a red read, and a "101."

- 13. The 0.5 mark on the RBC diluting pipet designates a volume equal to:
  - a. 0.5 ml.
  - b. 0.5 cu mm.
  - c. 1/20 of the dilution volume.
  - d. 1/200 of the dilution volume.
- 14. Which statement is correct concerning the use of the diluting pipets for cell counting?
  - a. Blood or fluid is drawn into the pipet to a predetermined point and diluted to the correct mark with diluting fluid. After proper mixing, the diluted substance is placed in the counting chamber and the cells are counted.
  - b. Blood or fluid is drawn into the pipet to an arbitrary point and diluted to the correct mark with diluting fluid. After proper mixing, the diluted substance is placed in the counting chamber and the cells are counted.
  - c. Blood or fluid is drawn into the pipet to a predetermined point and diluted to a guesstimated mark with diluting fluid. After gentle mixing, the diluted substance is placed in the counting chamber and the cells are counted.
- 15. Usually, the technique for diluting the blood specimen with a pipet calls for whole blood to be diluted only to the "11" or "101" mark with appropriate diluting fluid dependent upon the type of cell count. How is the dilution calculated?
  - a. The volume of fluid is measured in the stem and then calculated.
  - b. Dilution is calculated on the volume in the bulb.
  - c. Both a and b.
  - d. None of the above

- 16. In general, what is the proper procedure for puncturing the diaphragm using the Unopette System?
  - a. Using the protective shield on the capillary pipette, puncture the diaphragm of the reservoir.
  - b. Grasping the reservoir in one hand, take pipette assembly in other hand and pull tip of pipette shield firmly through diaphragm in neck of reservoir, then remove.
  - c. Grasping the reservoir in one hand, take pipette assembly in other hand and push tip of pipette shield firmly through diaphragm in neck of reservoir, then remove.
  - d. Remove shield from pipette assembly with a twist.
  - e. a, b, c, and d.
  - f. a, b, and d.
  - g. a, c, and d.
- 17. What phenomenon is used to fill the Unopette capillary with blood?
  - a. Gravity.
  - b. Capillary action.
  - c. Brownian movement.
  - d. Static electricity.
- 18. When using the Unopette system, at what point will the pipette automatically stop filling and be complete? When:
  - a. The blood reaches the top reservoir.
  - b. Negative pressure is applied,
  - c. The blood reaches the end of the capillary bore.

19.	the	w many times should the reservoir be squeezed and with what pressure to raise capillary bore and force the diluent up into, but not out of, the overflow amber using the Unopette system?
	a.	4 to 6; hard.
	b.	3 to 5; even.
	C.	2 to 3; gentle.
	d.	1 to 6; steady.
20.		ing the Unopette System, after the blood is thoroughly mixed with the diluent d left to stand for 10 minutes, the red cells will:
	a.	Overflow.
	b.	Settle out.
	C.	Hemolyze.
21.	Wr	nat are the outer dimensions of the Neubauer ruling?
	a.	0.20 by 0.20 mm.
	b.	0.25 by 0.25 mm.
	C.	1 by 1 mm.
	d.	3 by 3 mm.
22.		nat configuration does the most common type of hemacytometer look like for unting blood cells?
	a.	Two counting chambers separated by grooves or canals.
	b.	Three counting chambers separated by grooves or canals.
	C.	Two counting chambers separated by rough indentations.

d. Four counting chambers that are not separated.

23.	How many squares are used to count white cells, when the dimensions of the Neubauer ruling are further divided into 9 smaller squares, with dimensions of 1 mm by 1 mm?		
	a.	2.	
	b.	3.	
	C.	4.	
	d.	5.	
24.		ng the Neubauer ruling for the red cell count, which portion and how many as are used?	
	a.	Middle squares; 0 30 mm on a side.	
	b.	Outer 8 squares; 0.15 mm on a side.	
	C.	10 outer squares; 0 05mm on a side.	
	d.	Central square; 0.20 mm on a side.	
25.	Wh	en taking a blood cell count, the cover glass must be free of visible and optically on both sides.	
	a.	Defects; plane.	
	b.	Outer squares; plane.	
	C.	Stains; rough.	
	d.	Blood; clean.	

- 26. A method for cleaning Wintrobe sedimentation rate tubes, small bore tubes and pipets, and remaining blood clots in the bore of a pipet is to:
  - a. Attach a capillary pipet by a rubber hose, to a water type suction pump.
  - b. Attach the tube to the flat end of the pipet and hold the tube under water.
  - Draw the blood from the tube as water is drawn in.
  - d. Ensure the tube is clean by inverting it, remove the residual water by suction, and allow to dry.
  - e. All of the above.
  - f. None of the above.
- 27. What occasional method it used to clean small bore tubes and pipets?
  - Use dilute sodium hypochlorite (household bleach) to remove deposits of residual blood.
  - b. Use distilled water.
  - c. Insert a small rubber hose.
  - d. Use a suction tube.
  - e. Absorbent lens paper and water.
- 28. Blood dilution pipets can be washed by:
  - a. Soaking in calcium chloride.
  - b. Soaking in ammonia.
  - Flushing water and acetone through them.
  - d. Rinsing in distilled water.

29.	Glassware, which is used for coagulation studies, must becleaned in		
	a.	Gently; household bleach.	
	b.	Freely; distilled water.	
	C.	Scrupulously; non-organic detergent.	
	d.	Somewhat organic soapy water.	
	e.	Vigorously; in water with absorbent lens paper.	
30.		lect the correct items normally found on a modern microscope used in a matology laboratory.	
	a.	A darkening system, a substage condenser system, an objective system, a projector (eyepiece or ocular system), an iris diaphragm, nicol prisms, a tubular barrel (monocular or binocular bodies), and a mechanical stage.	
	b.	An illuminator system, a substage condenser system, an objective system, a projector (eyepiece or ocular system), an iris diaphragm, nicol prisms, a tubular barrel (monocular or binocular bodies), and a mechanical stage.	
	C.	An illuminator system, a substage condenser system, an objective system, a projector (eyepiece or ocular system), a round barrel (monocular or binocular bodies), and a mechanical stage.	
	d.	An illuminator system, a substage condenser system, an objective system, a projector (eyepiece or ocular system), an iris diaphragm, nicol prisms, a tubular barrel (monocular or binocular bodies), and a mechanical stage.	
31.	Wh	nat combination of lenses does a compound microscope use?	
	a.	Objective lens.	
	b.	Monocular lens.	
	C.	Aerial image magnifier.	
	d.	a and b.	
	e.	b and c.	

32.	32. Select the best explanation of an "aerial image".	
	a.	An image formed in the air.
	b.	An image formed in the air. The object is viewed through the projector or eyepiece that acts like a magnifier except that it magnifies an aerial object instead of an actual object.
	C.	An image formed in the air. The object is viewed through the magnifying glass except that it magnifies an aerial object instead of an actual object.
	d.	An image formed on a surface. The object is viewed through the projector or eyepiece that acts like a magnifier except that it magnifies an aerial object instead of an actual object.
33.		e ability of a microscope to render fine detail is dependent upon the numerical erture and proper adjustment of which lens (es)?
	a.	Ocular and objective.
	b.	Ocular and condenser.
	C.	Objective and condenser.
	d.	Objective only.
34.		nen rotation of a microscope's fine adjustment causes an object in the center of field to sway from side to side, the lighting is:
	a.	Central.
	b.	Oblique.
	C.	Too dim.
	d.	Too intense.

35.	Name one item the resolving power of the microscope is dependent upon during magnification?	
	a.	Focal lengths.
	b.	Binocular bodies.
	C.	Arial image.
	d.	N.A. of the objective.
36.	Pre	eliminary focusing of a microscope should be observed from the:
	a.	Ocular.
	b.	Objective.
	C.	Top of the microscope.
	d.	Side of the microscope.
37.		ce correct illumination of an object under study is an extremely important detail, at can incorrect lighting cause?
	a.	Inaccurate results and conclusions.
	b.	Inaccurate steps and timings.
	c.	Faulty conclusions and recommendations.
	d.	Changing positions and recommendations.
38.	What is the function of the substage condenser when illuminating slides under the microscopic?	
	a.	Indicates dark spots.
	b.	Reduces glare.
	C.	Directs a light beam.
	d.	Correct inconsistencies by changing the position left or right.

- 39. Which statement about the substage condenser is true?
  - a. The substage condenser functions to direct a hair beam of the desired numerical aperture (N.A.) and field size onto the specimen.
  - b. The size of the opening in the condenser together with its position up or down controls the light entering the system.
  - c. When the condenser is open to the stage, concentration of light is greater.
  - d. As the condenser is moved upward, less light passes downward through the object under observation.
- 40. Improper illumination is indicated when:
  - a. Light points appear on the outer edges of the slide.
  - b. The center of an object is bright on one side and dark on the other.
  - c. The object appears to be in dull light
  - d. Shadows appear in the field.
- 41. Which three parts of a microscope may be adjusted to control the illumination?
  - a. Mirror, iris diaphragm, and ocular.
  - b. Mirror, iris diaphragm, and objective.
  - c. Mirror, iris diaphragm, and condenser.
  - d. Ocular, objective, and condenser.

42.	When using a microscope the "virtual image" projected on the retina of the e is the image.	
	a.	Initial.
	b.	Intermediate.
	C.	Final.
	d.	Aperture.
43.	Wh	nich of the following is used to clean microscope lenses?
	a.	Xylene.
	b.	Acetone.
	C.	Household bleach.
	d.	Saturated sodium hydroxide.
44.	Wh	nich type of microscope is preferred for routine hematology?
	a.	Monocular.
	b.	Binocular.
	C.	Fluorescence.
	d.	Phase.
45.	Ве	cause of dust, the lens system should:
	a.	Be separated.
	b.	Never be separated.

	a.	Plastic sealing material.
	b.	Paper.
	C.	Glass.
	d.	Wax.
47.	Мо	st microhematocrit centrifuges have a speed of about:
	a.	500 rpm.
	b.	1,000 rpm.
	C.	5, 000 rpm.
	d.	10,000 rpm.
48.	Wh	ich centrifuge has the higher type of top speed for instruments?
	a.	Table top model.
	b.	Floor-mounted model.
49.	Onl	y which agent should be used to remove oil from a lens object?
	a.	An abrasive.
	a. b.	
	b.	A rag.

46. After a capillary centrifuge tube is filled with blood, it is sealed with:

- 50. Whenever centrifugation is required, which precaution must always be followed?
  - a. Careful attention must be given to balancing the units. This means that tubes must be placed exactly opposite each other, they must be of identical weight, and they must contain the same amount of fluid.
  - b. Centrifuges should be equipped with a microhematocrit so that speed may be checked and controlled.
  - c. Hematocrits require a regular force and can adapt for 6 to 19 tubes.
  - d. The heavier, floor-mounted models require rubber feet to absorb vibrations so as to accommodate a large number of tubes, which are housed in mounted receptacles on spokes from a distance hub.

**Check Your Answers on Next Page** 

# **SOLUTIONS TO EXERCISES, LESSON 2**

- 1. g (para 2-1)
- 2. d (para 2-1)
- 3. d (para 2-2)
- 4. d (para 2-2)
- 5. c (para 2-3a)
- 6. d (para 2-3b)
- 7. c (para 2-3c)
- 8. d (para 2-3f)
- 9. c (para 2-3g)
- 10. b (para 2-3i)
- 11. b (para 2-4a)
- 12. c (para 2-4a)
- 13. d (para 2-4c)
- 14. a (para 2-4b)
- 15. b (para 2-4c)
- 16. g (para 2-5b(1))
- 17. b (para 2-5b(2))
- 18. c (para 2-5b(2) (a))
- 19. c (para 2-5b(2)(f))
- 20. c (para 2-5b(2) (h))
- 21. d (para 2-6a)
- 22. a (para 2-6a)

- 23. c (para 2-6a)
- 24. d (para 2-6a)
- 25. a (para 2-6b)
- 26. e (para 2-7)
- 27. a (para 2-7)
- 28. c (para 2-7)
- 29. c (para 2-7)
- 30. d (para 2-8a)
- 31. d (para 2-8a)
- 32. b (para 2-8a)
- 33. c (para 2-8b(1))
- 34. b (para 2-8c(2))
- 35. d (para 2-8b(1))
- 36. d (para 2-8d(4) NOTE)
- 37. a (para 2-8c(1))
- 38. c (para 2-8ac(4))
- 39. b (para 2-8c(4))
- 40. d (para 2-8c(5))
- 41. c (para 2-8c(5))
- 42. c (para 2-8a)
- 43. a (para 2-8e(12))
- 44. b (para 2-8f(1))
- 45. b (para 2-8e(8))

- 46. a (para 2-9c)
- 47. d (para 2-9c)
- 48. b (para 2-9b)
- 49. c (para 2-8e(12))
- 50. a (para 2-9d)

**End of Lesson 2** 

#### **LESSON ASSIGNMENT**

LESSON 3

Collection of Blood and Preparation of Blood Smears.

**TEXT ASSIGNMENT** 

Paragraphs 3-1 through 3-8.

**LESSON OBJECTIVES** 

After completing this lesson, you should be able to:

- 3-1. Select the statement which best describes the requirements for selection and care of blood collection equipment.
- 3-2. Select the correct steps in collecting, processing, and recording blood specimens.
- 3-3. Select the names and functions of commonly used anticoagulants
- 3-4. Select the correct method and use the proper procedures in staining blood films for the type of blood cells found.
- 3-5. Select the correct function for each of Wright's component stain and buffer solutions.
- 3-6. Select the correct steps to prepare a blood smears.
- 3-7. Select the factors that affect the quality of stained blood smears.

SUGGESTION

After completing the assignment, complete the exercises at the end of this lesson. These exercises will help you to achieve the lesson objectives.

#### LESSON 3

## Section I. COLLECTION OF BLOOD SPECIMENS

## 3-1. INTRODUCTION

- a. Hematological laboratory procedures are based upon the examination of blood specimens. To obtain valid test results, specimens must be properly collected, processed, and recorded. Blood specimens are usually obtained by either venous or capillary puncture. The source of the specimen is determined chiefly by the quantity of blood required to perform the laboratory procedures and the age and condition of the patient.
- b. There is generally little difference in blood counts performed on venous or capillary blood if a free-flowing capillary blood specimen is obtained. Valid blood counts cannot be made when capillary specimens are not taken from a free-flowing sample or when they are obtained from cyanotic or calloused areas or areas of local stasis. White blood cell counts made on blood obtained from such sources can vary as much as 1000-1500 cells per cu mm from their real value. For general purposes, however, venous samples are preferable since they allow for multiple and repeated hematological examinations and provide a sufficient quantity of blood for performing any other required laboratory procedure. Further, with venous blood the chances of error are reduced because operations are made under ideal conditions and repeat operations are possible. In situations where there are limitations on the quantity of blood that can be obtained, that is, in small infants or extensive burn cases, microquantitative methods are satisfactory for performing an analysis on a specimen obtained by capillary puncture.

#### 3-2. VENIPUNCTURE

a. **Site**. To obtain blood by venipuncture, draw the specimen directly from a patient's vein with a sterile hypodermic needle and syringe or a vacuum blood sample device. In adults use the veins located in the proximal forearm or antecubital space as illustrated in figure 3-1. In infants employ the jugular or femoral vein for the venipuncture. The vein selected should be large, readily accessible, and sufficiently close to the surface to be seen and palpated. If venipuncture poses a problem due to the age of the patient, sclerotization due to repeated venipuncture, or any other unusual circumstance, the technician should consult a physician concerning the procedure. UNDER NO CIRCUMSTANCES SHOULD A TECHNICIAN WITHDRAW BLOOD FROM A SAGITTAL SINUS, JUGULAR VEIN, OR FEMORAL VEIN. This should be left to the discretion of the physician in charge. Occasionally, the best vein is found on the hand, leg, or foot. These areas are more sensitive, and the veins are not as firmly anchored as those of the arm. Veins can become distended and easier to enter by allowing the arm to hang down for 2 or 3 minutes, by massaging the blood vessel




Figure 3-1. Site of venipuncture.

toward the body, or by gently slapping the site of puncture. Young and vigorous persons usually have elastic veins well filled with blood. Elderly or debilitated persons can have sclerosed or fragile veins, which are hard to enter or which collapse easily.

- b. **Equipment**. All syringes, needles, lancets, and other instruments used for the collection of blood specimens must be sterile. Disposable syringes or blood collection sets with vacuum tubes are available through normal supply channels. These should be used whenever possible. Aseptic technique is necessary to prevent the possible transmission of homologous serum hepatitis. The following equipment is necessary to perform a venipuncture:
  - (1) Isopropyl alcohol, 70 percent, prep pads.
  - (2) Tourniquet.
  - (3) Sterile syringes or vacuum blood sample devices.
  - (4) Gauze pad, 2 x 2 inches.
  - (5) Needle, 1 to 1  $\frac{1}{2}$  inches long, 19–23 gauge.
  - (6) Suitable blood collection tubes and labels.
  - (7) Gloves, latex.

# c. Preparation.

- (1) Cleanse hands thoroughly with soap and water.
- (2) Place an identifying label on the blood collecting tube.
- (3) Assemble the sterile needle and syringe. If a vacuum system is used, screw the needle into the plastic holder. Always leave the cap over the needle when not in use.
- (4) Check to be sure that the syringe works smoothly. The syringe must be dry to avoid hemolysis of the red cells. The plunger must match the syringe and must be pushed firmly to the bottom of the cylinder to prevent injection of air into the vein. This can be fatal.

## d. Syringe Procedure.

(1) Place a tourniquet around the patient's arm above the elbow tightly enough to check venous circulation, but not so tightly as to stop arterial flow. (If latex tubing is used, place it approximately 2 inches above he proposed venipuncture site). Form a loop with the longer end and draw the loop under the shorter end so that the tails of the tubing are turned away from the proposed site (see figure 3-2a).

**CAUTION**: Do not allow the tourniquet to remain in place for more than 2 minutes. Check the pulse at the wrist to make sure that arterial circulation is not cut off.




Figure 3-2a. Venipuncture procedure: Locate the vein.

- (2) Instruct the patient to make a tight fist.
- (3) By inspection and palpation locate the desired vein, determine the direction of its course, and estimate its size and depth (see figure 3-2a venipuncture procedure, a through h).

- (4) Release tourniquet.
- (5) Cleanse the skin over the selected vein with prep pads in 70 percent isopropyl alcohol. Wipe off excess alcohol with a sterile dry gauze. Do not contaminate the area after cleaning (see figure 3-2b).





Figure 3-2b. Venipuncture procedure: Clean the puncture site.

- (6) Put on gloves.
- (7) Replace tourniquet on arm and have the patient straighten out the arm and make a fist.
- (8) Grasp the syringe in the right hand and place forefinger on the hub of the needle to guide it. Grasp the forearm with the left hand about 2 inches below the area to be punctured and hold the skin taut with the thumb (see figure 3-2c).



3-2c. Venipuncture procedure: Guide needle toward the vein.

(9) With the needle bevel up, parallel to, and alongside the vein, insert the needle quickly under the skin and then into the vein. The insertion into the skin and vein can be performed in one complete motion (see figure 3-2d). After entry into the vein, blood will appear in the needle hub. Do not probe or move the needle horizontally, as discomfort and possible nerve damage may result.




Figure 2-3d. Venipuncture procedure: Insert needle into the vein.

(10) Aspiration of the blood is accomplished by gently pulling upon the syringe plunger (see figure 3-2e). The syringe barrel should be held steady during this process. Withdraw the desired quantity.




Figure 3-2e. Venipuncture procedure: Aspirate the blood.

(11) Remove the tourniquet by pulling on the long, looped end of the tubing only after blood is drawn into the syringe (see figure 3-2f).




Figure 3-2f. Venipuncture procedure: Remove the tourniquet.

CAUTION: Do not remove the needle now. If the needle were remove prior to the Tourniquet being removed, blood would be forced out of the venipuncture site, resulting in blood loss and/or hematoma formation (tumor-like cluster of blood forming under the skin).

(12) Place a sterile gauze pad over the point where the needle entered the skin and deftly withdraw the needle simultaneously putting pressure on the site (see figure 3-2g).




Figure 3-2g. Venipuncture procedure: Place a sterile pad over the site and withdraw the needle.

(13) Have the patient extend the arm and maintain light pressure on the gauze pad over venipuncture site (see figure 3-2h).




Figure 3-2h. Venipuncture procedure: Have the patient extend the arm and maintain light pressure on the site.

#### e. Vacutainer Procedure.

(1) Place the Vacutainer tube in the holder until the rubber stopper reaches the guideline. The short needle should be embedded in the stopper, but the needle must not break the vacuum (see figure 3-3).




Figure 3-3. Vacutainer system.

- (2) Follow steps 1-6 in paragraph 3-2d.
- (3) Enter the vein with the needle parallel to and alongside the vein. Probing or horizontal movement of the needle while under the skin must be avoided.
- (4) After entry into the vein push the tube all the way into the holder; vacuum is broken, and blood flows freely into the tube. Release the tourniquet at this time by pulling the long, looped end of the tube.
- (5) If the multiple needle is used or more than one tube is required, release the tourniquet after the first tube is filled; remove the filled tube and insert the next one.

**CAUTION**: Ensure the needle is not moved while tubes are being changed.

- (6) Place a sterile gauze pad over the point where the needle enters the skin and deftly withdraw the needle, placing pressure on the site.
- (7) Have the patient extend the arm and maintain light pressure on the gauze pad over the venipuncture site.

#### f. Discussion.

- (1) Cleanliness is essential when performing a venipuncture.
- (2) It is most important that correct technique be practiced in order to avoid unnecessary pain to the patient, prevent tissue damage, secure a good representative blood specimen, and prevent contamination of the specimen or infection of the patient.
- (3) Syringes and needles must be thoroughly inspected for damage or malfunction

- (4) If difficulty is experienced in entering the vein or a hematoma begins to form, release the tourniquet and promptly withdraw the needle and apply pressure to the wound
- (5) Vigorous pulling on the plunger of the syringe can collapse the vein, produce hemolysis of the blood specimen, or cause air to enter the syringe.
- (6) When repeated venipunctures have to be performed on one patient, it is advisable to select different sites for blood withdrawal.
- (7) Remove the tourniquet as early as possible once a good flow of blood has been established. Prolonged application of the tourniquet results in partial stasis of blood and changes many quantitative values of blood components.
- (8) Blood drawn by venipuncture is often stored for a period of time before it is analyzed. For this reason, certain general precautions must be followed in order to ensure a valid analysis. Before withdrawing blood from its container, make sure the blood sample is thoroughly but gently mixed. Blood containers should Be tightly stoppered at all times to prevent drying or contamination. Store the blood specimen in the refrigerator. Blood counts should be done within 3 hours after collection. Under no circumstances should blood taken for hematological examinations be stored overnight.

#### 3-3. CAPILLARY PUNCTURE

a. **Site**. Several different sites are suitable for capillary puncture. Because it is the most accessible, the palmer or lateral surface of the tip of the finger (preferably ring finger) is the most common site in adults. However, certain problems can be encountered such as heavy calloused areas or excessive tissue fluids (edema) that tend to result in non-representative samples. The lobe of the ear can be used for capillary puncture. However, differences in cell concentration do occur when blood is obtained from this site, primarily because of higher lymphocyte concentrations in the ear lobe. Because of the small amount of tissue on the fingers of infants, preferred site is the heel or big toe. A modification of the normal technique that has proven quite satisfactory when working with the heel of infants is to make two incisions in a crisscross fashion or "T".

**NOTE**: To be a valid report, work done on capillary blood must be from a FREE-FLOWING puncture wound.

# b. Equipment.

- (1) Gauze pads 2 x 2 inches.
- (2) Blood lancet.
- (3) Glass slides, heparinized capillary tubes, and other devices to receive the specimen.

(4) Isopropyl alcohol, 70 percent, prep pads.

#### c. Procedure.

- (1) The puncture site should be warm to assure good circulation of blood. If it is cold, apply warm water (38°-40°C) for a few minutes. If blood is to be drawn from the ear, the edge of the lobe, not the flat side, should be punctured.
- (2) The site to be punctured is first rubbed with alcohol prep pads to remove dirt and epithelial debris, increase circulation, and render the area reasonably disinfected (see figure 3-4 capillary puncture procedure, a through d).




Figure 3-4a. Capillary puncture procedure: Clean the puncture site.

- (3) Allow sufficient time for the circulation to equalize.
- (4) While making a finger puncture, apply gentle pressure to the finger to hold the skin taut. Hold the finger in one hand and the lancet in the other. The puncture is made perpendicular to the lines of the fingerprints, which results in a more free-flowing wound (see figure 3-4b).




Figure 3-4b. Capillary puncture procedure: Puncture the finger.

(5) The first drop of blood that appears is wiped away before specimens are taken (see figure 3-4c).




Figure 3-4c. Capillary puncture procedure: Wipe away first drop of blood.

- (6) The blood must not be squeezed out since this dilutes it with fluid from the tissues, thus altering the ratio of cellular elements to fluid, as well as the ratio of cellular elements to each other.
- (7) After the desired specimens have been collected, have the patient hold a sterile dry gauze pad over the wound until bleeding stops (see figure 3-4d).




Figure 3-4d. Capillary puncture procedure: Apply pressure to the site.

## d. Discussion.

- (1) A disposable lancet is a very satisfactory instrument for puncture of the skin. However, if this is not available a sharp pointed surgical blade is quite suitable.
- (2) Do not use the finger on a hand, which has been hanging over the side of the bed as it is likely to be congested. Edematous or cyanotic areas should not be used.

- (3) The finger should be thoroughly dry prior to puncture; blood will not well up on a finger that is moist. Furthermore, the alcohol or other antiseptic used can coagulate the blood proteins causing cell clumping and erroneous values as well as dilute cell volumes. This will result in incorrect counts and differentials.
- (4) Finger punctures should be made along the lateral aspect of the fingertip. More nerve endings are located on the fingerprint area of the fingers; therefore, more pain results from punctures in this area. Scars can also form in these sensitive areas, and difficulty may be encountered in puncturing a callous. All of these difficulties are eliminated by drawing the blood from the lateral rather than the ventral aspect of the finger.

#### 3-4. ANTICOAGULANTS

- a. Anticoagulants are used to prevent the clotting of the blood specimens and the reagent employed should not bring about alteration of blood components. Unfortunately, many anticoagulants can alter cell structures as well as coagulation. The anticoagulants most often used are ethylene-diamine-tetra-acetate (EDTA), ammonium-potassium oxalate (Heller and Paul double oxalate), and heparin.
- b. The choice of anticoagulant will depend on the analysis to be made. Ethylene-diamine-tetra-acetate (EDTA) is the anticoagulant of choice for most hematological analyses. This anticoagulant causes a minimum of distortion to the cells and platelets. It does not dissolve quickly in blood, however, so the tube must be inverted four or five times after blood is added. The dipotassium salt is prepared as a 1 percent solution in distilled water, and a final concentration of 0.5 ml of anticoagulant for each 5 ml blood is used. Another common anticoagulant is arnmonium-potassium oxalate (Heller and Paul double oxalate). This combination of oxalates does not shrink or enlarge the red blood cells appreciably. It is essential, however, to add an optimal volume of blood to the oxalate, no less than 3.5 nor more than 6.0 ml.
- c. Heparin does not alter the size of cellular components. It is, in fact, the standard for comparison of anticoagulant distortion. Heparin is more expensive and dissolves less readily than double oxalate salts. Approximately 0.5 to 1.0 mg is required to anticoagulate 5 ml of blood for 72 hours. The quantity of anticoagulant noted above in each case is sufficient to prevent clotting of the blood specimen. On the other hand, an excess of anticoagulant should be avoided because too much will result in distortion of cells and hemolysis. Ideally, differential blood smears should not be prepared from blood that contains an anticoagulant.
- d. If oxalate is added to vials and dried in an oven, take great care to avoid temperatures above 80°C. Oxalates are converted to carbonates by prolonged exposure to elevated temperatures. Under normal circumstances, it should not be necessary to prepare your own oxalate solutions since prepared anticoagulant vacuum tubes are available from Federal medical supply sources.

- e. Sodium citrate is the anticoagulant of choice for coagulation studies. It is used in a concentration of 1 part 0.11 M sodium citrate to 9 parts whole blood. It prevents coagulation by binding the calcium of the blood in a soluble complex.
- f. Sodium oxalate is another anticoagulant widely used in coagulation studies. It is used in a concentration of 1 part 0.1 M sodium oxalate to 9 parts whole blood. The sodium oxalate combines with calcium in the blood to form insoluble calcium oxalate, thereby, preventing coagulation.

**NOTE**: Under normal circumstances, it should not be necessary to prepare an anticoagulant since prepared anticoagulant tubes are available through the Federal supply system.

g. A correctly anticoagulated blood sample is essential to the proper performance of a blood cell count. The cellular constituents must remain free in the plasma and should be as similar as is possible to those remaining in the patient's circulation.

### Section II. PREPARATION AND STAINING OF BLOOD SMEARS

#### 3-5. INTRODUCTION

- a. The type of blood cells found in the peripheral blood smears may be of diagnostic and prognostic importance. For this reason proper preparation and staining of blood films is essential for the identification and study of different kinds of leukocytes. The appearance of erythrocytes and thrombocytes will often give important clues that help distinguish between different types of diseases or other physical changes.
- b. There are two basic methods for the preparation of blood smears: the cover slip and the slide methods. The cover slip method has certain advantages over the slide method; distribution of cells is like that of the in vivo circulation. The principal disadvantage of the latter method is that covers lips are very fragile and easily broken during processing.
- c. The slides and cover glasses must be chemically clean and dry. New slides must first be washed with soapy water and rinsed thoroughly with distilled water. The slides are then placed in a beaker of 95 percent ethyl alcohol. As the slides are needed, dry them with a soft, lint free cloth. The slides may be reused by properly cleaning them and making sure they are not chipped or scratched.
- d. The foundation for the morphological study of blood was based on Ehrlich's investigations of the aniline dyes, dating back to 1877, while he was still a student. Originally, simple dyes were used in the clinical laboratory and tissues were stained successively if more than one color was desired. The majority of the aniline dyes are in the form of salts of acids and bases. During the process of staining, compounds are probably formed between the basic dyes and the acid nuclear substances of cells and

between the acid dyes (so called "neutral" dyes). In this way, the staining principles of the original components were preserved; and, in addition, new staining properties dependent upon the union of the component dyes were developed. These were, therefore, termed polychromic dyes.

- e. One modification of these polychromic stains is Wright's stain. This is the stain most used in Army laboratories today. Wright's stain is a methyl alcohol solution of an acid dye and a basic dye. The acid dye is known as eosin, which is red in color. The basic dye, methylene blue, is blue in color. The white cells are mostly identified by their preference for these dyes. In some cases the cells are even named for the dye that they prefer. For example, cells that prefer a mixture of the acid and basic dye are called neutrophils. In the staining process, a buffer solution is used to control the acid-base balance of the stain. This is a most important function. If the buffer solution is too acid it makes the acid dye too bright and the basic dye too faint. On the other hand, if the buffer solution is too basic it makes the basic dye too bright and the acid dye too faint. In either case, the result is a poorly stained slide. The acid-base balance of a solution is measured by its pH value. A buffer solution should have a pH value between 6.4 and 6.8. This allows the best color contrast between acid and basic dyes.
- f. When optimal staining conditions exist, Wright's stain is very satisfactory and easily differentiates cells. The eosin component stains cell cytoplasm, and the methylene blue component stains nuclear material, granules, and inclusions. Both stains oxidize rapidly because they are in alkaline solution. Giemsa, a purified polychrome stain, is added to compensate for this defect by maintaining the azurophilic staining property of the mixture.

#### 3-6. SLIDE METHOD

a. **Principle**. A small drop of blood is placed near one end of a clean glass slide. Using a second slide as the spreader, the blood is streaked into a thin film and allowed to dry. It is then fixed and stained with modified Wright's stain.

## b. Equipment.

- (1) Venipuncture or finger puncture material.
- (2) Clean glass slides.

## c. **Reagents**.

- (1) Wright's stain buffer.
- (a) Solution A. Dissolve 9.47 grams disodium phosphate (Na<sub>2</sub>HPO<sub>4</sub>) (dibasic), anhydrous, in a 1-liter volumetric flask containing about 750 ml of distilled water. Dilute to the mark with distilled water.

- (b) Solution B: Dissolve 9.08 grams potassium acid phosphate (KH<sub>2</sub>PO<sub>4</sub>) (monobasic), anhydrous, in a 1-liter volumetric flask containing about 750 ml of distilled water. Dilute to the mark with distilled water.
- (c) Mix 27 ml of the  $Na_2HPO_4$  solution (Solution A) with 73 ml of the  $KH_2PO_4$  solution (Solution B). quarts to 1 liter. Final pH should be 6.4.

# (2) Wright's stain solution.

- (a) Add 9.0 grams of Wright's powder stain, 1.0 gram of Giemsa powder stain, and 90 ml of glycerin to a mortar of suitable size. Triturate this mixture thoroughly for 15-30 minutes.
- (b) Transfer the mixture to a large open-mouth jar (a kitchen spoon is ideal for this purpose). Cover the jar and incubate the glycerin stain mixture overnight at 37°C.
- (c) Transfer the glycerin/stain mixture to a large brown bottle and add 2,910 ml of acetone-free methanol.
- (d) Age the Wright's stain solution approximately 2 weeks in the dark. Mix daily to assure that the Wright's stain powder is completely dissolved.

**NOTE**: Commercially available Wright's stain is available through the Federal supply system.

#### 3-7. SLIDE PREPARATION

- a. Make a finger puncture or venipuncture in the usual manner.
- b. Touch a drop of blood to a clean glass slide at a point midway between the sides of the slide and a short distance from one end. If a venipuncture is made, use a capillary tube to transfer a drop of blood from the tube to the slide. If a finger puncture is made, dispense the drop of blood from the puncture site after discarding the first drop.

**NOTE**: The drop of blood should be no larger than 1/8 to 3/16 inch in diameter (see figure 3-5, side method for preparation of blood films, a through c.

c. Lay the specimen slide on a flat surface and hold it securely. Place a smooth, clean edge of the spreader slide on the specimen slide at an angle of about 30<sup>0</sup> from the horizontal (see figure 3-5a).




Figure 3-5a. Side method for preparation of blood films: Place spreader slide at an angle of about 30<sup>0</sup> from the horizontal.

d. Pull the spreader slide toward the drop of blood until contact is made within the acute angle formed by the two slides as shown in figure 3-5b.



Figure 3-5b. Side method for preparation of blood films: Contact blood with spreader.

- e. Allow the blood to spread toward the sides of the slide.
- f. Push the spreader slide smoothly and lightly toward the opposite end of the specimen slide, drawing the blood behind it into a thin film (see figure 3-5c).




Figure 3-5c, Side method for preparation of blood films: Finished slide.

- g. Allow the blood film to air-dry completely. Do not blow on the slide in an effort to enhance drying.
- h. Using a lead pencil, write the name (or identification) of the patient in the thick area of the smear. Do not use a wax pencil as it dissolves during the staining process.

#### 3-8. SLIDE STAINING

- a. Cover the slide completely with Wright's stain and allow it to remain on the smear for about 2 minutes to fix the blood cells. The stain should cover the slide but should not be allowed to overflow the edges; the stain must be replenished should it begin to evaporate.
- b. Add an equal volume of Wright's stain buffer directly to the stain and blow the mixture gently to assure maximum mixing. Allow it to remain for about 4 minutes.
- **NOTE**: The times recommended for staining and buffering are approximate and should be adjusted with each fresh batch of stain to give the most satisfactory results.
- c. Using tap water, float off the mixture of stain and diluent from the slide to avoid the deposition of metallic scum on the smear. The scum appears after the addition of the buffer to the stain. Wash the slide thoroughly under cold, slowly-running water.
- d. Air dry the smear and wipe the excess stain from the under surface of the slide.
  - e. Discussion.
- (1) A properly prepared blood smear is margin-free; has no lines, ridges, or holes; is placed centrally on the slide; has an adequate thin area; and has a uniform distribution of leukocytes.
- (2) It is preferable that blood smears not be made from blood containing anticoagulants since the leukocytes change their staining characteristics, develop vacuoles, engulf oxalate crystals, and show nuclear deformities. However, satisfactory slides are made with blood anticoagulated with EDTA.
  - (3) Avoid the following errors:
- (a) Thick films made from an excess amount of blood placed on the slide.
  - (b) Delay in transferring the blood to the slide.

- (c) A spreader slide that has damaged or unpolished ends.
- (d) The use of dirty, dusty, greasy, or scratched slides.
- (4) If slides cannot be stained immediately, they should be dried and then fixed in methyl alcohol for 30 minutes.
- (5) In cases of marked leukopenia, smears can be prepared from the white cell layer ("huffy coat") obtained by centrifuging the blood slowly in a Wintrobe hematocrit tube at 500-800 rpm for 5 minutes.
- (6) It is important that the blood film be completely dried before staining; otherwise the wet areas will wash off the slide.
- (7) Protect blood slides from insects such as flies, cockroaches, etc. They can "clean" raw blood slides very rapidly.
- (8) Protect slides from areas of high humidity. Excessive moisture tends to hemolyze red blood cells.
- (9) Slides should be stained as soon as possible after preparation. White cells tend to become distorted and disintegrate very rapidly, thus causing considerable difficulty in identification.
- (10) Very little actual staining takes place during the fixation stage. Most of the staining actually occurs during the buffering stage.
- (11) During the buffering stage, it is important that only amounts of buffer equal to the stain be added, otherwise there is a tendency to over-dilute, causing the smear to stain weakly.
- (12) After the staining is complete, do not blot the smear but air-dry it. To speed up the drying process, the smear can be placed in the heat of the substage light. It is important that the slide not be heated too intensely or too long since overheating tends to darken the staining reaction.
- (13) A good quality smear should macroscopically pinkish-gray in hue. It should not be blue, green, or red. Microscopically, the red blood cells should be pink to orange and the white blood cells bluish if they display their true staining color.
- (14) If the RBCs are bluish or green, this indicates that the stain is too alkaline. With an alkaline stain, the WBCs stain heavily and generally display fair distinguishing characteristics. However, the heavy stain masks any abnormalities of the RBCs. Heavy staining can be caused by:
  - (a) Blood smears which are too thick.

- (b) Over-staining (prolonged buffer action).
- (c) Evaporation of the methanol in the stain.
- (d) Stain or diluent which is alkaline.
- (e) Alkaline fumes.
- (15) If the red blood cells are bright red, the stain is too acid. In this condition they stain well but the white blood cells (except eosinophilic granules) stain very poorly if at all. Thus, the stain is of no value for differential studies. "Tendency toward acid staining is caused by:
  - (a) Incomplete drying before staining.
  - (b) Insufficient staining (insufficient buffer action).
  - (c) Overdilution of the stain with buffer.
  - (d) Prolonged washing of the slide after staining.
  - (e) Stain or buffer which is acid.
  - (f) Acid tunes.
  - (16) The staining reactions of blood are as follows:

Type of blood cell or component	Good stain	Acid stain	Alkaline stain
Erythrocytes All neclei Eosinophilic granules	Pink to orange Purple-blue Granules red	Bright red Pale blue Brilliant red, distinct	Blue or green Dark blue Deep gray or blue
Neutrophilic Granules Lymphocyte Cytoplasm	Violet-pink Blue	Pale Pale blue	Dark, prominent

(17) The technician should strive for a staining reaction that is neither too alkaline nor too acid. Such a stain gives good distinguishing features for all the cells of the blood system.

- (18) If the staining reaction is excessively alkaline, this can be corrected by decreasing time of staining or neutralizing the stock stain solution with 1 percent acetic acid or 1 percent hydrochloric acid. Add the acid a drop at a time. Check the results after the addition of each drop of acid with trial slides.
- (19) If the staining reaction is excessively acid, this may be corrected by increasing the time of staining or neutralizing the stock stain solution with 1 percent potassium bicarbonate or a weak solution of ammonia water. Add the ammonia water or potassium bicarbonate one drop at a time. Check the results on trial slides after the addition of each drop of the neutralizer.
- (20) Staining reactions can also be varied by adjusting the proportions of disodium phosphate and potassium acid phosphate used in preparing the buffer. For example, increasing the proportion of disodium phosphate will make the buffer more alkaline; reducing it will make the buffer more acid.
- (21) A poorly stained smear can sometimes be saved by washing rapidly with 95 percent alcohol, washing quickly in water, then restraining.

**Continue with Exercises** 

# **EXERCISES, LESSON 3**

**INSTRUCTIONS**: Answer the following exercises by marking the lettered response that best answers the exercise, by completing the incomplete statement, or by writing the answer in the space provided at the end of the exercise.

After you have completed all of these items, turn to "Solutions to Exercises" at the end of the lesson and check your answers. For each exercise answered incorrectly, reread the material referenced with the solution. Some questions have more than one answer, so read them carefully.

ansv	ver,	so read them carefully.					
1.	То	To obtain valid blood test results, specimens must be properly:					
	a.	Collected.					
	b.	Processed					
	C.	Recorded.					
	d.	a and c.					
	e.	a, b, and c.					
2.	Blood counts on venous and capillary blood are nearly the same if the capillary puncture is:						
	a.	Shallow.					
	b.	Sterile.					
	C.	Free-flowing.					
	d.	Located on the finger.					

3.	Valid	blood	counts	cannot b	e.	made	when
$\sim$ .	v ana	21000	CCGIIC	OGI II IOL K	$\sim$	IIIGGG	**!!

- a. Capillary specimens are not taken from a free-flowing sample.
- b. When capillary specimens are obtained from cyanotic or calloused areas.
- c. a and b.
- d. When sources vary as much as 150 to 1550 cells per cu mm from the real value.
- 4. Venous samples are preferred over capillary samples because they:
  - a. Allow for several and repeated hematological examinations.
  - b. Provide a sufficient amount of blood to perform the various laboratory tests needed.
  - c. Provide for less chance of error because operations are made under better conditions and repeated operations are possible.
  - d. a, b, and c.
- 5. Which blood count method would be performed if blood from extensive burn victims was needed?
  - a. Venous.
  - b. Capillary (micro quantitative).
  - c. Neither method.
  - d. Both methods.

- 6. When collecting blood for white blood cell counts, blood obtained from free-flowing areas or areas of local stasis sources can vary as much as:
  - a. 400-1000 cells per cu mm from their real value.
  - b. 800-1200 cells per cu mm from their real value.
  - c. 1000-1300 cells per cu mm from their real value.
  - d. 1000-1500 cells per cu mm from their real value.
- 7. For adults, which veins should be used for venipuncture?
  - a. Veins located in the distal forearm or antecubital space.
  - b. Veins located in the proximal forearm or antecubital space.
  - c. The jugular vein.
  - d. The femoral vein or antecubital space.
- 8. For the elderly or debilitated persons, or those who may have sclerosed or fragile veins, what should you do for the venipuncture?
  - a. Consult with a physician concerning the procedure.
  - b. Take blood from the veins located in the proximal forearm or antecubital space.
  - c. Use the jugular vein.
  - d. Select the femoral vein.

- 9. If blood is needed from infants, which veins should be used for the venipuncture?
  - a. Sagittal sinus area.
  - b. Veins located in the proximal forearm or antecubital space.
  - c. The collapsed vein.
  - d. The jugular or femoral vein. The vein selected should be large, readily accessible, and sufficiently close to the surface to be seen and palpated.
- 10. If venipuncture poses a problem due to the age of the patient, sclerotization due to repeated venipuncture, or any other unusual circumstance, which procedure is to be followed?
  - a. Under some circumstances the technician should withdraw blood from a sagittal sinus, jugular vein, or femoral vein.
  - b. Under no circumstances should a technician withdraw blood from a sagittal sinus, jugular vein, or femoral vein.
  - c. The technician may withdraw blood from a sagittal sinus, jugular vein, or femoral vein.
  - d. Withdraw blood from a sagittal sinus, jugular vein, or femoral vein.
  - e. If the need arises, the technician should withdraw blood from, a sagittal sinus or jugular vein, but not from the femoral vein.
- 11. Which item of equipment is normally used for the collection of blood specimens?
  - a. Isopropyl alcohol, 40 percent.
  - b. Gauze pads, 6 x 6 inches.
  - c. Needle, 1 to 1 1/2 inches, 19-23 gauge.
  - d. Needle, large bevel.

12.	Ve	ins are made easier to enter if:
	a.	The site of puncture is gently slapped.
	b.	The vein is massaged toward the heart.
	C.	a and b.
	d.	The arm hangs down for 4-6 minutes.
	e.	All of the above.
13.		nerally speaking, veins from which group of people tend to collapse more sily; and, therefore, greater care may be needed to select and puncture the n?
	a.	Children.
	b.	Middle-aged people.
	C.	Elderly.
14.	Blo	ood collection instruments should be:
	a.	Glass and disposable.
	b.	Plastic and calibrated.
	C.	Sterile and disposable.
	d.	Aseptic and anticoagulated.
15.	Th	e syringe and needle for venipuncture must be dry to avoid of the red blood cells.
	a.	Hemolysis.
	b.	Coagulation.
	C.	Contamination.
	d.	Hemoglobin reduction.

16.	To prevent an injection of air into the vein, which can be fatal, what must the technician do?				
	a. Use a longer plunger than the syringe.				
	b.	Use a shorter plunger than the syringe.			
	C.	It makes no difference.			
	d.	The plunger must match the syringe.			
17.		atex tubing is used as a tourniquet, how far above the venipuncture site should e secured?			
	a.	1 inch.			
	b.	2 inches.			
	C.	3 inches.			
	d.	4 inches.			
18.	cor	blonged application of a tourniquet may change the concentration of many blood apponents. The maximum period over which a tourniquet should be applied for enipuncture is:			
	a.	1 minute.			
	b.	2 minutes.			
	C.	4 minutes.			
	d.	6 minutes.			

- 19. Besides inspecting and palpating to locate the desired vein for venipuncture, on what other two items should you focus?
  - a. Direction of vein course and estimate its size and depth.
  - b. Direction of vein course and estimate its length and color.
  - c. Direction of vein course and estimate its position and elasticity.
  - d. The vein's thickness, length, and size.
- 20. When preparing for the venipuncture, what should be done with the needle?
  - a. Keep the cap on until ready to draw.
  - b. Place it on a sterile pad.
  - c. Throw it away.
- 21. What are the reasons for inspecting a possible puncture site?
  - a. Estimate the size and depth of the vein (some may be too small or shallow).
  - b. Determine the direction of the vein's course (puncture with the grain, so to speak).
  - c. Palpate the vein (for resiliency).
  - d. All of the above.
- 22. What may not be done once the puncture area is cleansed and excess alcohol wiped off?
  - a. Grasp the forearm with the left hand.
  - b. Straighten the arm.
  - c. Contaminate the area.
  - d. Have the patient make a clenched fist.

23.	Puncture of the Vacutainer stopper is completed immediately:

- a. Before the needle enters the vein.
- After the needle enters the vein.
- c. Before withdrawal of the needle.
- d. After withdrawal of the needle.
- 24. Which way is the needle bevel supposed to be and how is it to be situated at time of entry?
  - a. Bevel side down, parallel with, and alongside the vein.
  - b. Bevel side up, adjunct with, and alongside the vein.
  - c. Bevel side perpendicular to, and close to the vein.
  - d. Bevel side up, parallel with, and alongside the vein.
- 25. After the needle for a venipuncture is withdrawn, what must the patient do?
  - a. Take an iron compound.
  - b. Lie down for 10 minutes.
  - c. Keep his fist clenched for 5 minutes.
  - d. Maintain light pressure on the gauze pad over the site.

- 26. What is an important vacutainer procedure?
  - a. The short needle should be embedded in the stopper, but the needle must not break the vacuum.
  - b. Any needle should be embedded in the stopper, but the needle must not break the vacuum.
  - c. The long needle should be embedded in the stopper, but the needle must not break the vacuum.
  - d. The first needle should be embedded in the stopper, but the needle must not break the vacuum.
- 27. If the multiple needle is used or more than one tube is required for venipuncture, what is to be followed?
  - a. Tighten the tourniquet after the first tube is filled; remove the filled tube and insert the next one.
  - b. Loosen the tourniquet after the first tube is filled; remove the filled tube.
  - c. Release the tourniquet after the first tube is filled; remove the filled tube and insert the next one.
  - d. Tighten the tourniquet after the first tube is filled and insert the next one.
- 28. Why must you be careful not to remove the needle while tubes are being changed?
  - a. The blood will continue to flow.
  - b. The skin may rip.
  - c. a and b may occur separately or at one tine.

29.	vvr	ny is it most important that correct venipuncture technique be practiced? To:
	a.	Avoid unnecessary pain to the patient.
	b.	Prevent tissue damage.
	C.	Secure a good representative blood specimen.
	d.	Prevent contamination of the specimen or infection of the patient.
	e.	All of the above.
	f.	a, b, and d.
30.		nat may occur to the donor if the tourniquet is not removed as early as possible ce the blood starts flowing well?
	a.	Coagulation of the blood.
	b.	Changes in the quantitative values of the blood components.
	C.	Hemolysis of the blood specimen.
	d.	All of the above.
31.		e blood count should be performed within once the blood collected.
	a.	30 minutes.
	b.	3 hours.
	C.	24 hours.
	d.	48 hours.
32.	Th	e site of a capillary puncture should be:
	a.	Warm.
	b.	Cold.

33.	. When is the tourniquet released and removed?			
	a.	When a hematoma begins to form, the first drop of blood appears, or when it is hard to enter the vein.		
	b.	When it is hard to enter the vein, a hematoma begins to form, or the first drop of blood appears, aspiration occurs.		
	C.	When a hematoma begins to form, or it is hard to enter the vein.		
	d.	All of the above.		
34.	Wh	nat will occur if blood is squeezed from a capillary puncture?		
	a.	Infections.		
	b.	Unnecessary pain.		
	C.	Free- flowing punctures.		
	d.	Inaccurate test results.		

35. Which aspect of the fingertip should be used as the site for a capillary puncture?

MD0853 3-31

a. Dorsal.

b. Ventral.

c. Frontal.

d. Lateral.

36.	Sodium	citrate is	a good	anticoagulant	for coagulation	studies because:
-----	--------	------------	--------	---------------	-----------------	------------------

- a. A concentration of one part 0.2 M sodium citrate is used to 9 parts of whole blood.
- b. It binds the calcium of the blood into a soluble complex to prevent coagulation.
- c. It combines cellular constituents in the plasma.
- d. A concentration of one part 0.1 M sodium citrate is used to 15 parts of whole blood.

# 37. Heparin:

- a. Does not alter the size of cellular components.
- b. Dissolves more rapidly than double oxalate salts.
- c. Is least expensive.
- d. Can be used in excessive amounts.
- 38. EDTA ammonium-potassium oxalate, and heparin are commonly used:
  - a. Stains.
  - b. Buffers.
  - c. Fixatives.
  - d. Anticoagulants.
- 39. What are the two basic methods for the preparation of blood smears?
  - a. Cover slip and polychromic stains.
  - b. Slide and acid dye.
  - c. Slide and cover slip.
  - d. Methylene blue and slide.

40.	Но	w many solutions are needed to perform a Wright's stain buffer?
	a.	1.
	b.	2.
	C.	3.
	d.	4.
41.	Wh	nat should be the final pH of the buffer used with Wright's stain?
	a.	6.0.
	b.	6.4.
	C.	6.8.
	d.	7.2.
42.	Wh	nere on a dried blood smear is the name or identification of the patient written?
	a.	Side.
	b.	Middle.
	C.	Thin area.
	d.	Thick area.
43.	bef	e Wright's stain should remain on a blood smear fore the buffer solution is added and should not be allowed to edges.
	a.	2 minutes; overflow.
	b.	4 minutes; extend beyond.
	C.	6 minutes; remain on.
	d.	8 minutes; reach out.

44.	the	blood smears for the differential leukocyte count cannot be stained immediately, by should be dried and then fixedin forin for nutes.
	a.	EDTA; 10.
	b.	Aniline dye; 15.
	C.	Methyl alcohol; 30.
	d.	Buffer solution; 40.
45.	wh	nen Smears for a differential leukocyte count contain a low concentration of ite blood cells, but marked leukopenia, they can be prepared from the layer by slowly centrifuging the blood specimen in a tube.
	a.	Top; volumetric.
	b.	Buffy coat; Wintrobe hematocrit.
	C.	Plasma layer; test.
	d.	Red blood cell layer; Vacutainer.
46.	If a	reas of a blood smear are still wet when staining is to begin, they will:
	a.	Hemolyze.
	b.	Wash away.
	C.	Stain well.
	d.	Stain too heavily.

47.	Increasing the proportion of disodium phosphate in the buffer for Wright's stain will make the buffer more:			
	a.	Red.		
	b.	Blue.		
	C.	Acid.		
	d.	Alkaline.		
48.	Wh	nat is indicated if when staining the slide the RBCs are bluish or green?		
	a.	The stain is too acidic.		
	b.	The WOC stains very lightly.		
	C.	An alkaline stain is observed.		
	d.	The film is too thick.		
49.	Wh	ny should slides be stained quickly after preparation?		
	a.	So buffers will appear unequal.		
	b.	WBC distort and disintegrate quickly.		
	C.	Lines and ridges will appear.		
	d.	Acid fumes will develop.		

50. Which errors should be avoided when staining slides?

- a. Routinely transferring of blood to the slide.
- b. Using an unpolished slide.
- c. Using clean, dust free, and smooth slides.
- d. Using thin films of blood and placing on slides.

**Check Your Answers on Next Page** 

# **SOLUTIONS TO EXECISES, LESSON 3**

- 1. e (para 3-1a)
- 2. c (para 3-1b)
- 3. c (para 3-1b)
- 4. d (para 3-1b)
- 5. b (para 3-1b)
- 6. d (para 3-1b)
- 7. b (para 3-2a)
- 8. a (para 3-2a)
- 9. d (para 3-2a)
- 10. b (para 3-2a)
- 11. c (para 3-2b)
- 12. c (para 3-2a)
- 13. c (para 3-2a)
- 14, c (para 3-2b)
- 15. a (para 3-2c(4))
- 16. d (para 3-2c(4))
- 17. b (para 3-2d(1))
- 18. b (para 3-2d(1), CAUTION)
- 19. a (para 3-2d (3))
- 20. a (para 3-2c(3))
- 21. d (para 3-2d(3))
- 22. c (para 3-2d(5)

- 23. b (para 3-2e(4))
- 24. d (para 3-c(9))
- 25. d (para 3-2d(11), e(7))
- 26. a (para 3-2e(l))
- 27. c (para 3-2e(5))
- 28. a, b, c (para 3-2d (11)CAUTION, e(5)CAUTION)
- 29. e (para 3-2f(2))
- 30. b (para 3-2f(7))
- 31. b (para 3-2f(8))
- 32. a (para 3-3c(1))
- 33. c (para 3-2f(4))
- 34. d (para 3-3a, c(6), d(3)
- 35. d (para 3-3d(4))
- 36. b (para 3-4e)
- 37. a (para 3-4c)
- 38. d (para 3-4a)
- 39. c (para 3-5b)
- 40. b (para 3-6c(1)
- 41. b (para 3-6c(1)(c)
- 42. d (para 3-7h)
- 43. a (para 3-8a)
- 44. c (para 3-8e(4))
- 45. b (para 3-8e(5))

- 46. b (para 3-8e(6))
- 47. d (para 3-8e(20))
- 48. c (para 3-8e(14))
- 49. b (para 3-8e(9))
- 50. b (para 3-8e(3))

**End of Lesson 3** 

## **LESSON ASSIGNMENT**

LESSON 4

Morphology of Blood Cells.

**TEXT ASSIGNMENT** 

Paragraphs 4-1 through 4-18.

**LESSON OBJECTIVES** 

After completing this lesson, you should be able to:

- 4-1. Select the statement that best describes a general rule of cell identification.
- 4-2. Select the stages of erythrocyte development, the order of maturation and series, and characteristics of each stage with variation.
- 4-3. Select the stages of leukocyte development, the order of maturation and series, and characteristics of each stage with variation.
- 4-4. Select the stages of thrombocyte development, the order of maturation and series, and characteristics of each stage.
- 4.5. Select the statement which best describes the occurrence of leukemia and its classifications and types. Select the appropriate signs and symptoms for the stage or type of leukemia.
- 4-6. Select the appropriate signs and symptoms for the stage or type of leukemia.

SUGGESTION

After completing the assignment, complete the exercises at the end of this lesson. These exercises will help you to achieve the lesson objectives.

MD0853 4-1

#### **LESSON 4**

#### MORPHOLOGY OF BLOOD CELLS

#### Section I. GENERAL INFORMATION

#### 4-1. BASIC CONCEPTS OF CELL MORPHOLOGY

- a. In this lesson, all normal and most commonly seen abnormal blood cells are morphologically described. Although general rules for identification are given along with representative photographs and drawings, it is important to realize that no biological entity fits the guidelines precisely.
- b. The present classification of blood cells is man's attempt at identifying stages of maturation by assigning artificial steps to a continuing process. The process is a smooth, continuing one, and therefore no one cell ever precisely fits the criteria for a specific stage. These stages are artificial classifications that exist to simplify identification.

#### 4-2. GENERAL RULES OF CELL IDENTIFICATION

Certain general rules are applied to all cell maturation (hemopoiesis) either in the erythrocyte, leukocyte, thrombocyte, or plasmocyte series. Although these rules are broken by individual cells, they are an aid to classifying cells.

- a. Immature cells are larger than mature cells and become smaller as they mature.
- b. The relative and absolute size of the nucleus decreases as the cell matures. In some cell series the nucleus disappears.
- c. The cytoplasm in an immature cell is quite blue in color and lightens as the cell matures.
  - d. The young nucleus is reddish and becomes bluer as the cell ages.
- e. Nuclear chromatin is fine and lacy (lacelike) in the immature cell It becomes coarse and clumped in the more mature cells.
  - f. If there is doubt in the identity of a cell, classify to the more mature form.

#### Section II. ERYTHROCYTES

#### 4-3. INTRODUCTION

a. In the normal development and maturation (erythropoiesis) of the erythrocytic series, the red blood cell undergoes a graduation of morphological changes. This cell

MD0853 4-2

development is a gradual transition (as noted in ASCP terminology—they are listed from the most immature to mature cells) and six different stages can be identified. The nomenclature used to describe red blood cells is recommended by the American Society of Clinical Pathologists and the American Medical Association. The terms with some of their synonyms are as follows:

# ASCP Terminology

# Rubriblast Prorubricyte Rubricyte Metarubricyte Diffusely basophilic erythrocyte Erythrocyte

## **Synonyms**

Pronormoblast
Basophilic normoblast
Polychromatophilic normoblast
Orthochromatic normoblast
Polychromatic erythrocyte
Normocyte

b. Erythropoiesis is regulated by the intake of substances to build the cells, the storage of these substances, and their proper utilization. When normal erythropoiesis occurs, both the cytoplasm and the nuclei of the cells grow at a synchronized rate. Individual differences in physiology and physical structure of the erythrocyte account for minor morphological changes so often encountered. In certain diseases, these morphological changes may vary to a greater extent. These variations occur in size, shape, staining, and inclusions in the erythrocyte.

## 4-4. ERYTHROCYTIC SERIES

a. **Rubriblast**. (See figure 4-1, Erythrocytes series, a through f).




Figure 4-1a. Erythrocytes series: Rubriblast.

- (1) Size. Fourteen to 19 microns in diameter.
- (2) <u>Nucleus</u>. This cell has a large round-to-oval purple nucleus that occupies most of the cell. The nuclear chromatin is arranged in a close mesh network forming a reticular appearance. There are 0-2 light blue nucleoli present within the nucleus.

(3) <u>Cytoplasm</u>. The cytoplasm is dark blue, granule-free, and limited to a thin rim around the nucleus. There is no evidence of hemoglobin formation.

# b. Prorubricyte.




Figure 4-1b. Erythrocytes series: Prorubricyte.

- (1) <u>Size</u>. Ten to 15 microns in diameter.
- (2) <u>Nucleus</u>. The nucleus is generally round, dark purple, am smaller than the nucleus of the rubriblast. The chromatin is coarse am clumped giving the nucleus a darker stain. Nucleoli are usually not present, but when they are, they appear more prominent than in the rubriblast.
- (3) <u>Cytoplasm</u>. The cytoplasm is royal blue and more radiant than in the rubriblast. Cytoplasmic granules are not present.

# c. Rubricyte.




Figure 4-1c. Erythrocytes series: Rubricyte.

(1) <u>Size</u>. Eight to 14 microns in diameter.

- (2) <u>Nucleus</u>. The nucleus is dark, round or oval, and smaller than the prorubricyte nucleus. The chromatin material is found in dense, irregular clumps. Nucleoli are not present.
- (3) <u>Cytoplasm</u>. The cytoplasm is more abundant than in the precursor cells. It is blue-pink (polychromatic), the pink resulting from the first visible appearance of hemoglobin. Cytoplasmic granules are absent.

# d. Metarubricyte.




Figure 4-1d. Erythrocytes series: Metrarubricyte.

- (1) Size. Seven to 10 microns in diameter.
- (2) <u>Nucleus</u>. This cell has a pyknotic nucleus (a homogeneous blue-black mass with no structure) that is round. This is the main difference between the rubricyte and the metarubricyte.
  - (3) Cytoplasm. The cytoplasm is abundant, bluish-buff to buff pink.
  - e. Reticulocyte (diffusely basophilic erythrocyte).




Figure 4-1e. Erythrocytes series: Reticulocyte.

- (1) Size. Seven to 9 microns in diameter.
- (2) <u>Nucleus</u>. The nucleus is absent.

(3) <u>Cytoplasm</u>. The cytoplasm stains a bluish-buff with Wright's stain and there is no central light pallor as in the erythrocyte. With supravital staining, this cell will show light blue reticulum strands in the cytoplasm.

# f. Erythrocyte.




Figure 4-1f. Erythrocytes series: Erythrocyte.

- (1) Size. Six to 8 microns in diameter.
- (2) Nucleus. The nucleus is absent.
- (3) <u>Cytoplasm</u>. The cytoplasm of the periphery is light orange with a central zone of pallor. The appearance of the central zone of pallor is due to the biconcave morphology of the cell, which allows more light through the center than through the margin areas.

# 4-5. VARIATIONS IN ERYTHROCYTES

# a. Size.

- (1) <u>Anisocytosis</u>. Anisocytosis is a variation in the size of erythrocytes beyond the normal limits. Cells of varying size are seen in the same fields.
- (2) <u>Macrocytes</u>. Macrocytes are erythrocytes larger than nine microns in diameter. These cells may be found in liver disease.
- (3) <u>Microcytes</u>. These erythrocytes are smaller than 6 microns in diameter. These cells are found in thalassemia and other anemias.

#### b. **Shape**.

(1) <u>Poikilocytosis</u>. This term describes a marked variation in the shape of erythrocytes. Poikilocytes can be pear-shaped, comma-shaped, oval- shaped, or various other bizarre forms. These cells are encountered in pernicious anemia and many other types of anemia.




Figure 4-2a. Variations in erythrocytes: Poikilocytosis: Sickle cells.

(2) <u>Sickle cell (Drepanocytes)</u>. Sickle cells are abnormal erythrocytes that assume a crescent or sickle-shaped appearance under conditions of reduced oxygen tension. The presence of sickle cells is an inherited abnormality due to the presence of hemoglobin S. Sickle cell anemia is encountered primarily in Blacks.




Figure 4-2b. Variations in erythrocytes:

a. Metarubricytes Anisocytosis
Polychromasia Spherocytosis

(3) <u>Spherocytes</u>. These are abnormal erythrocytes that are spherical in shape, having a diameter smaller than normal, and a darker stain (without central pallor) than normal erythrocytes. These cells are found in instances of hemolytic anemias and are particularly characteristic of congenital hemolytic jaundice, a hereditary disorder.




Figure 4-2c. Variations in erythrocytes:

Marked Poikilocytosis.

Anisocytosis and Target Cells.

(4) <u>Ovalocvtes</u>. These cells are abnormal erythrocytes that have an oval or "sausage" shape. They can be found in hereditary elliptocytosis.




Figure 4-2d. Variations in erythrocytes: Elliptocytes (Oval Erythrocytes).

- (5) <u>Target cells (leptocytes)</u>. Target cells are erythrocytes that have deeply stained (pink) centers and borders, separated by a pale ring, giving them a target-like appearance. They are associated with liver disease and certain hemoglobinopathies.
- (6) <u>Burr cells</u>. Burr cells are triangular or crescent-shaped erythrocytes with one or more spiny projections on the periphery. These cells are seen in uremia, acute blood loss, cancer of the stomach and pyruvate kinase deficiency.




Figure 4-2e. Variations in erythrocytes:

Crenated RBC Burr Cells
Acanthocytes 2 Leukocytes

(7) <u>Acanthocytes</u>. Acanthocytes are irregularity-shaped erythrocytes with long spiny projections. They are seen in a congenital abnormality characterized by serum concentration of low density (beta) lipoproteins.




Figure 4-2f. Variations in erythrocytes:

- a. Metarubricyte
- b. Target Cell
- c. Crenated RBC.
- (8) <u>Crenated erythrocytes</u>. This condition occurs when blood films dry too slowly and the surrounding plasma becomes hypertonic. There is no pathological significance when they are found in blood smears.
- (9) <u>Schistocytes</u>. These are red blood cell fragments. Frequently these cells have a hemispherical shape (helmet cells).
- (10) <u>Rouleaux formation</u>. This phenomenon is adherence of erythrocytes to one another presenting a stack-of-coins appearance. It occurs in conditions characterized by increased amounts of fibrinogen and globulin.

# c. Staining.

(1) <u>Hypochramia</u>. Hypochramia is a condition in which the normal central pallor is increased due to decreased hemoglobin content. This condition is characteristic of many anemias.




Figure 4-2g. Variations in erythrocytes: Hypochromic macrocytic erythrocytes.

(2) <u>Polychromatophilia</u>. This term describes non-nucleated erythrocytes that show bluish coloration instead of light pink. Polychromatophilia is due to the fact that the cytoplasm of these cells does not mature, resulting in the abnormal persistence of the basophilic cytoplasm of the earlier nucleated stages.

#### d. Inclusions.

(1) <u>Howell-Jolly bodies</u>. These are nuclear remnants found in the erythrocytes of the blood in various anemias. They are round, dark violet granules about one micron in diameter. Generally, only one Howell-Jolly body will be found in any one red cell. However, two or more may sometimes be present. Howell-Jolly bodies generally indicate absent or non-functioning spleen. They occur in megalobastic anemia and in other forms of nuclear maturation defects.




Figure 4-2h. Variations in erythrocytes:

- a. Metarubricyte
- b. Howell-Jolly Bodies
- (2) <u>Cabot's rings (ring bodies)</u>. These are bluish threadlike rings found in the red cells in the blood of patients with severe anemias. They are interpreted as remnants of the nuclear membrane and appear as ring or "figure-eight' structures. Usually only one such structure will be found in any one red cell.




Figure 4-2i. Variations in erythrocytes: a. Cabot Ring.

(3) <u>Basophilic stippling</u>. Round, small, blue-purple granules of varying size in the cytoplasm of the red cell represent a condensation of the immature basophilic substance (see poly-chromatophilia) that normally disappears with maturity. This is known as basophilic stippling. It can be demonstrated by standard staining techniques in contrast to reticulocyte filaments that require a special stain. Stippling occurs in anemias and heavy metal poisoning (lead, zinc, silver, mercury, bismuth) and denotes immaturity of the cell.




Figure 4-2j. Variations in erythrocytes:
a. Basophilic Stippled Erythrocyte.

- (4) <u>Heinz-Ehrlich bodies</u>. These are small inclusions found primarily in those hemolytic anemias induced by toxins. They are round, refractile bodies inside the erythrocyte and are visible only in unfixed smears. It is thought that they are proteins that have been dematured and that they are an indication of erythrocyte injury.
- (5) <u>Siderocytes</u>. These are erythrocytes containing iron deposits. These deposits indicate an incomplete reduction of the iron from ferric to the ferrous state that is normally found in hemoglobin. Prussian blue stain must be used to readily demonstrate these cells.




Figure 4-2k. Variations in erythrocytes:
a. Rubricytes (Pernicious Anemia).




Figure 4-2I. Variations in erythrocytes:
a. Metarubricyte (Pernicious Anemia)

e. **Megaloblastic Erythrocytes**. The development of megaloblastic cells is caused by a deficiency of vitamin  $B_{12}$  or folic acid. Pernicious anemia is a disease considered to be due to a deficiency in vitamin  $B_{12}$  and/or certain related growth factors. With this deficiency, the erythrocytes do not mature normally and are generally larger than normal. The most notable characteristic of this abnormal maturation is a difference in the rates of maturation of the cytoplasm and the nucleus. The development of the nucleus is slower than that of the cytoplasm, so that in the more mature of the nucleated forms a spongy nucleus as well as an exceptionally large size may be observed. Nuclear chromatin in the megaloblast is much finer and is without the clumps observed in the rubriblast. Such development is termed asynchronism. The mature cell is large (about 10 microns) and is termed a megalocyte. The younger cells of this series are named by adding the suffix "pernicious anemia type," that is rubricyte, pernicious anemia type, and so forth.

#### Section III. LEUKOCYTES

#### 4-6. GRANULOCYTIC SERIES

The stages in the normal maturation of the granulocytes are: myeloblast, promyelocyte, myelocyte (neutrophilic, eosinophilic, and basophilic), metamyelocyte (neutrophilic, eosinophilic, and basophilic), and segmented cell (neutrophilic, eosinophilic, and basophilic). As the granulocytes mature, the granules increase in number. These granules later become specific and differ in the affinity for various dyes. Neutrophilic granules do not stain intensely with either dye. Basophilic granules have an affinity for the basic or blue dye. Eosinophilic stain red with an affinity for the acid dye. The criteria for identification of the various stages of the granulocytic series are: size of cell, nucleus-cytoplasm ratio, nuclear shape, number of nucleoli, and the type and size of cytoplasmic granulation.

# a. Myeloblast.

- (1) Size. Fifteen to 20 microns in diameter.
- (2) <u>Nucleus</u>. The nucleus is round or ovoid and stains predominantly reddish-purple. The interlaced chromatin strands are delicate, well defined, and evenly stained. Two or more pale blue nucleoli are demonstrable. The nucleus occupies most of the cell with a nucleus-cytoplasm ratio of 6:1. It is separated from the cytoplasm by a definite nuclear membrane.
- (3) <u>Cytoplasm</u>. The cytoplasm is a narrow, deep blue, nongranular rim around the nucleus.




Figure 4-3a. Granulocytic series: Myeloblast.

# b. Promyelocyte.

- (1) Size. Fifteen to 21 microns in diameter.
- (2) <u>Nucleus</u>. The nucleus is round or ovoid with coarse-clumping, purple chromatin material. One to three oval, light-blue nucleoli are usually present. The nucleoli are less distinct than in the myeloblast. This cell has a nucleus-cytoplasm ratio of 4:1.
- (3) <u>Cytoplasm</u>. The cytoplasm is light purple and contains varying numbers and sizes of dark nonspecific granules that stain red to purplish-blue. The granules usually overlie the nucleus.




Figure 4-3b. Granulocytic series:

- a. Promyelocyte.
- b. Promyelocyte with Auer Body.
- c. **Myelocyte**. In the myelocytic stage the granules are definite and so numerous that frequently they obscure nuclear detail. While promyelocytes are sometimes distinguished as neutrophilic, eosinophilic, or basophilic, the differentiation is generally considered as first occurring in the myelocytic stage.




Figure 4-3c. Granulocytic series: Myelocyte.

# d. Neutrophilic Myelocyte.

- (1) <u>Size</u>. Twelve to 18 microns in diameter.
- (2) <u>Nucleus</u>. The nucleus is round, oval, or flattened on one side. The chromatin strands are light purple, unevenly stained, and thickened. Nucleoli are usually absent. The nucleus is smaller than the earlier cells of this series with a nucleus-cytoplasm ratio of 2:1.
- (3) <u>Cytoplasm</u>. The cytoplasm is bluish-pink and contains a small relatively light area of ill-defined, pink granules, which develop among the dark, nonspecific, azurophilic granules of the promyelocyte. As the myelocyte ages, the dark granules become less prominent and the light-pink-colored neutrophilic granules predominate.

# e. Neutrophilic Metamyelocyte.

- (1) <u>Size</u>. Ten to 18 microns in diameter.
- (2) <u>Nucleus</u>. The nucleus is indented or kidney-shaped. The nuclear chromatin stains dark purple and is condensed into irregular strands. Nucleoli are absent. The nucleus-cytoplasm ratio is approximately 1.5:1.
- (3) <u>Cytoplasm</u>. The cytoplasm is pinkish-blue and covered with many small, light pink granules.




Figure 4-3d. Granulocytic series:

- a. Metamyelocyte.
- b. Band Neutrophil.

# f. Neutrophilic Band.

- (1) Size. Ten to 16 microns in diameter.
- (2) <u>Nucleus</u>. The nucleus is shaped like a horseshoe with a dark pyknotic mass at each pole of the nucleus where the lobes develop. The nucleus-cytoplasm ratio is approximately 1: 2.

(3) <u>Cytoplasm</u>. The cytoplasm contains many small, evenly distributed light pink granules.

# g. Neutrophilic Segmented Cell.

- (1) Size. Ten to 16 microns in diameter.
- (2) <u>Nucleus</u>. The nucleus has definite lobes separated by a very narrow filament or strand. The nucleus-cytoplasm ratio is approximately 1: 3.
- (3) <u>Cytoplasm</u>. The cytoplasm is light pink and the small, numerous, and evenly distributed neutrophilic granules have a light pink color.
- h. **Development of the Eosinophilic Group**. Cells of the eosinophilic group are characterized by relatively large, spherical, cytoplasmic granules that have a particular affinity for the eosin stain. The earliest eosinophil (myelocyte) has a few dark spherical granules with reddish tints that develop among the dark, nonspecific granules. As the eosinophilic cells pass through their various developmental stages, these granules become less purplish-red and more reddish-orange. The dark blue, nonspecific granules, characteristic of the promyelocyte and the early myelocyte stages, disappear. Because the percentage of eosinophils is usually low in bone marrow peripheral blood smears, no useful clinical purpose is served by routinely separating the eosinophils into their various myelocyte, metamyelocyte, band, and segmented categories. On the other hand, in situations such as eosinophilic leukemia in which the eosinophils are greatly increased, an analysis of the incidence of the various stages would be useful in diagnosis.

# i. Mature Eosinophil.

- (I) Size. Ten to 15 microns in diameter.
- (2) <u>Nucleus</u>. The nucleus has definite lobes separated by a very narrow filament or strand. Seldom does an eosinophil have more than two lobes.
- (3) <u>Cytoplasm</u>. The cytoplasm contains bright reddish-orange, distinct granules. The granules are spherical, uniform in size, and evenly distributed throughout the cytoplasm, but rarely overlie the nucleus.




Figure 4-3e. Granulocytic series: Eosinophil.

j. **Development of the Basophilic Group**. These cells have round, indented, band, or lobulated nuclei and are classified according to the shape of the nuclei, as basophilic rnyelocytes, metamyelocytes, bands, and segmented forms. These cells are so few in peripheral blood and bone marrow that there is little clinical value in differentiation of the various maturation stages.




Figure 4-3f. Granulocytic series:

a. Basophil.

b. Neutrophil: Segmented.

# k. Mature Basophils.

- (1) Size. Ten to 15 microns in diameter.
- (2) <u>Nucleus</u>. The nucleus has definite lobes separated by a very narrow filament or strand. The nuclear details are obscured by the cytoplasmic granulation.
- (3) <u>Cytoplasm</u>. The cytoplasm is covered by many blue to black granules. These granules are unevenly distributed and vary in number, size, shape, and color.

#### 4-7. AGRANULOCYTES

Agranulocytes are leukocytes devoid of specific granulation. These cells generally originate in the lymphatic system, but can be found in normal bone marrow. Agranulocytes include the lymphocytic series, monocytic series, and plasmocytic series.

#### 4-8. LYMPHOCYTIC SERIES

The stages in the development of the lymphocytic series are: lymphoblast, prolymphocyte, and lymphocyte. These cells are fragile and can show shape variants. Lymphocytes usually have round contours, blue cytoplasm, and eccentrically located round nuclei. Cells of this series are differentiated on the basis of the nuclear chromatin.

# a. Lymphoblast.

- (1) Size. Ten to 18 microns in diameter.
- (2) <u>Nucleus</u>. The nucleus has an oval or round shape and stains reddishpurple. The nuclear chromatin is fine, well distributed, and coarser than in the myeloblast. Chromatin is condensed at the edges of the nucleus to form a definite nuclear membrane. One or more nucleoli are present. The nucleus is prominent with a nucleus-cytoplasm ratio of 6:1.
- (3) <u>Cytoplasm</u>. The cytoplasm is deep blue with a frequent perinuclear clear zone.




Figure 4-4a. Lymphocytic series:

- a. Lymphoblast.
- b. Lymphocyte.
- c. Smudge Cell.

# b. Prolymphocyte.

- (1) Size. Ten to 18 microns in diameter.
- (2) <u>Nucleus</u>. The nucleus is oval and slightly indented. The nuclear chromatin is coarse, slightly clumped, am dark purple. One light blue nucleolus is usually present. The nucleus- cytoplasm ratio is 5:1.
- (3) <u>Cytoplasm</u>. The cytoplasm varies from light blue to dark blue and it can show a few red-purple ( azurophilic) granules.




Figure 4-4b. Lymphocytic series: .Prolymphocyte.

# c. Lymphocyte.

- (1) <u>Size</u>. The mature cell of this series varies greatly in size. Small lymphocytes are 7 to 9 microns in diameter. The large lymphocytes are 10 to 18 microns in diameter.
- (2) <u>Nucleus</u>. The nucleus is round or oval and can be slightly indented. The nuclear chromatin is markedly condensed, dark purple-blue, and clumped. Nucleoli are absent and a definite nuclear membrane ids present. The nucleus-cytoplasm ratio is approximately 1.5:1.0.
- (3) <u>Cytoplasm</u>. The cytoplasm is light blue to blue with a perinuclear clear zone around the nucleus. A few azurophilic granules can be seen in the cytoplasm of larger lymphocytes.




Figure 4-4c. Lymphocytic series:
Lymphomcyte.
Azurophilic Granulation.

#### 4-9. MONOCYTIC SERIES

The stages in the development of the monocytic series are monoblast, promonocyte, and monocyte. Cells of the series are slightly larger than granulocytes. They are round with smooth margins and seldom show shape variants. The mature monocyte is differentiated from the lymphocyte and metamyelocyte by the very fine, light staining nucleus.

#### a. Monoblast.

- (1) <u>Size</u>. Twelve to 20 microns in diameter.
- (2) <u>Nucleus</u>. The nucleus is round or oval and red-purple in color. The nuclear chromatin is fine and well distributed. One to two nucleoli are present. The nucleus-cytoplasm ratio is 6:1.
- (3) <u>Cytoplasm</u>. The cytoplasm is a clear, deep blue and follow a thin rim around the nucleus.




Figure 4-5a. Monocytic series:

- a. Monoblast.
- b. Stem (Ferrata Cell)

# b. Promonocyte.

- (1) Size. Twelve to 20 microns in diameter.
- (2) <u>Nucleus</u>. The nucleus is oval or indented and light purple. The nuclear chromatin is fine and spongy. One to five nucleoli may represent. The nucleus-cytoplasm ratio is 5:1.
- (3) <u>Cytoplasm</u>. The cytoplasm is gray-blue with fine dust like azurophilic (red-purple) granules.




Figure 4-5b. Monocytic series: Promonocyte.

## c. Monocyte

- (1) <u>Size</u>. Fourteen to 20 microns in diameter.
- (2) <u>Nucleus</u>. The nucleus is round or kidney-shaped, but can be deeply indented or have two or more lobes. One of the most distinctive features of the monocyte is the presence of superimposed lobes, giving the nucleus the appearance of brain-like convolutions. Heavy lines marking the edges of the folds and grooves are

features that are not seen in other cells. Another feature of the nucleus, which is of value in diagnosis, is the tendency for the nuclear chromatin to be loose with light spaces in between the chromatin strands, giving a coarse linear pattern in contrast to the lymphocyte that has clumped chromatin. Nucleoli are absent. The nucleus-cytoplasm ratio is approximately 3:1.

(3) <u>Cytoplasm</u>. The cytoplasm of the monocyte is dull gray-blue while the cytoplasm of the neutrophils in the adjacent fields is definitely lighter and is pink rather than gray-blue. The nonspecific granules of the monocyte are usually fine and evenly distributed, giving to the cell a dull, opaque or ground-glass appearance. In addition to the background of evenly distributed nonspecific granules, there may be a few unevenly distributed larger azurophilic granules. Vacuoles are often demonstrable in the cytoplasm.

Figure 4-5c. Monocytic series:

- a. Neutrophil (Late Band)
- b. Monocyte.

# 4-10. PLASMOCYTIC SERIES

Plasmocytes constitute approximately one percent of the white cells in the normal bone marrow. These cells can represent in the peripheral blood in chronic infections, granulomatous and allergic diseases, and multiple myeloma. The stages of development are: plasmoblast, proplasmocyte, and plasmocyte.

#### a. Plasmoblast.

- (1) <u>Size</u>. Eighteen to 25 microns in diameter.
- (2) <u>Nucleus</u>. The nucleus is large, oval or round, and located off center in the cell. Nuclear chromatin is fine, purplish, and reticulated. There are multiple nuclei that may or may not be visible.




Figure 4-6a. Plasmocytic series: Plasmoblast.

# b. Proplasmocyte.

- (1) Size. Fourteen to 22 microns in diameter.
- (2) <u>Nucleus</u>. The nucleus is ovoid and located eccentrically. The chromatin is purple, coarser, and more clumped. One to two nucleoli are present. The nucleus-cytoplasm ratio is 2:1.
- (3) <u>Cytoplasm</u>. The cytoplasm is pale blue, nongranular. A lighter-staining area in the middle of the cell (in the cytoplasm, next to the nucleus) may become visible. This is termed a hof.




Figure 4-6b. Plasmocytic series: Proplasmocyte.

# c. Plasmocyte.

- (1) <u>Size</u>. Eight to 20 microns in diameter.
- (2) <u>Nucleus</u>. The nucleus is small, ovoid, and eccentrically located. The chromatin is coarse, lumpy, and purple. Nucleoli are not usually present. The nucleus-cytoplasm ratio is 1:2.
- (3) <u>Cytoplasm</u>. The cytoplasm adjacent to the nucleus is lightly stained in contrast to the periphery of the cell that has a high saturation of red and blue dyes. In same cells, the dark cytoplasm has a greenish or larkspur-blue color. The cytoplasm contains multiple small and relatively unstained globules embedded in a bluish-red filamentous matrix. It is the presence of these tapioca-like globules in the dark surrounding medium that gives the plasmocyte its characteristic mottled and foamy appearance and its brilliant translucency. In occasional cells, the globules can be quite prominent and take a red or bluish-red stain. Such globules are called Russell or fuchsin bodies, or eosinophilic globules. Vacuoles of various sizes are frequently demonstrable.




Figure 4-6c. Plasmocytic series: Plasmocyte.

#### 4-11. VARIATIONS OF LEUKOCYTES

Variations of leukocytes occur as a result of abnormal maturation of the nucleus and/or cytoplasm. These variations are induced by leukemic states, infectious diseases, and toxicity. Described below are the most frequently occurring variations.

- a. **Dohle Bodies**. Dohle bodies are light blue or blue-gray, small, round inclusions found in the cytoplasm of neutrophilic leukocytes. The variation may occur in toxic conditions such as severe infections, burns, poisoning, and following chemotherapy.
- b. **Auer Rods**. Auer rods are rods or spindle-shaped, cytoplasmic inclusions. They stain red-purple and are 1 to 6 microns long and less than 1.5 microns thick. They are frequently found in leukemias.
- c. **Toxic Granullation**. Toxic granulation occurs in the neutrophilic metamyelocyte, band, and segmented cells. These granules are distinguished from the normal granulation because they are coarser and stain a dark purple. The variations occur in toxic states, severe infections, and burns.




Figure 4-7a. Variations in leukocytes:

Band Neutrophil: Toxic Granulation.

- d. **Basket Cell**. A basket cell is a ruptured leukocyte that has a network appearance. These cells result from a partial breakdown of the immature and fragile leukocytes. Basket cells are found predominantly in diseases with an acute shift toward immature forms, for example, leukemias.
- e. **Vacuolated Cell**. A vacuolated cell is a degenerated cell with holes or vacuoles in the cytoplasm. Vacuolated cells can be seen in severe infections, poisoning, and leukemias, and in cells that have been in Heller & Paul oxalate too long.

f. **Hypersegmentation**. A normal neutrophilic segmented cell has a nucleus with an average of three lobes or segments. In a hypersegmented cell the nucleus is broken up into five to ten lobes or segments. This cell usually has a larger diameter than a normal neutrophilic segmented cell. Hypersegmentation is often seen in pernicious anemia and folic acid deficiency. They may also be found in chronic infections.




Figure 4-7b. Variations in leukocytes: Neutrophili: Hypersegmented.

- g. **Atypical Lymphocytes**. These lymphocytes are characteristic of infectious mononucleosis but they may also be seen in apparently healthy individuals and those with certain other diseases. Atypical lymphocytes are larger than normal and vary in appearance. Downey and McKinley described three types of atypical lymphocytes, but this classification has no real clinical purpose.
  - (1) Size. Large, up to 20 microns in diameter.
- (2) <u>Nucleus</u>. The nucleus is oval or kidney shaped with very coarse chromatin strands not as lumpy as a normal lymphocyte.
- (3) <u>Cytoplasm</u>. The cytoplasm is blue to dark blue. Often it is vacuolated which gives rise to a foamy appearance.




Figure 4-7c. Variations in leukocytes:
Atypical Lymphocytes:
Infectious Mononucleosis

#### h. L.E. Cells.

- (1) Persons having lupus erythematosus, one of the "collagen" diseases, have an abnormal; plasma protein that causes swelling and breakdown of certain blood cell nuclei in vitro. This degenerated nuclear material attracts phagocytic cells, particularly segmented neutrophils, which engulf this nuclear mass. The resulting phagocyte and inclusion material is termed an "L.E." cell.
- (2) The nucleus of an L.E. cell is adjacent to the peripheral outline of the inclusion material. The inclusion is smooth and silky or light purple and has no visible chromatin network.




Figure 4-7d. Variations in leukocytes: L.E. Cells.

- i. **Rosettes**. Rosette formation is the intermediate stage in the formation of an L.E. cell. A rosette formation consists of neutrophilic leukocytes surrounding free masses of lysed nuclear material.
- j. **Tart cells**. A tart cell, which may be confused with the L.E. cell, contains lysed nuclear material within its cytoplasm. It differs from an L.E. cell because the inclusion retains characteristic nuclear structure. This inclusion is not smooth and has a darker staining periphery. The significance of tart cells is not known but their presence in an L.E. preparation does not signify a positive test for systemic lupus erythematosus.

# Section IV. THROMBOCYTES

#### 4-12. INTRODUCTION

- a. The general pattern of thrombocyte maturation is slightly different from that of leukocyte maturation. The cells of the megakaryocytic series tend to grow larger as they mature until there is cytoplasmic fragmentation (or breaking off) to form the cytoplasmic thrombocytes seen in the peripheral blood.
- b. Azurophilic granulation begins to appear in the second stage *of* development and continues until it almost obscures the nuclear lobes. The nucleus develops from a fine single lobe to multiple ill-defined lobes. The stages in the normal maturation of the megakaryocytic series are: megakaryoblast, promegakaryocyte, megakaryocyte, and thrombocyte.

#### 4-13. MEGAKARYOCYTIC SERIES

# a. Megakaryoblast.

- (1) Twenty to 50 microns in diameter.
- (2) <u>Nucleus</u>. One to two large oval or kidney-shaped nuclei are present. There is a fine chromatin pattern. Multiple nucleoli may be present which stain blue. The nucleus cytoplasm ratio is approximately 10:1.
- (3) <u>Cytoplasm</u>. The cytoplasm is blue, nongranular, and may have small, blunt pseudopods. It is usually seen as a narrow band around the nucleus.




Figure 4-8a. Megakaryocytic series:

Megakaryoblast: Bone marrow.

# b. **Promegakaryocyte**.

- (1) Size. Twenty to 50 microns in diameter.
- (2) <u>Nucleus</u>. One-to-two indented round or oval nuclei are present. They may show slight lobulation. The nuclear chromatin is purple, coarse, and granular. Multiple nuclei are present but may be indistinct.

#### c. Megakaryocyte.

- (1) Size. Forty to 100 microns in diameter.
- (2) <u>Nucleus</u>. Two-to-sixteen nuclei may be visible or the nucleus may show multilobulation. No nucleoli are visible. The nuclear chromatin is purple, coarse, and granular.
- (3) <u>Cytoplasm</u>. The cytoplasm is pinkish-blue in color and very granular. Numerous blue-purple granules begin to aggregate into small bundles that bud off from the cell to become platelets.




Figure 4-8b. Megakaryocytic series: Megakaryocyte.

# d. Thrombocyte (Platelet)

- (1) Size. One to 4 microns in diameter.
- (2) Nucleus. None.
- (3) <u>Cytoplasm</u>. The cytoplasm is light blue to purple and very granular. It consists of 2 parts:
  - (a) The chromomere, which is granular and located centrally.
- (b) The hyalomere, which surrounds the chromomere and is nongranular and clear to light blue.




Figure 4-8c. Megakaryocytic series: Thrombocytes.

#### Section V. THE LEUKEMIAS

#### 4-14. INTRODUCTION

- a. Leukemia is an abnormal, uncontrolled proliferation of one or more of the white-cell-producing cells. It is a disease of the blood-forming tissues that always involves the bone marrow.
- b. The condition was first discovered microscopically by Donne in 1839 but was not recognized as a clinical entity until 1845. Leukemia was commonly called "white blood" because in many cases, the number of leukocytes in the peripheral blood was so great that upon separation of the blood elements, the buffy (or white cell layer) might be ten times greater than normal.

#### 4-15. OCCURRENCE

Leukemia occurs at any age. Often during the first five years of life, and until the age 21, the great majority of cases are acute. From this time until around 45 the tide turns and chronic granulocytic leukemia becomes the predominating form. It is much more difficult to ascertain the relative frequencies of the various forms of acute leukemia. One school of thought is that acute lymphoblastic leukemia predominates in childhood until about the age of puberty and approximately equals the sum of the cases of myeloblastic and monocytic; that monocytic leukemia never occurs in the aged; or that chronic leukemia is never found in children.

#### 4-16. CLASSIFICATION

- a. The leukemias are classified by the duration of the disease; by the number of white cells present in the peripheral blood; and by the predominant cell type found in the peripheral blood and bone marrow. On the basis of disease duration, leukemia may be described as:
- (1) Acute leukemia, a rapidly progressive disease that lasts several days to 6 months.
  - (2) Sabacute leukemia; 2 to 6 months.
- (3) Chronic leukemia; the length of this disease is variable, depending on the age of the patient and the type of cell involved. Most patients live a minimum of 1 to 2 years or more.
- b. On the basis of the leukocyte count in peripheral blood, the leukemias are classified as:
  - (1) Leukemic leukemia; white blood count is greater than 15,000.

- (2) Subleukemic leukemia; white blood count is less than 15,000 WBC per cu mm with immature or abnormal forms of white cells present in the peripheral blood.
- (3) Aleukemic leukemia; white blood count is less than 15,000 WBC cu mm with no immature or abnormal white cells present in the peripheral blood. The distinction between subleukemic and aleukemic leukemia will depend mostly on a thorough search for abnormal cells.
- c. Another classification is based on the type and the maturity of the cells found to be predominant in the blood and bone marrow. In most cases, classification by cell type makes it possible to distinguish between lymphocytic, neutrophilic, and other leukemias. In acute blast cell leukemia a pure population of blast cells is seen, and the cell type usually cannot be identified. Classification according to cellular maturation will often parallel the classification according to duration; acute leukemias show a large number of immature cells while the chronic cases have a significant number of well-differentiated cells. Morphologically, a "subacute" leukemia will usually resemble the acute stage more readily than the chronic condition.

#### 4-17. TYPES OF LEUKEMIA

Several of the more common types of leukemia are discussed below:

- a. **Acute Lymphocytic Leukemia**. At the time of diagnosis, the white count is usually noticeably elevated, with 60 percent or more lymphoblasts and some prolymphocytes present. Severe anemia and thrombocytopenia are characteristic. The lymphoblasts have a round or oval nucleus composed of coarse granular or stippled chromatin and usually contain one or two nucleoli.
- b. **Chronic Lymphocytic Leukemia**. The white count is usually 20,000 to 200,000 per cu mm, with the peripheral blood smear showing 60 to 95 percent lymphocytes. These cells are generally the small type of mature lymphocyte that often show a small cleft or indentation in the shape of the nucleus. Lymphoblasts are generally absent from the peripheral blood but a rare prolymphocyte may sometimes be found. These lymphocytes are somewhat more fragile than normal, resulting in many of the cells appearing to be entirely devoid of cytoplasm. A normocytic, normochromic anemia generally develops as the disease progresses. The platelet count is usually normal.
- c. **Acute Myelogenous Leukemia**. The white blood count usually shows moderate to marked elevation with 60 percent or more of the cells being myeloblasts. Auer rods may or may not be present in the cytoplasm of these cells. Some cases of this disease show micromyeloblasts, a much smaller myeloblast than normal. Severe normocytic, normochromic anemia develops, along with thrombocytopenia. The platelets that are present may be large and bizarre-looking.

- d. **Acute Promyelocytic Leukemia**. In this form of leukemia, the predominant cell in the bone marrow and blood is the promyelocyte. Often, the nucleus of this cell is more immature than usual and the cytoplasmic granules may be large and abnormal-appearing.
- e. **Chronic Myelogenous Leukemia**. The white count is usually 100,000 to 300,000 per cu nm at the time of diagnosis. Less than 10 percent myeloblasts are present in the peripheral blood and there are numerous immature granulocytes. Eosinophils and basophils are commonly increased, and the percentage of monocytes may also show an increase. Mild normochromic anemia is generally present. The platelet count is often increased and large forms of the platelets may be present.
- f. **Acute Monocytic Leukemia of Schilling**. The white cell count is usually moderately elevated with a predominant number of immature monocytes present. The nucleus is even more lacy than in typical monocytes and it may appear to be folded or segmented. One-to-five nucleoli may appear as fading or inconspicuous. The cytoplasm is variable in amount, generally has few to no visible granules, and may have a serrated-shaped border. Anemia and thrombocytopenia are usually present.
- g. **Acute Myelomonocytic Leukemia of Naegeli**. At diagnosis, the white count usually shows moderate to marked elevation. Anemia is commonly found and thrombocytopenia may also be present. The abnormal cell found most commonly in this disorder has characteristics of both the myeloblast and the monocyte. The nucleus is monocytoid with a fine chromatin pattern and appears convoluted or folded. Granules present in the cytoplasm show characteristics of granulocytes. All stages of these cells are present in the bone marrow and the peripheral blood. Auer rods may be present in the blast cell.

#### 4-18. CLINICAL SIGNS

- a. A patient with leukemia may appear to be in more or less perfect health, or, at the extreme, when the disease is well advanced, the picture of long-drawn, progressive ill health is noticeable; emaciation is usually apparent and pallor is well marked.
- b. The major symptoms of leukemia are fever, weight loss, and increased sweating. Enlargement of the liver, spleen, and lymph nodes may occur. The basal metabolic rate is often elevated, and there may be hemorrhagic tendencies if marked thrombocytopenia is present. When splenic enlargement exists, the contrast between protuberant abdomen and the general evidences of weight loss produces a striking and characteristic picture.

- c. In chronic myleocytic leukemia the spleen is usually markedly enlarged at the time the patient presents himself to the physician. The size of the spleen may be extreme and the organ may extend to the ilium and even to the right anterior superior spine. The size does not appear to be related to the total number of leukocytes. The entire course of development of the splenomegaly may be painless, but not infrequently there is a history of distress in the left upper abdominal quadrant or in the posterior portion of the body.
- d. Pain, limitations of movement, and other symptoms that suggest arthritis, gout, or osteomyelitis, (or those due to actual fracture of bone) may occur in acute or in chronic leukemia. These symptoms are common in children. In chronic leukemia, there may be no symptoms associated with the bones and joints. Bone tenderness can be elicited if firm pressure is made over the sternum. The patient is frequently unaware of any tenderness until pressure is applied at the proper point.

**Continue with Exercises** 

# **EXERCISES, LESSON 4**

**INSTRUCTIONS**: Answer the following exercises by marking the lettered response that best answers the exercise, by completing the incomplete statement, or by writing the answer in the space provided at the end of the exercise.

After you have completed all of the exercises, turn to "Solutions to Exercises" at

		of the lesson and check your answers. For each exercise answered incorrectly, ne material referenced with the solution.
1.	Th	e stages of blood cell maturation are:
	a.	Irrelevant categories.
	b.	Artificial classifications
	C.	Obsolete pigeonholes.
	d.	Strictly and universally applicable categories.
2.	As	a general role for cell identification, the cytoplasm in a mature cell is:
	a.	Green.
	b.	Dark green.
	C.	Blue.
	d.	Light orange.
3.		enerally speaking, what are the texture and consistency of the nuclear chromating an immature cell?
	a.	Fine and lacy.
	b.	Course and clumpy.
	C.	Rough and lacelike

MD0853 4-31

d. Fine and crumbled together.

Generally speaking, what is the size of the cell and the texture and consistency of the nuclear chromatin in a mature cell?			
a.	Larger than an immature cell; fine and lacy.		
b.	Smaller than an immature cell; course and clumpy.		
C.	Smaller than an immature cell; rough and lace-like.		
d.	Larger than an immature cell; fine and crumbled together.		
Wh	nich does NOT occur during the development of blood cells?		
a.	Nucleus disappears.		
b.	Nucleus reduces in size.		
C.	Cytoplasm lightens in color.		
d.	Nucleus becomes reddish in color.		
The	he most immature cell in the erythrocytic series is the:		
a.	Rubricyte.		
b.	Rubriblast.		
C.	Prorubricyte.		
d.	Metarubricyte.		
e.	Erythrocyte.		
f.	Diffusely basophilic erythrocyte.		
	the a. b. c. d. The a. b. c.		

7.	Th	The rubricyte cell has:	
	a.	Cytoplasm staining a bluish-buff and a purple nucleus.	
	b.	A large oval, homogeneous blue-black mass for a nucleus.	
	C.	Dense, irregular clumpy chromatin and a small nucleus.	
	d.	Light blue reticulum strands in the cytoplasm, but no nucleus.	
8.	Wh	nich cell has a nucleus and usually a few nucleoli?	
	a.	Rubricyte.	
	b.	Rubriblast.	
	C.	Megakaryocyte.	
	d.	Metarubricyte.	
9.	ln v	which stage of erythrocyte development does hemoglobin first become visible?	
	a.	Rubricyte.	
	b.	Rubriblast.	
	C.	Prorubricyte.	
	d.	Metarubricyte.	
10.	A c	diffusely basophilic erythrocyte is a/an:	
	a.	Erythrocyte.	
	b.	Rubricyte.	
	C.	Reticulocyte.	
	d.	Metarubricyte.	

11.	The	ne metarubricyte has a/an nuc	cleus.
	a.	Pyknotic.	
	b.	Absent.	
	C.	Round and small.	
	d.	Purple.	
12.	The	ne immediate precursor of the polychromatophilic r	normoblast erythrocyte is the:
	a.	Rubricyte.	
	b.	Rubriblast.	
	C.	Prorubricyte.	
	d.	Metarubricyte.	
13.	Wh	hich cell does not have a nucleus?	
	a.	Rubricyte.	
	b.	Reticulocyte.	
	C.	Prorubricyte.	
	d.	Metarubricyte.	
14.	The	ne term normocyte is synonymous with what blood	cell?
	a.	Rubricyte.	
	b.	Metarubricyte.	
	C.	Prorubricyte.	
	d.	Erythrocyte.	

	a.	3.8 Microns.
	b.	5.3 microns.
	C.	11.4 microns
	d.	18.8 microns
16.	Ab a. b.	normal variation in the size of erythrocytes is called:  Hypochromia.  Ovalocytosis
	C.	Anisocytosis.
	d.	Poikilocytosis
17.	Wł	nich cells are triangular in shape and are spiny looking?
	a.	Ovalocytes.
	b.	Sickle.
	C.	Acanthocytes.
	d.	Burr.
18.	In a.	microcytosis, the microcytes are erythrocyte variations that are:  Larger than normal.
	b.	Smaller than normal.
	C.	Abnormally varied in size.
	d.	Abnormally varied in shape.
	u.	Abhormany varieum shape.

15. The average diameter of a normal prorubricyte is:

19.	RBC fragments that are helmet shaped erythrocytes are called:		
	a.	Crenated erythrocytes.	
	b.	Schistocytes.	
	C.	Drepancytes.	
	d.	Poikilocytosis.	
20.		nich cell is particularly characteristic of congenital hemolytic anemia (called molytic jaundice in the text)?	
	a.	Target	
	b.	Crenated erythrocyte.	
	C.	Spherocyte.	
	d.	Siderocyte.	
21.	Wr	nich cell does have an irregular outline?	
	a.	Acanthrocytes.	
	b.	Burr cells.	
	C.	Target.	
	d.	Crenated erythrocytes.	
22. Which cell does NOT indicate a possible hereditary disorder?		nich cell does NOT indicate a possible hereditary disorder?	
	a.	Ovalocyte.	
	b.	Spherocyte.	
	C.	Sickle cell.	
	d.	Crenated erythrocyte.	

23.	An increase of globulin and fibrinogen presents a stack-of-coins appearance for erythrocytes. This is called a/an:		
	a.	Irregularly-shaped erythrocyte.	
	b.	Pale ring.	
	C.	"Sausage" shape.	
	d.	Rouleaux formation.	
24.	In hypochromic erythrocytes, the normal central pallor is increased as a result oflike many		
	a.	Cellular immaturity; nucleated stages.	
	b.	An increased hemoglobin content; sickle cell abnormalities.	
	C.	A decreased hemoglobin content; anemias.	
	d.	Basophilic cytoplasm; mature cells.	
25.	A megaloblastic cell is caused by what deficiency?		
	a.	Vitamin B <sub>6</sub> .	
	b.	Vitamin B <sub>12</sub> .	
	C.	Vitamin B <sub>1.</sub>	
26.	An	immediate precursor of the neutrophilic band cell is the:	
	a.	Myeloblast.	
	b.	Promyelocyte.	
	C.	Neutrophilic myelocyte.	
	d.	Neutrophilic metamyelocyte.	
	e.	Neutrophilic segmented cell.	

27.	In t	he granulocytic series, the immediate precursor of the promyelocyte is the:
	a.	Monoblast.
	b.	Myeloblast.
	C.	Neutrophilic myelocyte.
	d.	Neutrophilic metamyelocyte.
28.	The	e neutrophilic segmented cell belongs to which series?
	a.	Monocytic series.
	b.	Plasmocytic series.
	C.	Erythrocytic series.
	d.	Granulocytic series.
29.	Wh	nich cell is the least mature and stains unevenly?
	a.	Neutrophilic myelocyte.
	b.	Neutrophilic band cell.
	C.	Neutrophilic metamyelocyte.
	d.	Neutrophilic segmented cell.
30. Which cell has two or more blue nucleoli, no cytoplasmic granule nucleus occupying a ratio of 6:1 nucleus-cytoplasm?		nich cell has two or more blue nucleoli, no cytoplasmic granules, and the cleus occupying a ratio of 6:1 nucleus-cytoplasm?
	a.	Myeloblast.
	b.	Promyelocyte.
	C.	Neutrophilic myelocyte.
	d.	Neutrophilic metamyelocyte.

- 31. Which cell has only dark nonspecific granules within the cytoplasm, with the granules overlying the nucleus?
  - a. Myeloblast.
  - b. Promyelocyte.
  - c. Neutrophilic myelocyte.
  - d. Neutrophilic metamyelocyte.
- 32. A myeloblast cell has a:
  - a. Small, relatively light area of pink granules among dark azurophilic granules and has a nucleus that is round, oval, or flattened on one side?
  - b. Nucleus that is indented and a nucleus-cytoplasm ratio of about 1:5:1.
  - c. Narrow, deep blue, nongranular rim around the nucleus.
  - d. Nucleus with narrow filament that separates the nucleus lobes.
- 33. Which cell has a kidney-shaped nucleus and many small, light pink granules within the cytoplasm?
  - a. Myeloblast.
  - b. Promyelocyte.
  - c. Neutrophilic myelocyte.
  - d. Neutrophilic metamyelocyte.
- 34. The normal stages of granulocytes are:
  - a. Myeloblast, myelocyte, and neutrophilic.
  - b. Myeloblast, promyelocyte, and myelocyte.
  - c. Promyelocyte, myelocyte, and esoinophilic.
  - d. Neutrophilic, esoinophilic, and basophilic.

	C.	Plasmocytic.
	d.	Granulocytic.
37.	Th	e lymphocyte has:
	a.	No nucleus.
	b.	A segmented nucleus.
	C.	An indented, round or oval nucleus.
	d.	A spongy, sprawling nucleus.
38.	Wł	nat are the stages of the lymphocytic series?
	a.	Myeloblast, lymphoblast, and basophilic.
	b.	Lymphocyte, myeloblast, and lymphoblast.
	C.	Lymphoblast, lymphocyte, and monocyte.
	d.	Lymphoblast, prolymphocyte, and lymphocyte.
MD0853		4-40

35. Numerous blue to black granules obscure the nucleus of the:

Erythroblast.

b. Mature basophil.

c. Mature eosinophil.

Monocytic.

b. Lymphocytic.

d. Neutrophilic segmented cell.

36. The promonocyte is a part of what leukocyte series?

a.

a.

39.	Azurophilic (reddish-purple) granules may be found in the cytoplasm of:		
	a.	Lymphocytes.	
	b.	Erythrocytes.	
	C.	Mature eosinophils.	
•	d.	Neutrophilic segmented cells.	
40.	Au	er rods are frequently found in:	
	a.	Anemia.	
	b.	Leukemia.	
	C.	Multiple myeloma.	
	d.	Infectious mononucleosis.	
41.	To	xic granulation of neutrophilic cells occurs in:	
	a.	All of the below.	
	b.	Severe infections.	
	C.	Chemical poisoning.	
	d.	Burns.	
42.	Which leukocyte variation is often produced in blood that has been oxalated to long?		
	a.	Vacuoles.	
	b.	Auer rods.	
	C.	Hyposegmentation.	
	d.	Toxic granulation.	

43.	A h	ypersegmented neutrophilic cell has how many segments?
	a.	1-3.
	b.	3 - 4.
	C.	1 - 5.
	d.	5-10.
44.		cuolated cytoplasm is common in the atypicalaracteristic of infectious mononucleosis.
	a.	Monocyte.
	b.	Lymphocyte.
	C.	Plasmocyte.
	d.	Neutrophilic segmented cell.
45.	A segmented neutrophil that has phagocytized a homogeneous mass of nuclear material is called:	
	a.	A rosette.
	b.	An L.E. cell.
	C.	A tart cell.
	d.	A Dohle body.
46.	Pla	telets (thrombocytes) have a diameter of:
	a.	1-4 microns.
	b.	4-6 microns.
	C.	6-8 microns.
	d.	8-10 microns.

- 47. Leukemias, classified by duration, are described by the most progressive to the least rapid by:
  - a. Acute, chronic, and subacute.
  - b. Acute, subacute, and chronic.
  - c. Chronic, acute, and subacute.
- 48. Which type of leukemia has smaller than normal myeloblast, large and bizarre-looking platelets?
  - a. Acute rnyelomonocytic leukemia of Naegeli.
  - b. Acute promyelocytic leukemia
  - c. Chronic myelogenous leukemia.
  - d. Acute myelogenous leukemia.
- 49. With acute monocytic leukemia of Schilling, the patient has:
  - a. A moderately high WBC, a convoluted fine chromatin pattern, 1 to 5 nucleoli with immature granulocyte, and a serrated-shaped border.
  - b. Anemia, a moderately high WBC, a normal platelet count, convoluted fine chromatin patter, 1 to 4 nucleoli with immature monocytes, and a serrated-shaped border.
  - c. a moderately high WBC, a serrated-shaped border, mature lymphocytes, 1 to 5 nucleoli, thrombocytopenia, and immature granulocytes.
  - d. Cytoplasm with few visible granules; a moderately high WBC; 1 to 5 nucleoli with immature, lacy monocytes; and a serrated-shaped border.

- 50. In advanced leukemia the signs and symptoms include:
  - a. Enlarged lymph nodes.
  - b. Spleen extending to the ilium.
  - c. Fever.
  - d. Weight loss.
  - e. Increased sweating.
  - f. Elevated metabolic rate.
  - g. All of the above.

**Check Your Answers on Next Page** 

# **SOLUTIONS TO EXERCISES, LESSON 4.**

- 1. b (para 4-1b)
- 2. c. (para 4-2c)
- 3. a (para 4-2e)
- 4. b (para 4-2a, e)
- 5. d (para 4-2b-d)
- 6. b (para 4-3a; figure. 4-1)
- 7. c (para 4-4c)
- 8. b (para 4-4a(2))
- 9. a (para 4-4c(3))
- 10. c (para 4-4e)
- 11. a (para 4-4d(2))
- 12. c (para 4-3a)
- 13. b (para 4-4e(2))
- 14. d (para 4-3a)
- 15. c (para 4-4b(1)
- 16. c (para 4-5a(l))
- 17. d (para 4-5b(6))
- 18. b (para 4-5a(3))
- 19. b (para 4-5b(9))
- 20. c (para 4-5b(3))
- 21. c (para 4-5b(5))

- 22. d (para 4-5b(8))
- 23. d (para 4-5b(10))
- 24. c (para 4-5c(1))
- 25. b (para 4-5e)
- 26. c (para 4-6)
- 27. b (para 4-6)
- 28. d (para 4-6)
- 29. a (para 4-6)
- 30. a (para 4-6a(2), (3))
- 31. b (para 4-6b(3))
- 32. c (para 4-6a)
- 33. d (para 4-6e(2), (3))
- 34. b (para 4-6)
- 35. b (para 4-6k(2), (3))
- 36. a (para 4-9b)
- 37. c (para 4-8c(2))
- 38. d (para 4-8)
- 39. a (para 4-8c)(3))
- 40. b (para 4-11b)
- 41. a (para 4-11c))
- 42. a (para 4-11e)
- 43. d (para 4-11f)
- 44. b (para 4-11g)

- 45. b (para 4-11h(1))
- 46. a (para 4-13d(1))
- 47. b (para 4-16a)
- 48. d (para 4-17c)
- 49. d (para 4-17f)
- 50. g (para 4-18)

End of Lessson 4

# **LESSON ASSIGMMENT SHEET**

**LESSON 5** Manual Cell Counts.

**TEXT ASSIGNMENT** Paragraphs 5-1 through 5-6.

**LESSON OBJECTIVES**After completing this lesson, you should be able to:

be able to:

5-1. Select the statement that best describes the materials, procedures, and calculations needed to perform manual counts for blood specimens (WBC, total eosinophil, and reticulocyte cell counts).

5-2. Select the statement which best describes the materials, procedures, and calculations needed to perform manual counts for other body fluids (cerebrospinal fluid cell count and spermatozoa).

5-3. Select the statement which best describes the materials and procedures to perform a semen analysis.

**SUGGESTION** After completing the assignment, complete the

exercises of this lesson. These exercises will help you

to achieve the lesson objectives.

#### LESSON 5

#### MANUAL CELL COUNTS

# Section I. MANUAL COUNTS, BLOOD SPECIMENS

#### 5-1. INTRODUCTION

- a. Blood cells are subject to quantitative variations as well as the qualitative variations described in Lesson 4. Some diseases stimulate the production of blood cells while others prevent or diminish the production of blood cells. For this reason a cell count gives valuable information to the physician concerning his patient's condition. Furthermore, in the case of the leukocyte count, the total count is necessary to calculate absolute counts for each type of leukocyte. This is done by multiplying the total count by the percentage of the particular cell type.
- b. Cell counts can be performed by a variety of methods. Erythrocytes and leukocytes are counted by manual methods or automated methods. Other cell counts are performed only by manual methods. It is important when performing a cell count to maintain good quality control. Great care should be taken when performing any cell count.
- c. The following paragraphs outline procedures for white blood cell (WBC) count, total eosinophil count, and reticulocyte count. The WBC counts are routinely done; they are performed either by the hemacytometer method (manually) or by automated methods. Total eosinophil counts are performed by a hemacytometer method (manually) using special diluting fluids to accentuate these cells. Reticulocytes are demonstrated by using a supravital stain. Semen analysis and cerebrospinal fluid (CSF) counts use a hemacytometer to perform the procedure as well. They are included in the next section.

# 5-2. WHITE BLOOD CELL COUNT

- a. **Principle**. A sample of whole blood is mixed with a weak acid solution that lyses nonnucleated red blood cells. Following adequate mixing, the specimen is introduced into a counting chamber where the white blood cells (leukocytes) in a diluted volume are counted.
- b. **Reagent**. White-count diluting fluid. Either of the following diluting fluids may be used:
- (1) <u>Two percent acetic acid</u>. Add 2 ml glacial acetic acid to a 100 ml volumetric flask. Dilute to the mark with distilled water.
- (2) One percent hydrochloric acid. Add 1 ml hydrochloric acid to a 100 ml volumetric flask. Dilute to the mark with distilled water.

#### c. Procedure.

- (1) Draw well-mixed capillary or venous blood exactly to the 0.5 mark in a white blood cell diluting pipet. This blood column must be free of air bubbles.
- (2) Wipe the excess blood from the outside of the pipet to avoid transfer of cells to the diluting fluid. Take care not to touch the tip of the pipet with the gauze.
- (3) Immediately draw diluting fluid to the "11" mark while rotating the pipet between the thumb and forefinger to mix the specimen and diluent. Hold the pipet upright to prevent air bubbles in the bulb.
- (4) Mix the contents of the pipet for 3-5 minutes to ensure even distribution of cells. Expel unmixed and relatively cell-free fluid from the capillary portion of the pipet (usually 4 drops).
- (5) Place the forefinger over the top (short end) of the pipet, hold the pipet at a 45<sup>o</sup> angle, and touch the pipet tip to the junction of the cover glass and the counting chamber.
- (6) Allow the mixture to flow under the cover glass until the chamber is completely charged. Similarly, fill the opposite chamber of the hemacytometer.

**NOTE**: If the mixture overflows into the moat or air bubbles occur, clean and dry the chambers, remix the contents of the pipet, and refill both chambers.

- (7) Allow the cells to settle for about 3 minutes. Under low-power magnification and reduced light, focus on the ruled area and observe for even distribution of cells.
- (8) Count the white cells in the four 1 sq mm corner areas corresponding to those marked A, B, C, and D of Figure 5-1 in each of two chambers.
- (9) Count all the white cells lying within the square and those touching the upper and right-hand center lines. The white cells that touch the left-hand and bottom lines are not to be counted. In each of the four areas, conduct the count as indicated by the "snake-like" line in figure 5-1. A variation of more than 10 cells between any of the four areas counted or a variation of more than 20 cells between sides of the hemacytometer indicate uneven distribution and require that the procedure be repeated.




Figure 5-1. Hemacytometer counting chamber (WBCs). Areas marked A, B, C, and D are used to count white blood cells.

#### d. Calculation.

- (1) Routinely, blood is drawn to the 0.5 mark and diluted to the 11 mark with WBC diluting fluid. All the blood is washed into the bulb of the pipet (which has a volume of 10). Therefore, 0.5 volumes of blood are contained in 10 volumes of diluting fluid. The resulting dilution is 1:20. (These figures are arbitrary and refer strictly to dilution and not to specific volumetric measurements.)
- (2) The depth of the counting chamber is 0.1 mm and the area counted is 4 sq mm (4 squares are counted, each with an area of 1.0 sq mm therefore, 4 X 1.0 sq mm = a total of 4 sq mm). The volume counted is: area X depth = volume. Four sq mm X 0.1 mm = 0.4 cu mm.

(3) The formula is as follows:

Average number of chambers (2) WBCs counted X dilution (20)

WBCs per cu mm =

Volume (0.4)

(4) For example:

First Chamber	Second Chamber
Cells counted in	Cells counted
each square	in each square
35	45
40	37
44	36
_39_	44
158 WBCs counted	162 WBCs counted

Average of the 2 chamber counts. Total the amount and divide by 2:

Then:

### e. Sources of Error.

- (1) Improper collection of blood specimens causes variable results.
- (2) Wet or dirty pipets.
- (3) Poor condition or inaccurate calibration of pipets. Pipets must be in good condition and calibrated to have maximum error of + or -1percent.
- (4) Poor pipetting technique causes high or low counts. Poor pipetting technique includes:
  - (a) Undershooting desired line with blood or diluting fluid.
  - (b) Overshooting desired line with blood or diluting fluid.
  - (c) Air bubbles in the column on bulb.
  - (d) Failure to wipe tip free of blood.

- (e) Too slow manipulation following the withdrawal of the specimen thus, allowing some of the blood specimen to coagulate.
  - (f) Failure to mix the blood and diluent properly.
- (5) Failure to expel 2 or 3 drops in the pipet tips before charging the hemacytometer.
- (6) Overfilling the chamber of the hemacytometer, which causes erroneously high counts.
  - (7) Wet or dirty cover glasses and hemacytometers.
- (8) Uneven distribution of cells in the counting chamber causes erroneous results.
  - (9) Inaccuracy or carelessness in marking counts.
  - (10) Diluent that which is cloudy or contains debris.
  - (11) Failure to mix anticoagulated blood thoroughly before use.

#### f. Discussion.

- (1) The available error when four large squares are counted is  $\pm 20$  percent. Counting eight large squares decreases the error to  $\pm 15$  percent.
- (2) The importance of clean, dry diluting pipets cannot be stressed too much as the greatest source of error in the counting of WBC is the use of wet and/or dirty pipets.
- (3) The counting chamber must be scrupulously clean and free of debris that might be mistaken for cells.
- (4) The minimum blood sample recommended for performing routine white blood cell counts is that obtained using one pipet and counting two chambers as previously outlined.
- (5) In cases where the WBC count is exceptionally high, as in leukemia, the dilution should be made in the red blood cell diluting pipet. The blood is drawn to the "1.0" mark and the diluting fluid is drawn to the "101" mark. The resulting dilution is 1:100.
- (6) In cases of leukopenia, the white pipet should be filled to the "1.0" mark and diluted to the "11" mark with 2 percent acetic acid. The resulting dilution is 1:10.

(7) If nucleated erythrocytes are present, the count is corrected by the following formula:

The percent nucleated erythrocyte is obtained from the differential count, which is discussed in another subcourse.

# g. Normal Values.

- (1) Adults (both sexes): 4,500-11,500 WBCs per cu mm.
- (2) Childhood: 6,000-14,000 WBCs per cu mm.
- (3) Birth: 9,000-30,000 WBCs per cu mm.

### h. Unopette Procedure for White Blood Cell Count.

- (1) Follow procedure described in para 2-5 for blood dilution.
- (2) Prepare diluted specimen for count.
  - (a) Mix diluted blood by inverting reservoir to resuspend cells.
- (b) Convert to dropper assembly by withdrawing pipet from reservoir and reseating securely in reverse position.
- (c) To clean capillary more, invert reservoir, gently squeeze sides, and discard first three to four drops.
- (d) Carefully load hemacytometer with diluted blood by gently squeezing sides of reservoir to expel contents until chamber is properly filled.
- (e) Place hemacytometer in moisture chamber, let stand for 3 to 5 minutes (10 minutes for platelets) to allow cells to settle.
  - (3) Counting and calculation of leukocytes.
- (a) Under 100X (low power) magnification, count leukocytes in all nine large squares of the counting chamber.
- (b) Add 10 percent of count to total number of cells counted. This step simplifies the calculation that actually entails dividing the number of cells by the number of squares counted and multiplying by 10 to correct for the depth of the chamber.

(c) Multiply this figure by 1000 to get total leukocyte count.

### 5-3. TOTAL EOSINOPHIL COUNT

- a. **Principle**. A sample of blood is diluted with a solution that selectively stains the eosinophils and eliminates all other leukocytes and erythrocytes from view. Following mixing, the specimen is introduced into the counting chamber and the number of eosinophils in a known volume of blood is counted.
- b. **Reagent**. Pilot's Solution: Add 50 ml propylene glycol, 40 ml distilled water, and 1.0 ml sodium carbonate (10 percent aqueous) to a 150-ml beaker. Mix well and filter. Add 10 ml of phloxine (1 percent aqueous) just prior to use. The addition of 10 units of heparin to this solution prevents clumping of cells.

#### c. Procedure.

- (1) Draw capillary (or venous) blood to the "1.0" mark in each of two white cell diluting pipets.
  - (2) Draw Pilot's solution to the "11" mark of each pipet.
- (3) Mix by gently shaking the pipets for 30 seconds. Prolonged and harsh shaking will tend to cause rupturing of the eosinophils.
  - (4) Expel the cell-free liquid from the capillary portion of the pipets.
- (5) Using one pipet, charge both chambers of a hemacytometer and with the other pipet charge both chambers of the second hemacytometer.
- (6) Allow both hemacytometers to stand for 15 minutes to permit staining of the eosinophils. To prevent evaporation, the hemacytometers are placed on a damp towel and covered with Petri dish covers.
- (7) Under low-power magnification, count the red-stained eosinophils in the entire ruled area (9 sq mm) each of the four chambers (a total area of 36 sq mm). The chamber has a depth of 0.1 mm so the total vo1ume is 3.6 cu mm.

#### d. Calculations.

	Number of eosinophils counted X dilution (10)
Number of eosinophils =	
Per cu mm	volume (3.6)

e. Source of Error. See paragraph 5-2e.

#### f. Discussion.

- (1) Fuchs-Rosenthal or levy hemacytometer is preferable to a standard counting chamber since its greater volume (4.0 X 4.0 X 0.2 mm) allows for counting of more cells, thereby reducing the statistical error. Counting two of these chambers is equivalent in accuracy to seven standard chamber counts.
- (2) In eosinopenia, it is necessary to set up more chambers to provide an optional number of cells to be counted.
- (3) The eosinacetone diluting fluids are unsatisfactory and should not be used.
- (4) Estimation of eosinophils on a stained blood smear is too inaccurate for use because of poor cellular distribution.
- (5) The propylene glycol in Pilot's solution renders the erythrocytes invisible, and the sodium carbonate causes lysis of all the leukocytes except the eosinophils. The phloxine stains the eosinophils.
- (6) In the Thorn test an eosinophil count must be made prior to the initiation of the test proper. This establishes the patient's total eosinophil count, to which the response of the adrenal cortex to adrenocorticotropic hormone (ACTH) can be judged. The ACTH is then injected and, at an interval of 4 hours, another eosinophil count is made. The interpretation of this test is as follows:

Normal--approximately a 50 percent drop in eosinophils.

Cushing's disease (hyperadrenalism)--0-30 eosinophils per cu mm.

Addison's disease (hypoadrenalism) no change in eosinophil count.

- (7) Nasal smears are also submitted for eosinophil evaluation. These smears are stained with Wright's stain and examined for the presence of eosinophils.
  - g. **Normal Value**. Normal value: 150-300 eosinophils per cu mm.

### 5-4. RETICULOCYTE COUNT

a. **Principle**. Nonnucleated immature erythrocytes contain nuclear remnants of RNA and the cell is known as a reticulocyte. To detect the presence of this RNA, the red cells must be stained while they are still living. This process is called supravital staining. With supravital staining, the RNA appears as a reticulum within the red cell.

### b. **Reagent**.

- (1) New Methylene Blue Solution. Dissolve 0.5 grams of new methylene blue, 1.4 grams of potassium oxalate, and 0.8 grams of sodium chloride in distilled water. Dilute to 100 ml. Filter before use.
- (2) <u>Brilliant Cresyl Blue Solution</u>. Dissolve 1.0 grams of brilliant cresyl blue in 99 ml of .85 per cent sodium chloride. Filter before use.

#### c. Procedure.

- (1) Place 3 or 4 drops of new methylene blue and 3 or 4 drops of blood (venous or capillary) in a small test tube.
- (2) Mix the tube contents and allow to stand for a minimum of 15 minutes. This allows the reticulocytes adequate time to take up the stain.
  - (3) At the end of 15 minutes, mix the contents of the tube well.
- (4) Place a small drop of the mixture on a clean glass slide and prepare a thin smear.
  - (5) Counterstain with Wright's stain, if desired.
  - (6) Allow smear to air-dry.
- (7) Place the slide on the microscope stage and, using the low power objective, locate the thin portion of the smear in which the red cells are evenly distributed and are not touching each other.
- (8) Switch to oil immersion magnification and count the number of reticulocytes in 5 fields of 200 RBCs.

### d. Calculation.

Number of reticulocytes counted	
	= % reticulocytes
10	

### e. Sources of Error.

(1) Equal volumes of blood and stain give optimum staining conditions. An excess of blood causes the reticulum to understain. An excess of stain usually obscures the reticulum.

- (2) Crenated erythrocytes and rouleaux formation make an accurate count difficult to perform.
- (3) Stain precipitated on erythrocytes causes them to appear as reticulocytes
  - (4) Dirty slides cause uneven spreading.
- (5) The dye solution should have adequate time to penetrate the cell and stain the reticulum.

### f. Discussion.

- (1) Reticulocytes are nonnucleated erythrocytes that exhibit blue reticulum strands within their cytoplasm when stained supravitally. When stained only with Wright's stain, they are buff-pink in color and larger and darker than erythrocytes.
- (2) Reticulocytes serve as an index of the activity of the bone marrow in blood regeneration. As such, these counts are of value in following anti-anemia therapy. Satisfactory response to therapy is evidenced by an increase of reticulocytes in the peripheral blood. Increased reticulocyte counts also occur whenever there is rapid bone marrow activity as in leukemia or blood regeneration associated with hemorrhage or hemolysis. Decreased reticulocyte counts occur in conditions in which the bone marrow is not producing adequate red blood cells, such as aplastic anemia.
- (3) Several methods for staining and counting reticulocytes are in common use. Compared to the use of alcoholic solutions of dye, methods employing saline solutions of new methylene blue can give slightly higher values for reticulocytes. For comparative studies, the same method should be used throughout the work.
- (4) Precipitated stain is often confused with reticulum but can be recognized by its presence throughout the smear and apart from the red cells. Precipitation can be eliminated as a source of error by frequently filtering the stain.
- (5) An alternate method of counting reticulocytes utilizes the Miller disk that is placed inside the microscope eyepiece. This disc consists of 2 squares as shown below in figure 5-2. The area of the smaller square (B) is a tenth that of square A. Therefore, if there are 40 red cells in square A, there should be four red cells present in square B. When employing this method to count reticulocytes, the red cells in square B are counted in successive fields on the slide, until a total of 500 red cells have been counted. At the same time, the reticulocytes in square A are enumerated. At the completion of the count, theoretically, the reticulocytes obtained in this way are divided by 50, in order to obtain the percent reticulocytes present in the blood.




Figure 5-2. Miller disc.

# g. Normal Values.

- (1) <u>Birth</u>. Two and one-half to 6.0 percent, but falls to adult range by the end of second week of life.
  - (2) Adults (both sexes). Five-tenths to 1.5 percent.

# Section II. MANUAL COUNTS, OTHER BODY FLUIDS

### 5-5. CEREBROSPINAL FLUID CELL COUNTS

- a. **Principle**. Cerebrospinal fluid is delivered to a counting chamber and examined microscopically for blood cells. Normally, spinal fluid is clear. If it is cloudy, dilute before charging the counting chamber.
- b. **Reagent**. Cerebrospinal Fluid Diluting Fluid: Add 10 ml of glacial acetic acid and 0.2 grams of crystal violet to a 100-ml volumetric flask. Dilute to the mark with distilled water.

### c. Procedure.

**NOTE**: Set up cell counts on spinal fluids within 30 minutes after withdrawal of the specimen.

- (1) Clear spinal fluid is set up as follows:
- (a) With a capillary pipet introduce a drop of well-mixed spinal fluid into one counting chamber of a hemacytometer.

**CAUTION**: Avoid contamination by careful handling of spinal fluid.

If both leukocy (fresh or crena	(b) Examine the entire ruled area for the presence of cellular elements. rtes and erythrocytes are observed, note the condition of the red cells ated).
	(c) Count all cells in the entire ruled area (0.9 cu mm).
(2)	Turbid spinal fluid is set up as follows:
blood cell dilut	(a) Draw cerebrospinal fluid diluting fluid to "1.0" mark of the white ing pipet.
mark.	(b) Carefully draw a well-mixed specimen of spinal fluid to the "11"
	(c) Shake the pipet for two minutes to mix the specimen.
	(d) Discard the fluid in the capillary portion of the pipet.
minutes.	(e) Charge the counting chamber and allow the cells to settle for five
area (0.9 cu m	(f) Under low-power magnification count all cells in the entire ruled im).
pipet.	(g) Switch to high-power and perform a rough differential count diluting
(3)	For very cloudy spinal fluid a white blood cell dilution is made as follows:
pipet.	(a) Draw spinal fluid to the "0.5" mark in the white blood cell diluting
	(b) Draw cerebrospinal fluid diluting fluid to the 11 mark.
	(c) Perform steps (c) through (g) in paragraph 5-6c(2).
d. Cald	culations.
(1)	Clear spinal fluid:
	Number of cells counted
	volume (0.9) = cells per cu mm

(2) Turbid spinal fluid:

Number of cells counted X dilution (10)

volume (0.9 cu mm) = cells per cu mm

(3) Very clouded spinal fluid:

Number of cells counted X dilution (20)
volume (0.9 cu mm) = cells per cu mm

- e. Source of Error. See paragraph 5-2e.
- f. Discussion.
- (1) If more than 100 leukocytes per cu mm are present, centrifuge the undiluted specimen, make a smear, and stain with modified Wright's stain. Perform a routine differential count and also estimate the ratio of erythrocytes to leukocytes.

**NOTE**: It may be necessary to use egg albumin or cell-free serum to make the sediment adhere to the slide.

- (2) Normally the spinal fluid is water clear. It can be turbid if cell count is 500 or more cells per cu mm. If there is fresh blood with spontaneous clotting, the indications are those of a bloody tap. Xanthochromia develops after subarachnoid hemorrhage has been present for a few hours and is due to disintegration of blood pigments. Xanthochromia may also develop from tumors, abscesses, and inflammation.
- (3) Cell counts above ten are considered to be evidence of intracranial disease. The predominant cell in most viral infections, syphilis, and tuberculous meningitis is the lymphocyte. Bacterial infections due to meningococcus, pneumococcus, and so forth, usually result in a predominance of the neutrophil. Cerebral and extradural abscesses as well as subdural hemorrhages produce a neutrophilic response although bacteria are not demonstrated.
- (4) Biochemical, bacteriological, virological, serological, and hematological ands are all necessary to reflect the true condition of the cerebrospinal fluid. The current laboratory standing operating procedures should give guidance for the most efficient method to accomplish all the necessary ands.
  - g. **Normal Value**: Zero to five cells per cu mm (chiefly lymphocytes).

### 5-6. SEMEN ANALYSIS

- a. **Principle**. Semen analysis involves gross examination (volume, color, turbidity, viscosity, and pH) and microscopic examination (motility and spermatozoa count).
- b. **Reagent**. Spermatozoa Fixative Solution: Add 5 grams sodium bicarbonate (NaHCO<sub>3</sub>) and 1 ml of formalin to a 100 ml volumetric flask. Dilute to the mark with distilled water.
- c. **Collection Instructions**. A physician will usually give the instructions; however, the patient should be reminded of several critical points.
- (1) The patient may be required to abstain from intercourse for a period directed by the physician.
- (2) The specimen is collected in a clean container that has been prewarmed to body temperature.
  - (3) The specimen should be delivered to the laboratory within 30 minutes.
- (4) The specimen must be kept at body temperature (37°C) and not subjected to extremes of heat or cold.

#### d. Gross Examination.

- (1) Record. the time of collection and receipt of the specimen.
- (2) Measure and record the volume.
- (3) Observe and record the color (white, gray, yellow, and so forth), turbidity (clear, opalescent, opaque, and so forth), and viscosity (viscid, gelatin, liquid).
  - (4) Determine the pH with a pH reagent strip and record this.

### e. Motility Examination.

- (1) When the specimen becomes fluid (within 15 to 30 minutes after collection, the semen liquifies by the action of fibrinolysin), place 1 drop on a slide (prewarmed to 37°C) and place a cover slip on it.
- (2) Under high dry power, count motile and nonmotile spermatozoa in two or more areas until a total of at least 200 spermatozoa have been observed. It is necessary to focus through the entire depth of a given field so as to include nonmotile spermatozoa that may have settled to the bottom of the slide. Only those that move forward actively are considered motile. Record the percent of motile spermatozoa seen.

(3) Repeat this procedure in three hours and six hours, using a new drop from the original specimen each time.

# f. Spermatozoa Count.

- (1) Draw the liquefied semen to the "0.5" mark on white blood cell diluting pipet.
  - (2) Draw fixative solution to the "11" mark on the WBC pipet.
  - (3) Let the mixture stand until the mucus dissolves.
  - (4) Shake the pipet thoroughly and charge a hemacytometer.
- (5) Count the spermatozoa in the same manner as you would count white blood cells.
- (6) After counting the sperm, examine the morphology and report the percent of abnormal forms. Morphologically normal sperm are quite uniform in appearance. Any sperm with rounded, enlarged, small, or bilobed heads are abnormal. Abnormal tails are enlarged, small, irregular in length, absent, or multiple. See Figure 5-3 for morphology of spermatozoa.




Figure 5-3. Morphology of spermatozoa.

# g. Calculations.

- (1) Number of sperm

  counted X dilution (20)

  volume (0.4) = sperm per cu mm
- (2) Sperm per cu mm X 1000 = sperm per ml.

### h. Sources of Error.

- (1) Delay in analysis results in a lower percentage of motile forms and a lower count.
- I (2) Temperature extremes cause spermatozoa to die.
  - (3) See paragraph 5-2e for sources of error when counting.

#### i. Discussion.

- (1) Semen analyses are usually performed as part of infertility studies or following a vasectomy.
- (2) Semen analysis can be performed for medico-legal cases involving rape or to support or disprove a denial of paternity on the grounds of sterility.
- (3) Semen is derived from the following: testes, seminal vesicles, prostate, epididymides, vasa deferentia, bulourethral glands, and urethral glands.

### i. Normal Values.

(1) Volume: 1.5-5.0 ml.

(2) pH: 7.2-7.6.

(3) Motility: 60-90 percent.

(4) Spermatozoa Count: 60-150 million per ml.

**Continue with Exercises** 

# **EXERCISES, LESSON 5**

**INSTRUCTIONS**: Answer the following exercises by marking the lettered responses that best answer the exercise, by completing the incomplete statement, or by writing the answer in the space provided at the end of the exercise.

After you have completed all of the exercises, turn to "Solutions to Exercises" at the end of the lesson and check your answers. For each exercise answered incorrectly, reread the material referenced with the solution.

1.	For a WBC count, after drawing blood into the diluting pipet, it is necessary to wipe
	any excess blood from the outside of the pipet in order to avoid:

a.	W	astii	ng	cel	ls.
----	---	-------	----	-----	-----

- b. Dirtying the pipet.
- Dirtying the hemacytometer cover glass.
- d. Transferring cells to the diluting fluid.

2.	When doing a	WBC count.	to what mark	should the	dilutina flui	id be drawn?

- a. 3.
- b. 7.
- c. 11.
- d. 13.

3. When performing a WBC count, which reagents may be used as dilutants?

- a. Acetic acid or sulfuric acid.
- b. Acetic acid or citric acid.
- c. Acetic acid or hydrochloric acid.
- d. Acetic acid or trisphosphoric acid.

4.	For the WBC count, immediately after the contents of the pipet have been mixed for about three minutes, it is necessary to:		
	a.	Use a mechanical mixer.	
	b.	Expel three to four drops.	
	C.	Observe for even distribution of cells.	
	d.	Refill both chambers of the hemacytometer.	
5.	whi	er the WBCs have settled for about three minutes during a manual WBC count, ich powered magnification and lighting arrangements are used to focus on the ed area to observe for even distribution of WBC?	
	a.	Low-power (10X); reduced light.	
	b.	High-power (43X); bright light.	
	C.	Oil immersion (97X); bright light.	
	d.	High-power (43X); reduced light.	
6.	the side	nen counting WBCs, a variation of more than cells between any of four areas counted or a variation of more than cells between es of the hemacytometer indicate uneven distribution and require that the cedure be repeated.	
	a.	6; 12.	
	b.	7; 9.	
	C.	10; 20	
	d.	10; 18.	

	a.	To indicate a "snake-like" line.
	b.	To wash the blood.
	C.	To dilute the blood better.
	d.	To prevent air bubbles in the bulb.
8.	Wł	nich WBCs are counted?
	a.	Those touching the inner left-hand bottom lines.
	b.	All WBCs outside the squares.
	C.	All WBCs within the square and those touching the upper and right hand center lines.
	d.	All WBCs within the square and those touching the upper and left- hand center lines.
9.	Но	w many 1-sq-mm comer areas and chambers are used to count WBCs?
	a.	3; 4.
	b.	4; 2.
	C.	2; 6.
10.	Wh	nich chemical is mixed with whole blood when obtaining a WBC count?
	a.	Sodium chloride.
	b.	Weak acid.
	C.	Weak base.
	d.	Calcium carbonate.

5-20

MD0853

7. Why is the pipet held upright immediately after drawing the diluting fluid to the 11" mark and mixing it with the specimen?

11.		Blood is drawn to the mark and c VBC count.	liluted to the	_ mark for a
	a.	. 0.5, 10.		
	b.	. 0.1,20.		
	C.	. 0.5, 11.		
	d.	. 0.10, 20 .		
12.	Th	he usual blood dilution for the manual WBC	count is:	
	a.	. 1:10.		
	b.	. 1:20.		
	C.	. 1:100.		
	d.	. 1:200.		
13.	Th	he volume is the:		
	a.	. Area X width.		
	b.	. Width X length.		
	C.	. Width X depth.		
	d.	. Area X depth.		

col	unt is:
a.	Average number of WBCs counted X Dilution. = WBCs per cu mm
	Volume Volume
b.	Average number of WBCs counted X Dilution. = WBCs per sq in
	Volume Volume
C.	Average number of WBCs counted X Volume = WBCs per cu mm
	Dilution – WBCs per cu min
d.	Average number of WBCs per cu mm X WBCs counted

Volume

14. Using the hemacytometer counting chamber, the formula for calculating the WBC

- 15. If blood is drawn to the "0.5" mark on a WBC diluting pipet, and diluting fluid to the "11" mark, what is the WBC count of the patient if the average of two chamber counts is 163?
  - a. 3,260 WBCs per cu mm.
  - b. 8,530 WBCs per cu mm.
  - c. 8,150 WBCs per cu mm.
  - d. 10,320 WBCs per cu mm.
- 16. If blood is drawn to the "0.5" mark on a WBC diluting pipet, and diluting fluid to the "11" mark, what is the WBC count of the patient if the average of two chamber counts is 214?
  - a. 4,280 WBCs per cu mm.
  - b. 8,350 WBCs per cu mm.
  - b. 9,984 WBCs per cu mm.
  - c. 10,700 WBCs per cu mm.

"11" r		lood is drawn to the "0.5" mark on a WBC diluting pipet, and diluting fluid to the "mark, what is the WBC count of the patient if the average of two chamber ints is 198?	
	a.	3,960 WBCs per cu mm.	
	b.	9,900 WBCs per cu mm.	
	C.	9,984 WBCs per cu mm.	
	d.	10,540 WBCs per cu mm.	
		at is the maximum allowable error rate when using the four large nacytometer squares in the WBC count?	
	a.	10 percent.	
	b.	15 percent.	
	C.	20 percent.	
	d.	25 percent.	
19.		What is the maximum allowable error rate for the manual WBC count when 8 square areas are employed?	
	a.	10 percent.	
	b.	15 percent.	
	C.	20 percent.	
	d.	25 percent.	

- 20. When performing a WBC count, what is done when the white cell count is exceptionally reduced as in the case of leukopenia?
  - a. The dilution should be made in the red blood cell diluting pipet.
  - b. The blood is drawn to the "1.0" mark and the diluting fluid is drawn to the "11" mark. The resulting dilution is 1:100.
  - c. The white pipet should be filled to the "1.0" mark and diluted to the "11" mark with two percent acetic acid. The resulting dilution is 1:10.
  - d. The count is corrected calculating the observed count x 100 divided by 100 + the percent of nucleated erythrocytes.
- 21. When performing a WBC count, what is done when the white cell count is exceptionally high as in the case of leukemia?
  - a. The dilution should be made in the red blood cell diluting pipet. The blood is drawn to the "1.0" mark and the diluting fluid is drawn to the "101" mark. The resulting dilution is 1:100.
  - b. The white pipet should be filled to the "1.0" mark and diluted to the "11" mark with 2 percent acetic acid. The resulting dilution is 1:10.
  - c. The count is corrected calculating the observed count x 100 divided by 100 + the percent of nucleated erythrocytes.
  - d. Add 10 percent of the count to the total number of cells counted + the percent of nucleated erythrocytes.
- 22. If blood for a WBC count is drawn to the "1.0" mark on a RBC diluting pipet, and diluting fluid to the "101" mark, what is the WPC count of the patient if the average of two chamber counts is 356?
  - a. 22,250 WBCs per cu mm.
  - b. 44,500 WBCs per cu mm.
  - c. 89,000 WBCs per cu mm.
  - d. 178,000 WBCs per cu mm.

		lood for a WBC count is drawn to the "1.0" mark on an RBC diluting pipet, and uting fluid is drawn to the "101" mark, what is the WBC count if the average of chamber counts is 290?
	a.	72,500 WBCS per cu mm.
	b.	145,000 WBCs per cu mm.
	C.	217,500 WBCs per cu mm.
	d.	362,500 WBCs per cu mm.

- 24. If blood is drawn to the "1.0" mark on a WBC diluting pipet, and diluting fluid to the "11" mark, what is the WBC count of the patient if the average of two chamber counts is 153?
  - a. 3,825 WBCs per cu mm.
  - b. 7,650 WBCs per cu mm.
  - c. 38,250 WBCs per cu mm.
  - d. 76,500 WBCS per cu mm.
- 25. If the WBC count is 10,210 and the differential indicates there are 19 nucleated RBCs per 100 WBCs, what is the corrected WBC count?
  - a. 8,650.
  - b. 8,580.
  - c. 9,580.
  - d. 1,021,000.

26.		ne WBC count is 9,640 and the differential indicates there are 14 nucleated Cs per 100 WBCs, what is the corrected WBC count?
	a.	7,390.
	b.	8,256.
	C.	8,456.
	d.	946,000.
27.	Wh	nat is the normal range of the WBC count in adults?
	a.	4,500 -11,500 WBCs per cu mm.
	b.	6,000 -14,000 WBCs per cu mm.
	C.	9,000 -30,000 WBCs per cu mm.
	d.	4.2 -5.4 million WBCs per cu mm.
28.		8 eosinophils are counted in a 36-sq mm area of a hemacytometer using a 1:10 ution, what is the eosinophil count?
	a.	20 eosinophils per cu mm.
	b.	150 eosinophils per cu mm.
	C.	244 eosinophils per cu mm.
	d.	300 eosinophils per cu mm.
29.	Δn	unchanged eosinophil count 4 hours after the injection of ACTH is indicative of:
	,	
	a.	Hypoadrenalism.
	b.	Hyperadrenalism.
	C.	Cushing's disease.
	d.	Normal adrenocortical function.

30.	Which stain is used to evaluate eosinophil nasal smears?		
	a.	Pink.	
	b.	Orange.	
	C.	Supravital.	
	d.	Wright's.	
31.	The reticulocyte is an immature:		
	a.	Rubriblast.	
	b.	Erythroyte.	
	C.	Prorubricyte.	
	d.	Metarubricyte.	
32.		5 reticulocytes are counted in a total of 1,000 erythrocytes, what percentage of iculocytes should be reported?	
	a.	0.015 percent.	
	b.	1.15 percent.	
	C.	2.5 percent.	
	d.	1.5 percent.	
33.	If 86 reticulocytes are counted in a total of 1,000 erythrocytes, what percentage of reticulocytes should be reported?		
	a.	8.6 percent.	
	b.	4.6 percent.	
	C.	66 percent.	
	d.	86 percent.	

34. In tuberculous meningitis, the predominant WBC type usually found in fluid is the:		
	a.	Monocyte.
	b.	Neutrophil.
	C.	Eosinophil.
	d.	Lymphocyte.
	e.	Bone marrow activity.
	f.	Hemolysis.
	g.	All of the above.
	h.	Hemacytometer placement.
35.	If 1	98 cells are counted in an undiluted spinal fluid, what is the cell count?
	a.	2.2 per cu mm.
	b.	22 per cu mm.
	C.	220 per cu mm.
	d.	2,200 per cu mm.
36.	If 4	7 cells are counted in a spinal fluid diluted 1:10, what is the cell count?
	a.	4.7 per cu mm.
	b.	52.2 per cu mm.
	C.	470 per cu mm.
	d.	522 per cu mm.

37.		nite blood cell counts on spinal fluid that are above are ually considered indicative of some type of intracranial disease.
	a.	10 per cu mm.
	b.	15 per cu mm.
	C.	20 per cu mm.
	d.	25 per cu lm1.
38.	ln i	most viral infections, the predominant cell usually found in the spinal fluid is the:
	a.	Neutrophil.
	b.	Basophil.
	C.	Eosinophil.
	d.	Lymphocyte.
39.	In :	subdural hemorrhages, the predominant cell type found in spinal fluid is usually ::
	a.	Lymphocyte.
	b.	Neutrophil.
	C.	Monocyte.
	d.	Segmented lymphocyte.
40.	Th	e neutrophil cell is predominant in which disease or infection?
	a.	Tuberculous meningitis.
	b.	Syphilis.
	C.	Bacterial infections.
	d.	Hodgkin's disease.

41.		Except for the diluting fluid used, the spermatozoa count is almost identical in procedure to the:		
	a.	RBC count.		
	b.	WBC count.		
	C.	Reticulocyte count.		
	d.	Total eosinophil count.		
42.		at is the patient required to abstain from prior to having a Semen analysis lected?		
	a.	Water intake.		
	b.	Using a pre-warmed container.		
	C.	Intercourse.		
	d.	Using the Miller disk.		
43.		at three factors should be observed and recorded during gross and of the nen specimen?		
	a.	Color, viscosity, and temperature.		
	b.	Color, amount of blood, and viscosity.		
	C.	Mucus dissolved, temperature, and color.		
	d.	Viscosity, color, and turbidity.		
44.	Du	ring a motility and of spermatozoa, which cells are considered to be motile?		
	a.	The entire mixture.		
	b.	The entire depth of the field.		
	c.	All active ones moving forward.		
	d.	Only those that are floating.		

	spe	ecimens?		
	a.	Every 15 minutes.		
	b.	In 3 hours and 6 hours.		
	C.	Within 30 minutes of collection.		
46.	Fib	orinolysin causes what type of change to the semen?		
	a.	It solidifies.		
	b.	It is modified with Wright's stain.		
	C.	It causes it to liquify.		
	d.	Nothing.		
47.		men analysis can be performed for cases involving be or in support or denial of paternity on the grounds of		
	a.	Medico-legal; sterility.		
	b.	Medico-vasa deferentia; sterility.		
	C.	Abnormal; epididymides.		
48.	Fro	om which of the following male body parts is semen derived?		
	a.	Testes, seminal vesicles, prostate, epididymides, vasa deferentia, reticulocytes, and urethral glands.		
	b.	Testes, seminal vesicles, prostate, erythrocytes, vasa deferentia, bulbourethra glands, and urethral glands.		

c. Testes, motility spermatozoa, seminal vesicles, prostate, epididymides, vasa

d. Testes, seminal vesicles, prostate, epididymides, deferentia, bulbourethral

deferentia, bulbourethral glands, and urethral glands

45. When should the motility procedure be repeated when examining spermatoza

MD0853 5-31

glands, and urethral glands.

- 49. Which normal value of spermatoza is correct?
  - a. Volume: 0.5-5.0 ml.
  - b. pH: 7.4-7.6.
  - c. Motility: 60-90 percent.
  - d. Spermatozoa Count: 25-150 million per ml.
- 50. Which normal value of spermatoza is correct?
  - a. Volume: 1.5-5.2 ml.
  - b. pH: 7.2-7.6.
  - c. Motility: 60-94 percent.
  - d. Spermatozoa Count: 60-170 million per ml.

**Check Your Answers on Next Page** 

# **SOLUTIONS TO EXERCISES, LESSON 5**

- 1. d (para 5-2c(2))
- 2. c (para 5-2c(3))
- 3. c (para 5-2b)
- 4. b (para 5-2c(4))
- 5. a (para 5-2c(7))
- 6. c (para 5-2c(9))
- 7. d (para 5-2c(3))
- 8. c (para 5-2c(9))
- 9. b (para 5-2c(8))
- 10. b (para 5-2a)
- 11. c (para 5-2c(1), (3); 5-2d(1))
- 12. b (para 5-2d(1))
- 13. d (para 5-2d(2))
- 14. a (para 5-2d(2), (3))
- 15. c (para 5-2d(2), (3))
- 16. d (para 5-2d(2), (3))
- 17. b (para 5-2b(2), (3))
- 18. c (para 5-2f(1))
- 19. b (para 5-2f(1))
- 20. c (para 5-2f(6))
- 21. a (para 5-2f(5))
- 22. c (para 5-2f(5))

- 23. a (para 5-2d(3), (4))
- 24. a (para 5-2d(3), (4))
- 25. b (para 5-2f(7))
- 26. c (para 5-2f(7))
- 27. a (para 5-2g(1))
- 28. c (para 5-3d)
- 29. a (para 5-3f(6))
- 30. d (para 5-3f(7))
- 31. b (para 5-4a)
- 32. d (para 5-4d)
- 33. a (para 5-4d)
- 34. d (para 5-5f(3))
- 35. c (para 5-5d(1))
- 36. d (para 5-5d(2))
- 37. a (para 5-5f(3))
- 38. d (para 5-5f(3))
- 39. b (para 5-5f(3))
- 40. c (para 5-5f(3))
- 41. b (para 5-6f(5))
- 42. c (para 5-6c(1))
- 43. d (para 5-6d(3))
- 44. c (para 5-6e(2))
- 45. b (para 5-6e(3))

- 46. c (para 5-6e(1))
- 47. a (para 5-6i(2))
- 48. d (para 5-6i(3))
- 49. c (para 5-6j(3))
- 50. b (para 5-6j(2))

End of Lesson 5

### LESSON ASSIGMMENT

LESSON 6

Hematocrit, Erythrocyte Sedimentation Rate, and Hemoglobin.

**TEXT ASSIGNMENT** 

Paragraphs 6-1 through 6-14.

**TASK OBJECTIVES** 

After completing this lesson you should be able to:

- 6-1. Select the statement that best describes the procedures used to perform the microhematocrit test.
- 6-2. Select the statement that best describes the procedures used to perform the erythrocyte sedimentation using the Wintrobe-Landsberg and modified Westergren techniques.
- 6-3. Select the statement that best describes hemoglobin, the compounds *of* hemoglobin, and the variations *of* hemoglobin.
- 6-4. State the four basic ways to measure hemoglobin concentrations and correctly describe the cyanmethemoglobin method.
- 6-5. Select the procedures required to perform the detection of hemoglobin S and to demonstrate the sickle cell phenomenon.
- 6-6. Select the statement that correctly describes the purpose of hemoglobin electrophoresis.
- 6-7. Select the statement that correctly describes the procedures to perform a fetal hemoglobin test.

SUGGESTION

After completing the assignment, complete the exercises of this lesson. These exercises will help you to achieve the lesson objectives.

#### **LESSON 6**

# HEMATOCRIT, ERYTHROCYTE SEDIMENTATION RATE, AND HEMOGLOBIN

### Section I. HEMATOCRIT

#### 6-1. INTRODUCTION

- a. The hematocrit or the packed-cell volume is the percentage of the total volume of red blood cells in relation to the total volume of whole blood. The term "hematocrit" is derived from two Greek words: <u>Hemato</u>, meaning blood; and <u>Krites</u>, meaning to judge (then reduced to <u>crit</u>, meaning to separate). The procedure is performed by filling a capillary tube with blood and centrifuging at a constant speed for a constant period of time. The packet cell volume is then measured. The hematocrit can also be determined by automated sequential analyzers but is usually a calculated value.
- b. The hematocrit is the most useful single index for determining the degree of anemia or polycythemia. It can be the most accurate (2-4 percent error) of all hematological determinations. In contrast, the direct red blood cell chamber count has a percent error of 8-10 percent. The hematocrit is, therefore, preferable to the red blood cell count as a screening test for anemia. The values for the hematocrit closely parallel the values for the hemoglobin and red blood cell count.

### 6-2. MICROHEMATOCRIT

a. **Principle**. A capillary tube is filled with whole blood by capillary action to within 1 to 2 cm of the end. The unfilled end is sealed and the tube is centrifuged. After centrifugation, the capillary tube is placed in a reading device and the hematocrit value determined.

### b. Procedure.

- (1) If anticoagulated venous blood is the specimen, fill a plain capillary tube with blood. If blood without anticoagulant is used, fill a heparinized capillary tube with the blood specimen. A heparinized capillary tube is identified by red line on the tube specimen. A heparinized capillary tube is identified by a red line on the tube.
- (2) Allow blood to enter two capillary tubes until they are approximately 2/3 filled with blood. (Air bubbles denote poor technique, but do not affect the results of the test.
  - (3) Seal the unfilled ends of the tubes with clay.
- (4) Place the two hematocrit tubes in the radial grooves of the centrifuge head exactly opposite each other, with the sealed end away from the center of the centrifuge.

- (5) Screw the flat centrifuge head cover in place.
- (6) Centrifuge at 10,000 rpm for 5 minutes.
- (7) Remove the hematocrit tubes as soon as the centrifuge has stopped spinning. Determine the hematocrit values with the aid of a microhematocrit reader. Results should agree within +2 percent. If they do not, repeat the procedure.

**NOTE**: Since there are a variety of readers available, it is necessary that the technician carefully follow the directions of the manufacturer for the particular device utilized.

#### c. Sources of Error.

- (1) Inadequate mixing *of* the blood prior to sampling.
- (2) Improper sealing of the capillary tube causes the blood to blow out of the capillary tube during centrifugation.
- (3) Capillary tubes must be properly identified. Numbered holders for capillary tubes are available. Place the tubes in slots on the holder and record the numbers on the laboratory request slip.
- (4) Improper centrifugation leads to varied results. For good quality control, maintain prescribed centrifuge speed and time.
- (5) Misreading the red cell level by including the buffy coat causes elevated values.

# d. Discussion.

- (1) The microhematocrit technique is advantageous because of speed, and because only a small quantity of blood is necessary for the determination. An additional advantage is the ease with which this procedure is adapted to infants and small children. The microhematocrit technique requires only a simple capillary puncture whereas in the Wintrobe method venous blood oust be used. Another advantage is the use of disposable capillary tubes.
- (2) If the microhematocrits cannot be read promptly, the capillary tubes must be properly identified and placed in a vertical position. Slanting of the cell layer will occur if tubes are left in a horizontal position for more than 30 minutes.
- (3) Hematocrit results may also re obtained or computed through use of an electronic cell counter (Coulter Models). Results of computed hematocrits are generally 1.3 to 3.0 percent higher due to their method of derivation. An accurate erythrocyte count (three simultaneous electronic counts averaged) is multiplied by an accurate MCV

measurement resulting in the area that the red cells would occupy if packed. The computation makes no allowance for the trapped plasma that always occurs to sate degree in manual packing procedures.

#### e. Normal Values.

(1) Birth: 50-62 percent.

(2) Age 1: 31--39 percent.

(3) Adult males: 42–52 percent.

(4) Adult females: 36–46 percent.

#### Section II. ERYGHROCYTE SEDIMENTATION RATE

#### 6-3. ERYTHROCYTE SEDIMENTATION RATE

- a. **Introduction**. When anticoagulated whole blood is allowed to stand for a period of time, the red blood cells settle out frail the plasma. The rate at which the red cells fall is known as the erythrocyte sedimentation rate (ESR). The ESR is affected mainly by three factors: erythrocytes, plasma, and mechanical and technical factors.
- b. **Erythrocytes**. Size and shape of erythrocytes cause the ESR to fluctuate. Microcytes tend to settle slower than normal cells, while macrocytes fall more rapidly. An increase in spherocytes and/or bizarrely-shaped red cells retards the sedimentation rate. With a decreased hematocrit there is less retardation of sedimentation by the erythrocytes themselves and they tend to settle faster. Corrections for anemic blood are available; however, most experts consider this correction useless and invalid. An increased hematocrit (above 55 percent) retards the sedimentation rate.
- c. **Plasma Composition**. In certain diseases plasma proteins, namely fibrinogen and globulin, may be altered, causing rouleaux formation. The speed of the sedimentation corresponds to the length of the rouleaux formation. Increases in fibrinogen or globulin will produce long rouleaux that are difficult to disperse. This leads to a larger mass and an increased sedimentation velocity.
- d. **Mechanical Technical Factors**. It is important that the ESR tube be exactly perpendicular. A tilt of 30 can cause errors up to 30 percent. Vibration or movement of the ESR tube or holding rack can affect the ESR as can large changes in temperature.
- e. **ERYTHROCYTE SEDIMENTATION RATE**. Determination of the ESR is performed by many methods of the more common methods; the Wintrobe-Landsberg and modified Westergen methods are stated in this lesson.

# 6-4. DETERMINATION OF SEDIMENTATION RATE (WINTROBE-LANDSBERG)

a. **Principle**. Anticoagulated blood is placed in a narrow tube. The blood cells settle out of the suspension, leaving clear plasma above them.

The distance that the erythrocytes fall within a given interval of time is measured.

### b. Procedure.

- (1) Draw 5 ml of blood by venipuncture and place in a test tube containing ethylenediaminetetraacetic acid (EDTA) (lavender top vacuum tube).
- (2) Thoroughly mix the blood and anticoagulant by gently inverting the tube several times, being careful not to cause bubbles.
- (3) Draw the blood into the capillary pipet and fill the Wintrobe tube to the "0" mark. This is done by inserting a capillary pipet to the bottom of the Wintrobe tube, while holding it at an angle of 45°. As the tube is filled, slowly withdraw the pipet so that the tip is always just below the level of the blood. The blood count must be free of bubbles.
- (4) Place the filled Wintrobe tube in a rack in an exactly vertical position and note the time and roan temperature.
- (5) At the end of exactly 1 hour, read the level to which the red cells have settled on the descending scale etched on the tube. Each mark equals 1.0 mm while each numbered mark equals 10 mm (1 cm). The figure obtained is reported in mm per hour as the "uncorrected" erythrocyte sedimentation rate.
- (6) If a "corrected" sedimentation is requested, perform a hematocrit. Calculators to correct for sedimentation rate are available in the Federal Supply Catalog.

### c. Sources of Error.

- (1) The blood specimen must be properly mixed with the proper anticoagulant to obtain an undiluted representative sample.
- (2) The test should be set up within two hours after the blood sample is collected to avoid a false low sedimentation rate.
- (3) Increase in temperature accelerates the rate. Desirable temperature range is 22°C to 27° C.
- (4) The tube must be vertical. A 3° variation from the vertical rate accelerates the rate by 30 percent.

- (5) Dirty Wintrobe tubes or capillary pipets can decrease the rate.
- (6) Tubes should be placed free from vibration or disturbance.

### d. Normal View.

(1) Males: Zero to nine mm per hour.

(2) Females: Zero to 20 mm per hour.

(3) Children: Zero to 13 mm per hour.

### 6-5. DETERMINATION OF SEDIMENTATION RATE

a. **Principle**. Well-mixed, whole blood is diluted with 0.85 percent sodium chloride placed in a Westergren pipet, and allowed to stand for exactly one hour. The number of millimeters the red cells fall during this tie period constitutes the ESR result.

#### b. **Procedure**.

- (1) Mix the whole blood and the anticoagulant by gently inverting the tube several times or by placing on a rotator for two minutes.
- (2) Place 0.5 ml of 0.85 percent sodium chloride in a plain 13X 100 mm test tube.
  - (3) Add 2.0 ml of well-mixed whole blood to the test tube.
  - (4) Mix the tube for two minutes.
  - (5) Make certain that the Westergren ESR rack is exactly level.
- (6) Fill the Westergren pipet to exactly the "0" mark, making certain there are no air bubbles in the blood.
- (7) Place the pipet in the rack. Be certain the pipet fits snugly and evenly into the grooves provided for it.
  - (8) Allow the pipet to stand for exactly 60 minutes.
- (9) At the end of 60 minutes, record the number of millimeters that the red cells have fallen. This result is the erythrocyte sedimentation rate in millimeters per hour.
  - c. Sources of Error. See para 6-4c.

#### d. Normal Values.

(1) Males: Zero to 15 mm per hour.

(2) Females: Zero to 20 mm per hour.

(3) Children: Zero to 10 mm per hour.

## e. Discussion.

- (1) The erythrocyte sedimentation rate is a nonspecific test that suggests the possibility of a disease process and tissue damage in the body. It is not diagnostic but is extremely useful in following the course of some diseases.
- (2) The rate is usually increased in inflammatory infections, toxemia, cell or tissue destruction, severe anemia, active tuberculosis, syphilis, acute coronary thrombosis, rheumatoid arthritis, and malignant processes.
- (3) Sickle cell anemia, polycythemia, hypofibrinogenemia, and certain drugs usually decrease the rate.

## Section III. HEMOGLOBIN

### 6-6. GENERAL INFORMATION

Hemoglobin is a conjugated protein composed of the basic protein globin linked to 4 heme molecules. Ninety-eight percent of all the iron found in the blood is contained in hemoglobin. Hemoglobin transports oxygen and carbon dioxide. This important substance reacts with oxygen to form oxyhemoglobin. In the tissues, oxygen is released and reduced hemoglobin is formed. Hemoglobin can react with acids, bases, and oxidizing and reducing agents. It also can exist in a variety of forms. These hemoglobin compounds and variants are discussed briefly in the following paragraphs. For more detailed information, refer to the standard hematological texts.

### 6-7. COMPOUNDS OF HEMOGLOBIN

- a. **Oxyhemoglobin**. Oxygen combines loosely with iron (ferrous state) in hemoglobin. The loosely attached oxygen diffuses into the tissues for oxidative processes. The hemoglobin then binds carbon dioxide and exists as reduced hemoglobin.
- b. **Carboxyhemoglobin**. Hemoglobin combines with carbon monoxide to form carboxyhemoglobin. Carbon monoxide has an affinity 200 times greater for hemoglobin than oxygen does. Hemoglobin in this combination is incapable of oxygen transport.

- c. **Methamoglobin**. This compound is formed when the ferrous state of the heme is oxidized to the ferric state. This compound is incapable of oxygen transport.
- d. **Sulfhemoglobin**. This compound results from the combination of inorganic sulfides and hemoglobin. This compound is incapable of oxygen transport. This is an irreversible reaction.
- e. **Cyanmethemoglobin**. This compound results when methemoglobin combines with the cyanide radical. This compound is used in hemoglobinometry.

### 6-8. VARIATIONS OF HEMOGLOBIN

The variations of hemoglobin occur due to structural differences in the globin protein. These differences are genetically controlled. The normal hemoglobin components are hemoglobin A (HbA), hemoglobin A<sub>2</sub> (HbA<sub>2</sub>), and fetal hemoglobin (HbF). HbA constitutes most of the hemoglobin of a normal adult while HbA<sub>2</sub> constitutes a much smaller amount. HbF is present during the first 4 to 6 months of life and not normally present in adults. Hemoglobin S and hemoglobin C are the most commonly occurring abnormal hemoglobins. Others (D, E, H, etc.) are found in rare occurrences associated with several types of anemia. The various types of hemoglobin are separated by electrophoresis.

### 6-9. HEMOBLOGINOMETRY

The hemoglobin concentration is directly proportional to the oxygen-combining capacity of blood. Therefore, the measurement of the hemoglobin concentration in the blood is important as a screening test for diseases associated with anemia and for following the response of these diseases to treatment. There are four basic ways to measure the hemoglobin concentration: (1) measurement of the oxygen-combining capacity of blood (gasometric), (2) measurement of the iron content (chemical method), (3) colorimetric measurement of specific gravity (gravimetric method), and (4) method is the most widely used. The cyanmethemoglobin method is the method of choice and is recommended by the Technical Subcommittee on Hemoglobinometry of the International Committee for Standardization in Hematology.

## 6-10. CHYANMETHEMOGLOBIN METHOD

a. **Principle**. Blood is diluted with a dilute solution of potassium ferricyanide and potassium cyanide at a slightly alkaline pH. The ferricyanide converts the hemoglobin to methemoglobin. The cyanide then reacts with the methemoglobin to form the stable cyanmethemoglobin. The color intensity is measured in a spectrophotometer at a wavelength of 540 mm. The optical density is proportional to the concentration of hemoglobin.

#### b. Discussion.

- (1) Cyanmethemoglobin is the most stable of the various hemoglobin pigments showing no evidence of deterioration after 6 years of storage in a refrigerator. The availability of prepared standards is a distinct advantage of this technique. All hemoglobin derivatives are converted to cyanmethemoglobin with the exception of sulfhemoglobin.
- (2) This method is highly accurate and is the most direct analysis available for total hemin or hemoglobin iron. Its disadvantage is the use of cyanide compounds, which, if handled carefully, should present little hazard.
- (3) For accuracy in hemoglobin determinations, it is absolutely necessary that the spectrophotometer and Sahli pipets be accurately calibrated.
  - (4) Venous samples give more constant values than capillary samples.
- (5) If the procedure is performed properly, the degree of accuracy is  $\pm 2$  to 3 percent.

### c. Normal Values.

- (1) Infants at birth: 17–23 g hemoglobin per deciliter.
- (2) Childhood: 12–14 g hemoglobin per deciliter.
- (3) Adult males: 14–17 g hemoglobin per deciliter.
- (4) Adult females: 13–15 g hemoglobin per deciliter.

### 6-11. DETECTION OF HEMOGLOBIN S AND NON S SICKLING HEMOGLOBINS

a. **Principle**. Erythrocytes are introduced into a phosphate buffer solution containing a reducing agent and hytic agent. The red cells are lysed and the hemoglobin is reduced. Reduced sickling types of hemoglobin are insoluble in phosphate buffer and turbidity results. On addition of urea, hemoglobin S dissolves.

# b. Reagents.

- (1) Stock phosphate buffer solution. Dissolve 160.48 grams anhydrous potassium dihydrogen phosphate ( $KH_2PO_4$ ), AR, and 28.88 grams anhydrous potassium monohydrogen phosphate ( $K_2HPO_4$ ) and AR in a 1-liter volumetric flask containing 500 ml of distilled water. Dilute to 1-liter with distilled water.
- (2) Dithionite reagent. Add 20.0 grams dithionite ( $Ha_2S_4O_4$   $2H_2O$ ) and 0.25 grams saponin to a 100-ml volumetric flask. Add 80 ml of stock phosphate buffer

solution. Mix well. Dilute to the mark with distilled water. This reagent remains stable under refrigeration at 4°C for approximately 1 month. See figure 6-1.




Figure 6-1. Dithionite tube test interpretation.

WARNING: Sodium dithionite reagent may ignite if allowed to become damp. Use only clean dry utensils in handling.

(3) Urea (USP).

# c. Procedure.

- (1) Pipet 2 ml of dithionite reagent in 12 x 75 mm test tube.
- (2) Add 0.02 ml of well-mixed anticoagulated blood (collected in EDTA).
- (3) Mix the contents and al1a.v to stand at roan temperature for 5 minutes (see figure 6-1).
- (4) After 5 minutes examine the tube for turbidity against a lined reader. Hemoglobin 5, if present, produces turbidity in the tube .

## d. Sources of Error.

- (1) The use of 10 x 75 mm test tubes could cause a false negative result.
- (2) The dithionite reagent has a limited stability. The freshness of this reagent must be checked with positive and negative controls. The test should show the blue-pink color of reduced hemoglobin and adequate lysis of erythrocytes.

- (3) Unstoppered tubes containing dithionite reagent decompose when left out at room temperature.
- (4) False negative results could occur if the blood sample for testing is drawn within 4 months of transfusion.

#### e. Discussion.

- (1) Hemoglobin S is an inherited type of hemoglobin found primarily in blacks and people from Mediterranean areas.
- (2) The degree of erythrocyte sickling is dependent on the concentration of hemoglobin. SS, SC, and SD cells sickle more rapidly than AS cells. Newborns with sickle cell anemia have erythrocytes more resistant to sickling due to the presence of hemoglobin F.
- (3) The dithionite test also detects other sickling types of hemoglobin. Urea causes hemoglobin S (and structural variants of hemoglobin S) to dissolve. Other hemoglobins remain turbid in the presence of urea.
- (4) This test is a rapid screening test for hemoglobin S. All positive dithionite tests should be electrophoresed for confirmation.
- f. **Interpretation**. Hemoglobin S causes turbidity in the tube. Hemoglobin A is soluble in the phosphate buffer.

### 6-12. DEMONST4RATION OF THE SICKLE CELL PHENOMENON

- a. **Principle**. Erythrocytes of persons with sickle cell anemia or trait will assume a sickle shape when the oxygen tension is lowered. This may be demonstrated by mixing a drop of blood with a reducing agent such as sodium metabisulfite.
- b. **Reagent**. Sodium metabisulfite, 2 percent. Add 2 g of sodium metabisulfite to a 100-ml volumetric flask. Dilute to the mark with distilled water. This solution, if stored at 3° or 4°C remains effective for about 1 week.

#### c. Procedure.

- (1) Place 1 drop of capillary (or venous without anticoagulants) blood on a clean glass slide.
  - (2) Add 1 or 2 drops of 2 percent aqueous sodium metabisulfite and mix.
- (3) Place a cover glass on the preparation and express the excess blood by gently pressing the cover glass. The gentle pressure will produce a film thin enough to permit examination of individual red cells. It is not necessary to seal the preparation.

(4) Observe immediately and at intervals of 15 to 30 minutes after preparation for signs of sickle cell information.

#### d. Discussion.

- (1) The sickling phenomenon is a consequence of an inherited, abnormal type of hemoglobin (hemoglobin S). The severity of the anemia as compared to the trait will vary with the proportion of defective to normal hemoglobin. This abnormal hemoglobin can be identified electrophoretically.
- (2) When a fresh solution of sodium rnetabisulfite is used, positive cases should show 10 to 75 percent sickling within 15 minutes. There is a positive correlation between the degree of sickling and the severity of the disease.
  - (3) Normal blood treated in this way shows no sickling.

# 6-13. HEMOGLOBIN ELECTROHORESIS (CELLULOSE ACETATE)

a. **Principle**. Hemoglobin fractions are separated by the rate of their protein migration in an electrical medium. The fractions are stained with ponceau S and quantitated on a densitometer. The order of mobility from the cathode toward the anode is  $A_3 > A_1 > F > S-D>C-A_2$ .

### b. Discussion.

- (1)  $A_2$  hemoglobin migrates identically to hemoglobin C. They are distinguished by the quantity present. If this band is 40 percent or more of the total hemoglobin, it is C.  $A_2$  hemoglobin should always be less than 20 percent.
- (2) Two slow-moving, nonhemoglobin components are see using this technique. These fractions are carbonic anhydrases I and II ( $CA_I$  and  $CA_{II}$ )
  - (3) Hemoglobin A<sub>2</sub> is elevated in thalassemia minor.
  - (4) Genotype SS is found in patients with sickle cell anemia.
  - (5) Genotype AS is found in patients with sickle cell trait.
- (6) This method separates hemoglobin  $A_2$  in the presence of hemoglobin S in patients manifesting sickle-thalassemia disease.
- (7) Hemoglobin F is quantitated by the alkali denaturation test because it migrates close to the hemoglobin A<sub>1</sub> fraction on the electrophoretic pattern.
  - (8) Include known A, S, and C controls in each analysis.

#### c. Normal Values

- (1) A<sub>3</sub> Hemoglobin: One and three tenths to 3 percent.
- (2) Genotype: AA.
- (3) F Hemoglobin: Zero to 2 percent (except in infants).

# 6-14. FETAL HEMOGLOBIN (ALKALI DENATURATION TEST)

a. **Principle**. Fetal hemoglobin (HbF) is more resistant to denaturation in alkaline solution than adult hemoglobin (HbA). Alkali converts HbA to alkaline hematin. Alkaline hematin is insoluble and precipitates. HbF is quantitated by measuring the hemoglobin concentration before and after denaturation.

# b. Reagents.

- (1) <u>Drabkin's reagent</u>. Available in pellets or powder in the Federal Supply Catalog.
- (2) <u>Sodium hydroxide (NaOH) 0.10 N</u>. Dilute 5.55 saturated NaOH to 1-liter with distilled water. Mix well and store in a polyethylene bottle at room temperature.
- (3) Ammonium sulfate, saturated. Add 400 g ammonium sulfate, AR, to 500 ml of distilled water. Stir mechanically for 10 minutes. Heat until all the salt has dissolved and filter rapidly through a Whatman number one filter paper. Store in polyethylene bottle.

# c. **Preparation of Hemolysate**. The procedure is as follows:

- (1) Dilute 0.2 ml of hemolysate with 1.8 of Drabkin's reagent and mix gently.
- (2) Place 2.5 ml of 0.1 N NaOH in a 13 X 100 mm test tube. At zero time, add. 0.5 ml of diluted hemolysate and mix rapidly.
- (3) Exactly 2 minutes after addition of the hemolysate, add 2.0 ml of saturated ammonium sulfate. Mix by inverting six times. DO NOT SHAKE!
- (4) Filter through a 7 or 9 cm Whatman number 42 filter paper. If the filtrate is not crystal clear, filter again. This is the fetal hemoglobin solution.
- (5) Prepare a total hemoglobin solution by adding 0.05 ml of diluted hemolysate to 5.0 ml distilled water. Mix well by inversion. DO NOT SHAKE!
- (6) Prepare a fetal hemoglobin blank by adding 0.5 ml of Drabkin's reagent to 4.5 ml of distilled water. Mix as above.

- (7) Prepare a total hemoglobin control by adding 0 05 ml of Drabkin's reagent to 5.0 ml of distilled water. Mix as above.
- (8) Measure the absorbances of the fetal and total hemoglobin solutions against the appropriate blanks.

### d. Calculations.

Absorbance fetal hemoglobin solution

Absorbance total hemoglobin solution

Report percent to nearest tenth (0.1).

### e. Sources of Error.

- (1) If the hemoglobin concentration of the hemolysate is less than 5 g per the sensitivity of the procedure is decreased.
- (2) Poor pipetting technique causes a significant error due to the small volumes involved.
- (3) If the absorbance of the fetal hemoglobin solution is greater than 0.700, dilute 1.0 ml of the fetal hemoglobin solution with 9.0 ml distilled water. Read the absorbance against the total hemoglobin blank. Multiply the calculated result by 10.

### f. Discussion.

- (1) The alkali denaturation test can be performed on bloody rectal discharges from infants. This test is used to determine whether the discharge is due to ingested maternal blood or a gastrointestinal lesion. The procedure is as follows:
  - (a) Add a small amount of bloody discharge to 10 ml of distilled water.
  - (b) Add 0.2 ml 10 percent sodium hydroxide solution.
- (c) If the solution turns pink, it is fetal blood. A muddy brown solution indicates blood of maternal origin.
- (2) Fetal hemoglobin constitutes approximately 80 percent of a newborn infant's hemoglobin. This decreases to approximately 5 percent at 6 months.

(3) Fetal hemoglobin is increased in adults with thalassemia major and minor, hereditary persistence of HbF, erythroleukemia, aplastic anemia, and spherocytosis.

# g. Normal Values.

- (1) Infants (less than 6 months): Up to 80 percent.
- (2) Adults: Zero to two percent.

**Continue with Exercises** 

# **EXERCISES, LESSON 6**

**INSTRUCTIONS**: Answer the following exercises by marking the lettered responses that best answer the exercise, by completing the incomplete statement, or by writing the answer in the space provided at the end of the exercise.

After you have completed all of the exercises, turn to "Solutions to Exercises" at the end of the lesson and check your answers. For each exercise answered incorrectly, reread the material referenced with the solution.

- 1. When blood is centrifuged, the volume percentage occupied by the packed red cells is known as the:
  - a. Erythrocyte sedimentation rate (ESR).
  - b. Hematocrit.
  - c. Hemoglobin concentration.
  - d. Mean corpuscular volume.
- 2. The procedure for determining hematocrit is performed by:
  - a. Filling a capillary tube capillary with blood.
  - b. Centrifuging at constant speed for a constant period of time.
  - c. Measuring the packed-cell volume.
  - d. Automated sequential analyzers (the results are usually a calculated value).
  - e. All of the above.
  - f. a, c, and d.

•	т.	and the second and the first second of the first sector fields of the second second fields and the second s			
3.	In	ne hematocrit is the most useful single index in determining the degree of:			
	a.	Anemia.			
	b.	Hypochromia or anemia.			
	C.	Leukopenia.			
	d.	Thrombocytopenia or thrombocvtosis.			
4.	The	e hematocrit error rate for determining the degree of polycythemia is:			
	a.	1-3 percent.			
	b.	2-5 percent.			
	C.	2-4 percent.			
	d.	3-6 percent.			
5.		contrast to hematological determinations, what is the percent error rate for the ect red blood cell chamber count?			
	a.	2-4 percent.			
	b.	5-9 percent.			
	C.	6-8 percent.			
	d.	8-10 percent.			
6.	The	e hematocrit values closely parallel the values for the:			
	a.	Packed WBCs and hemoglobin.			
	b.	Packed-cell blood count and reagent.			
	C.	WBC count and hemoglobin.			
	d.	Hemoglobin and RBC count.			

7.	Which	microhe	ematocrit	princi	ple is	correct?

- a. A capillary tube is filled with plasma by capillary action to within 1 to 2 cm of the end. The unfilled end is sealed and the tube is centrifuged. After centrifugation, the capillary tube is placed in a reading device and the hematocrit value determined.
- b. A capillary tube is filled with whole blood by capillary action to within 1 to 2 cm of the end. The unfilled end is sealed and the tube is centrifuged. After centrifugation, the capillary tube is placed in a reading device and the hematocrit value determined.
- c. A capillary tube is filled with whole blood by capillary action to within 2 to 4 cm of the end. The unfilled end is sealed and the tube is centrifuged. After centrifugation, the capillary tube is placed in a reading device and the hematocrit value determined.
- d. A capillary tube is filled with whole blood by capillary action to within 1 to 2 cm of the end. The filled end is sealed and the tube is centrifuged. After centrifugation, the capillary tube is read and recorded.

		continugation, the capitally tabe to read and recorded.
8.	Се	ntrifugation for the microhematocrit lasts:
	a.	30 seconds.
	b.	30 minutes.
	C.	1 minute.
	d.	5 minutes.
9.		ring the microhematocrit test, blood without anticoagulant is identified by a parinized capillary tube with a line.
	a.	Green.
	b.	Red.
	C.	Yellow.
	d.	Pink.

- 10. When performing the microhematocrit test, if blood contains anticoagulant, how far up should the capillary tube re filled with blood?
  - a. One forth.
  - b. Halfway.
  - c. Three fourths.
  - d. Completely.
- 11. Measuring the microhematocrit test, when blood is allowed to enter the two capillary tubes to approximately 2/3 full and air bubbles appear, what does this signify?
  - a. A poor technique was used but it does not affect the results of the test.
  - b. The heparinized capillary tube line was passed.
  - c. The tubes were dirty.
  - d. The seal was broken.
- 12. At what rpm and for how long are the 2 hematocrit tubes centrifuged for the microhematocrit test?
  - a. 4,000 rpm; 2 minutes.
  - b. 5,000 rpm; 4 minutes.
  - c. 10,000 rpm; 5 minutes.
  - d. 15,000 rpm; 7 minutes.

13.		When using the microhematocrit reader, the results should agree within If they do not, then what should occur?			
	a.	1 percent; nothing.			
	b.	2 percent; nothing.			
	C.	2 percent; repeat the procedure.			
	d.	3 percent; repeat the procedure.			
14.	He	matocrit results computed with an electronic cell counter are generally than a manual hematocrit.			
	a.	1 per cent lower.			
	b.	4-6 percent higher.			
	C.	1-4 percent lower.			
	d.	1-3 percent higher.			
15.		e rate at which red blood cells fall when anticoagulated whole blood is allowed stand is known as:			
15.		<u> </u>			
15.	to	stand is known as:			
15.	to s	etand is known as:  Plasma composition.			
15.	to s a. b.	Plasma composition.  Erythrocyte sedimentation.			
15. 16.	to: a. b. c. d.	Plasma composition.  Erythrocyte sedimentation.  Coulter models.			
	to: a. b. c. d.	Plasma composition.  Erythrocyte sedimentation.  Coulter models.  Spherocytosis.			
	to s a. b. c. d.	Plasma composition.  Erythrocyte sedimentation.  Coulter models.  Spherocytosis.  Atthrocyte sedimentation is retarded when the hematocrit exceeds:			
	to s a. b. c. d.	Plasma composition.  Erythrocyte sedimentation.  Coulter models.  Spherocytosis.  Atthrocyte sedimentation is retarded when the hematocrit exceeds:  35 percent.			

17.		The normal hematocrit readings for adult males and adult females are respectively:			
	a.	38-47 percent and 34-41 percent.			
	b.	44-64 percent and 34-41 percent.			
	C.	44-64 percent and 40-54 percent.			
	d.	42-52 percent and 36-46 percent.			
18.		e and shape of the erythrocytes cause the erythrocyte sedimentation rate SR) to:			
	a.	Increase.			
	b.	Remain the same.			
	C.	Decrease.			
	d.	Fluctuate.			
19.	ma	ring erythrocyte sedimentation and with certain diseases, what kind of formation by form in the plasma protein if the plasma protein fibrinogen and globulin are ered?			
	a.	Round.			
	b.	Long.			
	C.	Rouleaux.			
	d.	Short.			
20.		nat happens to the mass of plasma and to the sedimentation rate when the sma protein of fibrinogen and globulin are altered?			
	a.	Mass enlarges; rate increases.			
	b.	Mass decreases: rate decreases.			

c. Mass shrivels; rate increases.

d. Mass enlarges; rate decreases.

- 21. Keeping in mind the mechanical and technical factors, why is it important that the ESR tube be exactly perpendicular?
  - a. A tilt of 50 can cause errors up to 55 percent.
  - b. A tilt of 30 can cause errors up to 30 percent.
  - c. A tilt of 10 can cause errors up to o. 05 percent.
  - d. There are no other factors that effect the ESR tube.
- 22. What other mechanical and technical factors are important and why when working with the ESR tube or holding rack?
  - a. Spilling the ESR tube or tilting holding rack can affect the ESR, as can large changes in temperature.
  - b. Static electricity or movement of the ESR tube or holding rack can affect the ESR, as can large changes in temperatube.
  - c. Vibration or movement of the ESR tube or holding rack can affect the ESR as can large changes in temperature.
  - d. There are no other factors.
- 23. When determining the ESR, using the Wintrobe-Landsberg method, what happens to the anticoagulated blood and how is this procedure measured?
  - a. The anticoagulated blood is placed in a narrow tube/
  - b. The blood cells settle out of the suspension, leaving clear plasma above them.
  - c. The distance that the erythrocytes fall within a given interval of time is measured.
  - d. a, b, and c all happen.
  - e. a and d happen.
  - None of the above.

- 24. With the Wintrobe-Landsberg method, which tube is used to draw 5 ml of blood by venipuncture?
  - a. Green top vacuum tube.
  - b. Lavender top vacuum tube.
  - c. Red lined tube.
  - d. Blue lined tube.
- 25. After the Wintrobe tube is placed in a rack in an exactly vertical position and the time and room temperature are noted, when is a reading taken and what is observed?
  - a. At the end of exactly 1 hour, read the level to which the red cells have settled on the descending scale etched on the tube.
  - b. At the end of exactly 2 hour, read the level to which the red cells have settled on the descending scale etched on the tube.
  - c. Within 15 minute, read the level to which the red cells have settled on the ascending scale etched on the tube.
  - d. In 5 minutes, read the level to which the red cells have settled on the descending scale etched on the tube.
- 26. If measurement of the ESR is delayed more than 2 hours after blood collection, the reading may be inaccurate because of a:
  - a. False varied sedimentation rate.
  - b. False low sedimentation on rate.
  - c. False high sedimentation rate.
  - d. Varied sedimentation rate.

27.			andsberg method, to determine _ mark while holding it at		
	a.	0, 30 degrees.			
	b.	0, 45 degrees.			
	C.	5, 10 degrees.			
	d.	10,50 degrees.			
28.			variation from vertical this is a by perce		
	a.	30 percent.			
	b.	30 percent.			
	C.	40 percent.			
	d.	50 percent.			
29.	Wh		dified Westergren method, who sodium chloride and mixed for		
	a.	0.85 percent, 2.			
	b.	0.90 percent, 3.			
	C.	0.95 percent, 4.			
	d.	0.80 percent, 2.			
30.		ing the modified \ ldren?	Westergren method, what is th	e normal value ESR for	
	a.	0-15 mm per ho	our.		
	b.	0-20 mm per ho	our.		
	C.	0-10 mm per ho	our.		
	d.	0-25 mm per ho	our.		

31.	On	ce hemoglobin gives up its oxygen to the tissues, it is known as:
	a.	Methemoglobin.
	b.	Carboxyhemoglobin.
	C.	Cyanmethemoglobin.
	d.	Reduced hemoglobin.
32.	Не	moglobin reacts with oxygen to form:
	a.	Oxyhemoglobin.
	b.	Methemoglobin.
	c.	Cyanmethemoglobin.
	d.	Carboxyhemoglobin.
33.	Wh	ich compound results when methemoglobin combines with the cyanide radical?
	a.	Oxyhemoglobin.
	b.	Sulfhemogobin.
	C.	Cyanmethemoglobin.
	d.	Carboxyhemoglobin.
34.		ferrous iron in hemoglobin is oxidized to the ferric state, which of the following produced?
	a.	Methemoglobin.
	b.	Carboxyhemoglobin.
	C.	Carbaminohemglobin.
	d.	Reduced hemoglobin.

35.	Wh	nich constitutes most of the hemoglobin of a normal adult?
	a.	Hemoglobin F.
	b.	Hemoglobin A <sub>2</sub> .
	C.	Hemoglobin A.
	d.	Hemoglobin S.
36.		nich is normally present in infants of less than 6 months but not normally present adults?
	a.	Hemoglobin A.
	b.	Hemoglobin A <sub>2</sub> .
	C.	Hemoglobin F.
	d.	Hemoglobin S.
37.	Wh	nen hemoglobin combines with oxygen, its iron must be in what state?
	a.	Ferrous.
	b.	Globulin.
	C. =	= Anemic.
	d.	Active.
38.	Но	w many basic ways are there to measure the hemoglobin concentration?
	a.	2.
	b.	3.
	C.	4.
	d.	5.

39.		nich method is the most widely used to measure the hemoglobin concentration blood?
	a.	Gasometric.
	b.	Cyanmethemoglobin.
	C.	Chemical.
	d.	Specific gravity.
40.		nat does the spectrophotometer's 540 mm wavelength measure during the moglobin reaction using the cyanmethemoglobin method?
	a.	Specific gravity.
	b.	Proportionalism.
	C.	Color intensity.
	d.	Concentration.
41.		though the cyanmethemoglobin method is accurate, what is a disadvantage of ng it?
	a.	It is not the most direct method.
	b.	If the cyanide compounds are handled incorrectly, they can be hazardous.
	C.	Venous samples give erratic values.
	d.	Its hemoglobin pigments are not stable.
42.	The	e normal concentration of hemoglobin in blood of the adult male is:
	a.	10-15 g/dl.
	b.	12-16 g/dl.
	C.	14-17 g/dl.
	d.	18-27 g/dl.

43.	Reduced hemoglobin S is insoluble in:		
	a.	Urea.	
	b.	Anticoagulated blood.	
	C.	Cyanmethemoglobin standards.	
	d.	Phosphate buffer solution.	
44.	What may happen to tile sodium dithionite reagent if it becomes damp?		
	a.	Dissolve.	
	b.	Enlarge.	
	C.	Ignite.	
	d.	Remain the same.	
45. Which		ich cells sickle more rapidly than AS cells?	
	a.	SS cells.	
	b.	SC cells.	
	C.	SD cells.	
	d.	a, b, and c.	
	e.	a and c.	
46.	Erythrocytes of persons with sickle cell anemia or trait will assume a sickle shape when:		
	a.	The oxygen tension is lowered.	
	b.	The oxygen tension is raised.	
	C.	An electrophoretic pattern is run.	
	d.	Highly oxygenated blood is observed.	

6-28

MD0853

47.	Sic	ickle cell anemia is caused by:	
	a.	Endocrine disorders.	
	b.	Massive hemorrhage.	
	C.	Chronic hemorrhage.	
	d.	An inherited protein abnormality of hemoglobin.	
48.		Which reagent is added and mixed with the blood to determine if sickle cells are present?	
	a.	Sodium rnetabisulfite, 2 percent.	
	b.	Diathionite.	
	C.	Sodium metabisulfite, 9 percent.	
	d.	Acetic acid. 4 percent.	
49.		In hemoglobin electrophoresis, how is hemoglobin $A_2$ differentiated from hemoglobin $C$ ?	
	a.	Quantity present.	
	b.	Color after staining.	
	C.	Degree of mobility.	
	d.	Symptoms of patient.	
50.	Dith	Dithionite is:	
	a.	Totally stable.	
	b.	Very stable.	
	C.	Not stable.	
	d.	Limited in stability.	
Check Your Answers on Next Page			

## **SOLUTIONS TO EXCERCISES, LESSON 6**

- 1. b (para 6-1a)
- 2. e (para 6-1a)
- 3. a (para 6-lb)
- 4. c (para 6-lb)
- 5. d (para 6-lb)
- 6. d (para 6-lb)
- 7. b (para 6-2a)
- 8. d (para 6-2b(6))
- 9. b (para 6-2b(1))
- 10. d (para 6-2b(1))
- 11. a (para 6-2b(2))
- 12. c (para 6-2b(6))
- 13. c (para 6-2b(7))
- 14. d (para 6-2d(3))
- 15. b (para 6-3a)
- 16. d (para 6-3b)
- 17. d (para 6-2e)
- 18. d (para 6-3b)
- 19. c (para 6-3c)
- 20. a (para 6-3c)
- 21. b (para 6-3d)
- 22. c (para 6-3d)

- 23. d (para 6-4a)
- 24. b (para 6-4b(1))
- 25. a (para 6-4b(5))
- 26. b (para 6-4c(2))
- 27. b (para 6-4b(3))
- 28. a (para 6-4c(4))
- 29. a (para 6-5a,b(4))
- 30. c (para 6-Sd(3))
- 31. d (~as 6-6,6-7a)
- 32. a (para 6-7a)
- 33. c (para 6-7e)
- 34. a (para 6-7c)
- 35. c (para 6-8)
- 36. c (para 6-8)
- 37. a (para 6-7a)
- 38. c (para 6-9)
- 39. b (para 6-9)
- 40. c (para 6-10a)
- 41. b (para 6-10b(2))
- 42. c (para 6-10c(3))
- 43. d (para 6-11a)
- 44. c (para 6 -11b (2) WARNING)
- 45. d (para 6-11e(2))

- 46. a (para 6-12a)
- 47. d (paras 6-8,12d(1))
- 48. a (para 6-12b)
- 49. a (para 6-13b(1))
- 50. d (para 6-1lb(2))

**End of Lesson 6** 

## **APPENDIX**

## **GLOSSARY OF TERMS**

**Agranulocyte**: A leukocyte without definite cytoplasmic granules.

**Agranulocytosis**: Complete or nearly complete absence of the granular leukocytes from the blood and bone marrow.

**Aleukemic Leukemia**: A fatal condition of the blood-forming tissues, characterized by marked proliferation of immature cells in the bone marrow, without their presence, in any great numbers, in the blood steam.

**Anemia**: A condition in which the blood is deficient in quantity or quality of erythrocytes.

**Anisocytosis**: Variation in size of the erythrocytes.

**Anomaly**: Abnormality.

**Anoxemia**: Lack of normal proportion of oxygen in the blood.

**Antecubital Space**: The area on the forearm frontal to the elbow.

**Anticoagulant**: A substance that prevents the coagulation of blood. Commonly used ones are potassium oxalate, sodium oxalate, sodium citrate, EDTA and heparin.

**Aplasia**: Incomplete or defective blood development; cessation of blood cell formation.

**Aplastic Anemia**: Anemia characterized b incomplete or effective blood development.

**Asynchronous**: Uncoordinated development as in abnormal cell development.

**Azurophilic Granule**: Rounded, discrete, reddish-purple granule, smaller than the granules of neutrophils; 1-10 are common in lymphocytes, and they are very numerous, and smaller, in the cytoplasm of monocytes

**Band Form**: In the Schilling classification, a neutrophil with the nucleus unsegmented and ribbonlike; also stab, staff, nonfilamented.

**Basket Cell**: A degenerated primitive cell which has ruptured and in which the cell nucleus appears as a pale staining smear without prescribed form or shape.

**Basopenia**: An abnormal decrease in the number of basophils.

**Basophil**: A granular leukocyte, the granules of which have affinities for the basic dye of Wright stain (methylene blue). The granules are large, irregular and blue-black in color.

Basophilia: An abnormal increase in the number of basophils.

**Basophilic**: Staining readily with basic dyes, for example, blue with Ramanovsky type stains.

Binary Fission: Simple cell division.

**Bleeding Time**: The time required for a small standardized wound, made in the capillary bed of the finger or ear lobe, to stop bleeding.

**Blood Dyscrasia**: A disease of the blood or blood-forming organs.

**Buffy Coat**: The layer of leukocytes that collects immediately above the erythrocytes in sedimented or centrifuged whole blood.

**Cabot's Rings**: Lines in the form of loops or figures-of-eight seen in erythrocytes in severe anemias.

**Centriole**: A minute cell organoid within the centrosome.

**Centrosome**: An area of condensed cytoplasm active in mitosis.

**Chemotaxis**: the phenomenon of movement of leukocytes caused by a chemical influence.

**Chromatin**: The more stainable portion of the cell nucleus contains genetic materials.

**Clot Retraction**: The rate and degree of contraction of the blood clot.

**Coagulation Time**: The time required for venous blood, in the absence of all tissue factors, to clot in glass tubes under controlled conditions.

**Cocatalyst**: A substance that works in tandem with another group of chemicals to accelerate a reaction velocity without being used up in the reaction.

**Color Index**: The ratio between the amount of hemoglobin and the number of red blood cells.

**Complete Blood Count**: A hematology study which consists of a red cell count, white cell count, hematocrit, hemoglobin, and blood smear study including differential white cell count.

**Congenital**: Born with a person; existing at or before birth.

**Cooley's Anemia (Mediterranean Disease or Thalassemia)**: A chronic progressive anemia commencing early in life and characterized by many normoblasts in the blood, unusual facies, splenomegaly and familial and racial incidence. Target type red blood cells are often present in the peripheral blood.

**Crenation**: The scalloped or notched appearance of the periphery of erythrocytes found when the cells are suspended in a hypertonic solution. Also found in smears, caused by dirty glassware, slow drying, and poor smearing technique.

**Cytoplasm**: Protoplasm of a cell excluding the nucleus.

**DNA**: Deoxyribonucletic acid.

**Differential Count**: An en1.meration of the types of white blood cells seen on a stained blood smear

**Discrete**: Separate.

**Dyscrasia**: Abnormality.

**Ecchymosis**: Subcutaneous extravastion *of* blood covering a large area.

**Endothelial Leukocyte**: Monocyte.

**Eosinopenia**: An abnormal decrease in eosinophils.

**Eosinophil**: A granular leukocyte, the granules of which have an affinity for the acid dye of Wright's stain (eosin). The granules are large, round, uniform in size, red-orange in color and are shiny and refractile.

**Eosinophilia**: A relative or absolute leukocytosis in which the main increase is in eosinophils.

**Eosinophilic**: Readily stained with eaosin, red-orange stain.

**Epigastric**: Pertaining to the upper middle portion of the abdomen.

**Erythremia**: A disease marked by persistent polycythemia and increased blood volume also polychythemia vera.

**Erythrocyte**: Red blood cell.

**Erythrocytosis**: An increase in the total number of erythrocytes.

**Erythrogenic**: Producing erythrocytes.

**Erythroleukemia**: An abnormal condition characterized by proliferation of erythroblastic and myeloblastic cells.

**Erythropenia**: A decrease in the number of red cells in the blood.

**Erythropoiesis**: The production of erythrocytes.

**Etiology**: The theory of the causation of a disease.

**Extravascular**: Occurring outside the blood vessels.

**Extrinsic**: Originating outside of the particular area.

**Fibril**: A microscopic filament often composed of fibrin.

**Fibrin**: The end product of the clotting mechanism that forms a network of fibers that enmesh the formed elements of blood.

**Fibrinogen**: The precursor of fibrin that is present normally in the plasma and produced by the liver.

**Fragility Test (Osmotic)**: A test devised to measure the resistance of the erythrocytes to break down (hemolyze) when subjected to varying concentrations of hypotonic salt solutions.

**Fulminating**: Sudden and severe.

**Golgi Apparatus**: A meshwork of lipid containing fibrils within the cytoplasmic portion of a cell.

**Granulocyte**: A white blood cell that contains specific cytoplasmic granules (neutrophils, eosinophils, and basophils); these granules are peroxidase positive.

**Granulocytosis**: The presence of increased numbers of granulocytes in the blood.

**Granulocytopenia (Granulopenia)**: A decrease in the number of granulocytes in the blood.

**Granulopoiesis**: The production of granulocytes.

**Hemacytometer**: A calibrated chamber in which blood cells are counted.

**Hematin**: A brown or blue-black amorphous iron substance that unites with globin and forms hemoglobin.

**Hematocrit**: The packed cell volume (PVC) of red blood cells obtained by globin and forms hemoglobin.

**Hematology**: The branch of medicine that deals with the study of blood cells, blood producing organs and the manner in which these cells and organs are affected in disease

**Hematoma**: Subcutaneous effusion of blood with resulting swelling, pain, and discoloration, forming a tumorlike mass.

Hematopoietic (Hemopoietic): Blood forming.

**Hemoglobin**: The coloring matter of the red blood cells. A complex iron-bearing pigment that carries oxygen and carbon dioxide.

**Hemoglobinuria**: The presence of free hemoglobin in the urine.

**Hemogram**: The blood picture.

**Hemolysis**: The dissolution or dissolving of the erythrocytes.

**Hemolytic Anemia**: That type of anemia characterized by excessive intra-vascular destruction of red cells.

**Hemophilia**: A hereditary disease characterized by a prolonged coagulation time and repeated hemorrhages, occurring only in males and transmitted only by females and affected males. The cause is a deficiency in a plasma factor (antihemophilic globulin or thromboplastinogen) resulting in a defect in thromboplastic activity.

**Hemoptysis**: The spitting of blood; coughing up blood.

**Hemostasis**: The checking of the flow of blood, especially from a vessel.

**Hepatic**: Originating from the liver.

**Heterozygous**: Derived from germ cells unlike in respect to one or more factors.

**Hemeostasis**: Stability in normal body states.

**Homozygous**: Derived from germ cells that are alike.

**Howell-Jolly Bodies**: Small basophilic particles sometimes found in erythrocytes, remnants of nuclear material.

**Hygroscopic**: Readily taking up and retaining water.

**Hyperplasia**: An increase in cell formation.

**Hypertonic**: Greater than isotonic concentration.

**Hypertrophy**: Enlargement of an organ or part due to increase in the size of the constituent cells.

**Hypochromia**: A decrease in color of the erythrocytes, hence a decrease in their hemoglobin content.

**Hypoplasia**: A decrease in cell formation.

**Hypotonic**: Less than isotonic concentration.

Idiopathic: Disease of unknown cause.

**Inclusion**: Usually lifeless, an accumulation of fats, proteins, crystals pigments or secretory granules within a cell cytoplasm.

**Inhibitor**: A substance, directed against a coagulation factor or factors, which interferes with the coagulation process.

**Intravascular**: Occurring within the blood vessels.

**Intrinsic**: Situated within the particular part.

**In Vitro**: Within a test tube (glass, etc.).

**In Vivo**: Within the living organism, as in life.

**Isotonic**: Solutions with the same osmotic pressure.

**Jaundice**: Yellow mass of the skin and eyes due to the presence of blood pigments in the blood; follows excessive destruction of the blood, obstruction of the bile passage, diffuse liver disease, certain infections, toxic chemical agents and drugs.

**Juvenile Cell**: In the Schilling classification, the cell between the myelocyte and band forms; also metamyelocyte.

**Karyolysis**: Apparent destruction of the nucleus of a cell.

**Karyorrhexis**: Fragmentation of the nucleus; a degenerative process usually followed by karyolysis.

**L.E. Cell**: A large segmented neutrophil or eosinophil that contains ingested autolyzed nuclear fragments in its cytoplasm.

**Leukemia**: An ultimately fatal disease of the blood-forming organs characterized by increased numbers of leukocytes and associated anemia.

**Leukemoid Crisis or Reaction**: A temporary appearance of immature leukocytes in the blood stream, with a marked increase in the total white count. In the laboratory sometimes temporarily indistinguishable from leukemia.

Leukocyte: White blood cell.

**Leukocytosis**: An increase in leukocytes in the blood.

**Leukopenia**: A reduction in the number of leukocytes in the blood.

**Leukopoiesis**: Leukocyte formation.

**Lymphoblast**: The parent cell of the lymphocytic series.

**Lymphocyte**: A white blood cell having a round or oval nucleus and sky blue cytoplasm. The nuclear chromatin is densely clumped but separated by many clear areas giving a "hill and valley" effect. A few red-purple (azurophilic) granules may be present in the cytoplasm.

**Lymphocytosis**: A relative or absolute increase in the number of circulating lymphocytes.

**Lymphopenia**: An abnormal decrease in the number of lymphocytes.

**Lysis**: Destruction by a specific agent.

**Macrocyte**: An erythrocyte larger than normal.

**Macrocytosis**: An increase in the number of macrocytes.

Mast Cell: A basophil or a true tissue cell.

**Maturation Factor**: A substance that will cause cells to ripen and care to maturity.

**Mean Corpuscular Hemoglobin (MCH)**: The average amount of hemoglobin in the red blood cell.

**Mean Corpuscular Hemoglobin Concentration (MHC)**: The average percent hemoglobin saturation in the red blood cell.

**Mean Corpuscular Volume (MCV)**: The volume of the average red blood cell.

**Megakaryoblast**: The parent cell of the megakaryocytic series.

**Megakaryocyte**: An extremely large cell with an irregular lobed, ring or doughnut-shaped nucleus that stains blue-purple. The cytoplasm is abundant, light blue and is packed with fine azurophilic granules. This cell gives rise to thrombocytes.

**Megaloblast**: The type of red cell precursor found in pernicious anemia. This differs from the normal erythrocyte precursor (normoblast) in that the megaloblast is larger and the nuclear chromatin has a fine meshwork or scroll design.

**M.E. Ratio**: The ratio of myeloid to erythroid cells in the bone marrow.

**Mesentery**: The fold of peritoneum that attaches the intestine to the posterior abdominal wall.

Metamyelocyte: Juvenile cell of Schilling.

**Metarubricyte:** An erythrocyte with a pyknotic, contracted nucleus. Also called orthochromatophilic normoblast.

**Methemoglobin**: A spectroscopically detected compound of hemoglobin found in nitrobenzol, and other poisonings. The blood is a chocolate brC1.oln color to the eye.

**Microcyte**: An erythrocyte smaller than normal.

**Microcytosis**: An increase in the number of microcytes

**Micron**: One-thousandth of a millimeter, the common unit of microscopic measure.

**Mitochondria**: Granular components of a cell cytoplasm active in oxidative processes.

**Mitosis**: A series of changes through which the nucleus passes in indirect cell division. A tissue showing many cells in mitosis indicates rapid gro.&1th of that tissue.

**Monoblast**: The parent cell of the monocytic series.

**Monocyte**: A large white blood cell with a pale blue-gray cytoplasm containing fine azurophilic granules. The nucleus is spongy and lobulated.

**Monocytosis**: A relative or absolute increase in the number of circulating monocytes.

Mucosa: Mucous membrane.

**Myeloblast**: The parent cell of the granulocytic or myelocytic series.

**Myelocyte**: The stage in development of the granulocytic series that is characterized by the first appearance of specific granules (eosinophilic, neutrophilic or basophilic) and a round nucleus.

**Myeloid Cells**: The granular leukocytes and their stem cells.

**Myelopoiesis**: Formation of bone marrow and the blood cells that originate in the bone marrow.

**Myeloproliferative**: Rapid production of bone marrow constituents.

**Necrosis**: The death of a circumscribed portion of tissue. Simple necrosis is degeneration of the cytoplasm and nucleus without change in the gross appearance of the tissue.

**Neutropenia**: A decrease in the number of neutrophils in the blood.

**Neutrophil (Polymorphonuclear Neutrophil or Segmented Neutrophil)**: A granulocyte having fine neutrophilic (pink-violet) granules in the cytoplasm. The nucleus is divided into two or more lobes; each lobe is usually connected by a filament.

**Neutrophilia**: An increase in neutrophils.

**Normoblast**. The nucleated precursor of the normal red blood cell. Also called a rubriblast.

Normocyte (Erythrocyte): A red blood cell of normal size.

**NRBC**: Nucleated red cell, usually a metarubricyte when seen in the peripheral blood smear.

**Nucleolus**: An intranuclear pale blue body, surrounded: by a dense condensation of chromatin.

**Occult Blood**: The presence of blood that cannot be detected except by special chemical tests.

Oligochromemia: A decrease in hemoglobin.

Oligocythemia: A decrease in the number of erythrocytes.

**Organoid**: Structures present in cells resembling organs.

**Ovalocyte**: An elliptical erythrocyte.

**Oxyhemoglobin**: The bright red hemoglobin that is loosely combined with oxygen and found in arterial blood.

**Pancytopenia**: A reduction in all three formed elements of the blood, namely, the erythrocytes, leukocytes and thrombocytes.

**Pathologic Increase (Or Decrease)**: Due to abnormal function or disease, as contrasted to physiological (due to normal body function).

**Pernicious Anemia**: A chronic, macrocytic anemia caused by a defect in production of "intrinsic factor" by the stomach. There is accompanying megaloblastic erythropoiesis, poikilocytosis, granulocytic hypersegmentation, achlorhydia, and neurological disturbances.

**Petichiae**: Small spots on the skin formed by subcutaneous effusion of blood (also purpura and ecchymoses).

**Phagocytosis**: The destruction of organisms and extraneous matter by a process of envelopment and absorption.

**Plasma**: The fluid portion of the blood composed of serum and fibrinogen, obtained when an anticoagulant is used.

**Plasma Cell**: A lymphocyte-like cell with an eccentrically placed deep-staining nucleus. The nuclear chromatin is distributed in a "wheel-spoke" fashion. The cytoplasm is deep blue with a lighter halo about the nucleus.

Platelet: Thrombocyte.

Poikilocyte: A red blood cell having abnormal shape (pear-shape, sickle-shape, etc.).

**Poikilocytosis**: Increased number of abnormally shaped erythrocytes.

**Polychromasia**: Diffuse basophilia of the erythrocytes.

**Polychromatophilia**: The presence in the stained blood smear of immature, nonnucleated, bluish-staining red blood cells.

**Polycythemia**: An increase in the total number of erythrocytes. (See erythremia.).

**Precursor**: A substance from which another substance is formed.

**Promyelocyte**: The precursor of the myelocyte having nonspecific azurophilic (red-purple) cytoplasmic granules.

**Prorubricyte**: The second stage of development of the red cell.

**Prothrombin**: The inactive precursor of thrombin that is formed in the liver and present normally in the plasma. Its formation depends upon adequate vitamin K.

**Punctate Basophilia**: Small basophilic aggregates in the erythrocytes that stain blue with the basic dye of Wright's stain; also basophilic stippling.

**Purpura**: Small spots on the skin formed by subcutaneous effusion of blood.

**Pyknosis**: A condensation and reduction in size of the cell and its nucleus.

**Reduced Hemoglobin**: A combination of hemoglobin and carbon dioxide that is found in venous blood.

**Reticulocyte**: A red blood cell showing a reticulum or network when stained with vital dyes (for example, brilliant cresyl blue). The stage between the nucleated red cell and the mature erythrocyte.

**Reticulocytosis**: An increase above normal values of reticulocytes *in* peripheral blood.

RNA: Ribonucleic acid.

**Rouleaux Formation**: The arrangement of red cells with their flat surfaces facing, in which they appear as figures resembling stacks of coins.

**Rubricyte**: Polychromatophilic normoblast.

**Sedimentation Rate, Erythrocyte (ESR)**: The rate at which red cells will settle out in their own plasma in a given time under controlled conditions.

**Serum**: The fluid portion of the blood, after clot formation.

**Shift to the Left**: A term used to designate that condition in which the immature forms of the neutrophils are increased above their normal number.

**Shift to the Right**: Increase in mature, pyknotic, and hypersegmented neutrophils.

**Sickle Cell**: A sickle- or crescent-shaped erythrocyte.

**Sickle Cell Anemia**: This is a hereditary and familial form of chronic, hemolytic anemia essentially peculiar to Negroes. It is characterized clinically by symptoms of anemia, joint pains, leg ulcers, acute attacks of abdominal pain and is distinguished hematologically by the presence of distinct hemoglobin, peculiar sickle-shaped and oatshaped red corpuscles, and signs of excessive blood destruction and active blood formation.

**Smudge Cell**: A ruptured white cell; also basket cell, or degenerated cell.

**Spherocyte**. A red blood cell that is more spherical, smaller, darker, and more fragile than normal.

**Stasis**: A stoppage of blood flow.

**Supravital Stain**: A stain of low toxicity that will not cause death to living cells or tissues.

**Synchronous**: Occurring at the sane time and in a regular pattern.

**Target Cell (Leptocyte)**: An abnormal, thin erythrocyte characteristic of Cooley's or Mediterranean anemia.

**Triturate**: To grind together.

**Thrombin**: This is an enzyme formed from prothrombin that converts fibrinogen to fibrin. This is not present in circulating blood.

**Thrombocyte**: A blood platelet.

**Thrombocytopenia**: A decrease in blood platelets; also thrombopenia.

**Thrombocytosis**: An increase in blood platelets.

**Thromboplastin**: The substance that initiates the process of blood clotting. It is released from injured tissue and/or formed by the disintegration of platelets in combination with several plasma factors.

**Thrombopoiesis**: The production of thrombocytes.

**Thrombosis**: Formation of a thrombus, or blood clot.

**Vacuole**: A space or cavity formed in the protoplasm of a cell.

**Venipuncture**: The act of puncturing a vein in order to remove a sample of blood.

**Viscous Metamorphosis**: Friction between molecules resulting in a structural change.

**Vitamin K**: A vitamin constituent of the normal diet requiring bile salts for absorption. The liver in the production of prothrombin utilizes this vitamin.

**Xanthochromia**: A yellowish discoloration, usually associated with spinal fluid.

**End of Appendix**